

6.1.

Sequent™ Meniscal Repair Device

Sequent™ Meniscal Repair Device, Curved Needle, 4 Implants • Dispositif de réparation méniscale Sequent™, aiguille courbe, 4 implants • Dispositivo di riparazione menisco Sequent™, ago curvo, 4 impianti • Sequent™-Meniskus-Refixationsgerät, gebogene Nadel, 4 Implantate • Dispositivo de reparación del menisco Sequent™, aguja curva, 4 implantes • Συσκευή αποκατάστασης μηνίσκου Sequent™, καμπή βελόνα, 4 εμφυτεύματα • Dispositivo de reparação do menisco Sequent™, agulha curva, 4 implantes • Dispositivo de reparo de menisco Sequent™, agulha curva, 4 implantes • Sequent™-meniscusreparatiehulpmiddel, gebogen naald, 4 implantaten • Sequent™ meniskreparationsanordning, buet nål, 4 implantater • Sequent™ anordning för reparation av menisk, böjd nål, 4 implantat • Nivelkierukan Sequent™-korjauslaite, käyrä neula, 4 implanttia • Sequent™ reparasjonsenhet for menisk, buet nål, 4 implantater • Sequent™ 半月板修复设备, 弯针, 4个植入物 • Urządzenie do naprawy łagodki Sequent™, zakrzywiona igła, 4 implanty • Sequent™ Meniskai Onarim Cihazı, Kivnk İğne, 4 Implant • Dispozitiv de reparare meniscală Sequent™, ac curbat, 4 implanturi • جهاز إصلاح الغضروف الهلالي Sequent™ - بابرة منحنية، ذو 4 غرسات • Устройство для реконструкции мениска Sequent™, изогнутая игла, 4 имплантата • Proizvod za korekciju meniska Sequent™, zakrivljena igla, 4 implantata • Prostředek pro reparaci menisku Sequent™, zakřivená jehla, 4 implantáty • Sequent™ meniscusrekonstrukciós eszköz ívelt tűvel és 4 implantáttal • „Sequent™“ menisko gydymo priemonė, lenktą adatą, 4 implantai • Pomůcka na opravu menisku Sequent™, zakřivená ihla, 4 implantáty • Pripomoček za popravilo meniskusa Sequent™, ukrivljena igla, štirijs vsadki

LOT 1374816 2029-03-04 2024-03-05



http://elfu.conmed.com MD STERILE EO

PEEK-CP/IMA Made by Medtronic CE 2797 EC REP MDSS GmbH Schiffgraben 41 D-30175 Hannover, Germany

PRDT™

Made In USA
CONMED
 CONMED CORPORATION
 525 French Road
 Utica, NY 13502-5994 USA

(01)10845854036587
 (17)290304

LMR004C0001 AK I 1R2 A

Sequent™ Meniscal Repair Device



KNEE
PRESERVATION
SYSTEM

6.3 Sequent Meniscal Repair Device is an All-Inside, Stay-Inside system comprised of implants and suture that can create up to six continuous stitches for a knotless repair in any configuration. This technology provides surgeons a meniscal repair technique with the opportunity for a faster repair. It also reduces meniscal trauma and the risk of chondral damage.

Specifications	
Sequent Devices	MR003C * 3 implants, 2 stitches MR004C * 4 implants, 3 stitches MR007C * 7 implants, 6 stitches <i>6.1.</i>
Instrument	SCEC047 Sequent™ Disposable Kit (Suture Cutter, Entry Cannula) SC047D Sequent™ Disposable Suture Cutter EC047 Sequent™ Reusable Entry Cannula
Pull-out Strength	96.66N *(after 500 cycles)
Suture	#0 Hi-Fi®, 28 inches <i>6.1.</i>
Implant	PEEK, 1.3mm diameter x 5.1mm Length <i>6.2.</i>

DEMO PRODUCT INFORMATION

Catalog #	Description
MR007C	Sequent Meniscal Repair Device, 7 implant
1703-350*	Sequent Model
1402-9-8*	Sequent Puck

* Ordered from Sawbones Website: www.sawbones.com Phone number: (206) 463-5551

Sequent Model



Sequent Puck





Knee Preservation System™

Multiple continuous stitches for repair of any size meniscal tear.

Sequent™ Meniscal Repair Device



- 6 Stitch
- 3 Stitch
- 2 Stitch



MULTIPLE
CONTINUOUS STITCHES



ALL-INSIDE
STAY-INSIDE



KNOTLESS
TECHNIQUE



To learn more about these and other innovative products, call 800-448-6506 or visit ConMed.com.