

**EC Design-Examination Certificate**  
**Directive 93/42/EEC Annex II, Section 4**  
**Medical Devices**

**Registration No.:** ID 60149181 0001

**Report No.:** 21257349 003

**Manufacturer:** TERUMO EUROPE N.V.  
Interleuvenlaan 40  
3001 Leuven  
Belgium

**Product Identification:** Catheter, Angiography  
RADIFOCUS OPTITORQUE

(see attachment for products included)  
Replaces certificate, Registration No.: ID 60114894 0001

The Notified Body hereby declares that an examination of the design dossier relating to the listed products has been performed according to Annex II, section 4 of the directive 93/42/EEC and that the design of the devices conforms to the requirements of the abovementioned directive.

**Expiry Date:** 2024-05-26

**Effective Date:** 2020-05-07

**Date:** 2020-05-07



**TÜV Rheinland LGA Products GmbH - Tillystraße 2 - 90431 Nürnberg**

TÜV Rheinland LGA Products GmbH is a Notified Body according to Directive 93/42/EEC concerning medical devices with the identification number 0197.

**TÜV Rheinland  
LGA Products GmbH  
Tillystraße 2, 90431 Nürnberg**

**Attachment to certificate**

**Registration No.:** ID 60149181 0001  
**Report No.:** 21257349 003

**Manufacturer:** TERUMO EUROPE N.V.  
Interleuvenlaan 40  
3001 Leuven  
Belgium

**Radifocus Optitorque**

**Product code system**

(1) Thoracic use

R H - □ □ □ □ □ □ □ □ □ □  
1 2 3 4 5 6 7 8 9 10 11

Position	Indication & Meaning
1-2	Product name RH: Radifocus® Optitorque™
3	Manufacturing site -: TERUMO Europe N.V.
4	Outer diameter of catheter Indication: 4, 5 6 Size (Fr.): 4 (1.40 mm); 5 (1.70 mm); 6 (2.00 mm)
5-8	Tip shape: 4 alphanumerical digits
9	Number of side holes Indication: 0 ~ 9 Number of side holes: 0 - 9
10	Catheter Length Indication : 6 G 7 8 9 0 1 Length (cm): 60 65 70 80 90 100 110
11	Language used for the indications M: Multi-language

Date, 2020-05-07



**TÜV Rheinland  
LGA Products GmbH  
Tillystraße 2, 90431 Nürnberg**

**Attachment to certificate**

**Registration No.:** ID 60149181 0001  
**Report No.:** 21257349 003

**Manufacturer:** TERUMO EUROPE N.V.  
Interleuvenlaan 40  
3001 Leuven  
Belgium

**Radifocus Optitorque**

(2) Visceral & Cerebral use

R H -

1 2 3 4 5 6 7 8 9 10 11

Position	Indication & Meaning
1-2	Product name RH: Radifocus® Optitorque™
3	Production site -: TERUMO Europe N.V.
4	Indication of catheter use A : visceral use B : cerebral use
5-6	Tip shape, including number of side holes: 2 Alphanumerical digits
7	Outer diameter of catheter Indication: 4, 5 Size (Fr.): 4 (1.40 mm); 5 (1.70 mm)
8	Availability of stopcock 1 : Without stopcock
9-10	Catheter Length Indication: 02 0G 07 08 09 10 11 Length (cm): 20 65 70 80 90 100 110
11	Language used for the indications M: Multi-language



**Date, 2020-05-07**

