

Embedding cassettes

In histology, tissue sections are placed in embedding cassettes to be dehydrated and impregnated by Paraffin.


- made of high-quality, technical plastic (POM)
 - resistant against the solvents usually used in laboratories
 - different colours offer the possibility of colour coding (for example for different users, priorities, etc.)
 - with one 45° inclined and 2 vertical marking areas
 - the roughness of the marking areas is optimized to allow labelling by different systems: pencil, permanent markers (like Marienfeld Superior laboratory marker Cat. No. 61 306 03), inkjet or laser printers
 - suitable for all common clamps
 - ready for use. There is no need to break the lids off the cassettes' bases. Upper and base parts can also be supplied separated on demand if ordered in large quantities. Thus all types are also available for printing on cassette printers.
 - the hinge of our cassettes prevents the cassettes from unintended opening – even when they fall down accidentally. Furthermore, our cassettes can be repeatedly opened and securely closed again.
 - approved for in-vitro diagnostic (IVD) applications according to IVD directive 98/79/EC, with CE-mark and batch number for comprehensive information and traceability
- without lid
 - with accurately fitting neps at the back to fasten standardised stainless steel lids
 - entirely perforated bases increase the adhesion of paraffin blocks to the cassettes
 - 500 pieces in a bag, 20 bags in a carton



Cat. No.	Colour	Unit
Standard		
Perforation: 2.3 mm diameter		
52 525 00	light blue	10'000
▶ 52 525 01	yellow	10'000
▶ 52 525 02	light green	10'000
52 525 03	grey	10'000
52 525 04	white	10'000
52 525 07	pink	10'000
52 525 08	beige	10'000
52 525 09	orange	10'000


Embedding cassettes

- base and upper parts are connected by a hinge and supplied assembled
- 500 pieces in a box, 12 boxes in a carton



Cat. No.	Colour	Unit
Universal		
Perforation: 1 x 5 mm		
▶ 52 530 00	light blue	6000
▶ 52 530 01	yellow	6000
▶ 52 530 02	light green	6000
52 530 03	grey	6000
▶ 52 530 04	white	6000
▶ 52 530 07	pink	6000
52 530 08	beige	6000
52 530 09	orange	6000

Cat. No.	Colour	Unit
Biopsy		
Perforation: 1 x 1 mm		
▶ 52 526 00	light blue	6000
▶ 52 526 01	yellow	6000
52 526 02	light green	6000
▶ 52 526 03	grey	6000
▶ 52 526 04	white	6000
52 526 07	pink	6000
52 526 08	beige	6000
52 526 09	orange	6000



Cat. No.	Colour	Unit
BioNet		
Perforation: 0.35 x 0.35 mm		
It is not necessary to insert pads, because no tissue sections will be lost. This saves time and money.		
52 531 00	light blue	6000
52 531 01	yellow	6000
52 531 02	light green	6000
52 531 03	grey	6000
52 531 04	white	6000
52 531 07	pink	6000
52 531 08	beige	6000
52 531 09	orange	6000



- ▶ This arrow indicates products that are usually available on short notice.