

March 17, 2023

UAB VITROLAB  
Balty Pr. 36 -11,  
KAUNAS, 48196  
LITHUANIA

Dear Mr Raimundas Kaleda,

Thank you for contacting Ortho Clinical Diagnostics regarding your question about Quality Control (QC) that can be used in combination with ORTHO™ Sera Anti-D (DVI) Reagent (Product Code 6904623). You would like to confirm that ORTHO CONFIDENCE™ Quality Control Kit (Product Code 6842785) can be used as a positive and negative QC. This information is needed for a public tender in Lithuania. Your inquiry was registered under the reference number 32203013.

As per product Instructions for Use, ORTHO Sera Anti-D (DVI) will detect D category DVI cells.

ORTHO CONFIDENCE WB Quality Control Kit is prepared from red blood cells collected from blood donors. As per product Instructions for Use:

- Vial 1 group A1B R1R2
- Vial 2 group O R1R1
- Vial 3 group A R2R2
- Vial 4 group B rr

ORTHO Sera Anti-D (DVI) will react positively with ORTHO CONFIDENCE WB Quality Control Kit red blood cells from vial1, 2, 3 (positive QC) and negatively with ORTHO CONFIDENCE WB Quality Control Kit vial 4 red blood cells (negative QC).

This Quality Control WB is suitable for usage on ORTHO VISION®, ORTHO VISION® Max, and BioVue® manual method (ORTHO® Workstation).

It is our highest priority to maintain customer satisfaction. If you have additional questions, please contact Ortho Care™ Technical Solutions Center at +33 3 88 65 47 21.

Sincerely,

Malgorzata Wierzbicka  
Service Engineer Transfusion Medicine