

# Angio-Seal™

Vascular Closure Device

76 pirkimo dalis


## INDICATIONS TO USE

The Angio-Seal device is indicated for use in closing and reducing time to hemostasis at the femoral arterial puncture site in patients who have undergone diagnostic angiography procedures or interventional procedures using an 8 French or smaller procedural sheath for the 8F Angio-Seal device and a 6 French or smaller procedural sheath for the 6F Angio-Seal device.


The Angio-Seal device is also indicated for use to allow patients who have undergone diagnostic angiography and/or interventional procedure to ambulate safely as soon as possible after sheath removal and device placement.

## PRODUCT SPECIFICATIONS

- Insertion Sheath
- Arteriotomy Locator
- 6F (2.0 mm) - 0.035 in. (0.89 mm) Guidewire with J-Straightener or 8F (2.7 mm) - 0.038 in. (0.96 mm) Guidewire with J-Straightener
- Angio-Seal consists of 3 bioabsorbable components: anchor, collagen, and suture
- 10 units per box
- Single use, one year shelf life after sterilization
- Sterilized by gamma radiation. Do not re-sterilize
- Available in select markets



Angio-Seal™ STS Plus & VIP



Angio-Seal™ Evolution

Now, by **TERUMO**  
INTERVENTIONAL  
SYSTEMS