


Rapid Rhino anatomy

Platelet aggregation and tamponade working together



Rapid Rhino with Carboxymethylcellulose (CMC) combines platelet aggregation with tamponade.

Available in a range of sizes:

Catalog no.	Epistaxis solution	Size	Units per box
Rapid Rhino Anterior Device Ultra low profile provides for unsurpassed ease of insertion			
RR 450	for small anatomy/pediatrics	4.5 cm	10
RR 550	for anterior adult epistaxis	5.5 cm	10
RR 551	for COPC patients (Airway included)	5.5 cm	10
Rapid Rhino Posterior/Anterior Device Cuff will conform to nasal anatomy and provide gentle tamponade			
12. RR 750	for challenging posterior and anterior epistaxis	7.5 cm	10
RR 751	for posterior and anterior epistaxis including COPD patients (Airway included)	7.5 cm	10
13. RR 900	Rapid Rhino Posterior Device Two independent balloons for severe posterior epistaxis	9 cm	10
Rapid Rhino Bilateral Device Offers simultaneous inflation via single Pilot Cuff			
RR 552	for anterior/bilateral epistaxis ideal for counter packing	5.5 cm	10
RR 752	for anterior/posterior epistaxis	7.5 cm	10



International
ArthroCare ENT
 Skeppsbron 2, 111 30 Stockholm,
 Sweden
 Tel: +46 8 546 172 00
 Fax: +46 8 546 172 39
 info@arthrocare.se

Australia
ArthroCare (Australasia Pty Ltd)
 PO Box 2186, Normanhurst NSW 2076
 Australia
 Tel: +61 2 9479 9200 Fax: +61 2 9479 9299
 Email: service@arthrocare.net.au

Austria
Atlantech Med. Prod. Vertriebs-AG
 Kaufmannstrasse 15, 6020 Innsbruck, Austria
 Tel: +43 512 361 538 Fax: +43 512 361 548

France
ArthroCare France
 9, rue Noirot, 70000 Vesoul, France
 Tel: +33 84 76 81 38 Fax: +33 3 84 76 74 63

Germany
ArthroCare (Deutschland) AG
 Berghauser Str. 62, 42859 Remscheid, Germany
 Tel: +49 2191 93342 0 Fax: +49 2191 93342 49
 info@arthrocare.de

Scandinavia
ArthroCare Scandinavia ApS
 Science Center Skejby
 Brendstrupgårdsvej 21
 8200 Århus N.
 Denmark
 Tel: +45 7023 2850
 Fax: +45 7023 2851
 Email: info@arthrocare.dk

UK
ArthroCare UK Ltd
 St James Business Park
 9, Grimbald Crag Court
 Knaresborough
 North Yorkshire HG5 8QB
 United Kingdom
 Tel: +44 (0) 1423 888806
 Fax: +44 (0) 1423 888816
 Email: info@arthrocareuk.com

USA
ArthroCare ENT
 680 Vaqueros Avenue
 Sunnyvale, CA 94085-3523
 Tel: +1 408-736-0224

Rapid Rhino products are latex free.

Not for distribution in the United States

© 2009 ArthroCare Corporation. ArthroCare, Rapid-Pac and Rapid Rhino are registered trademarks of ArthroCare Corporation.