


# *Run & Run*


Micro Guidewire

Dedicated for Targeted Chemoembolization




PIOLAX

# "Reshape" is a key for successful procedures.



## Shapability

- Same technology in Japanese swords.
- High Shapability by "Ribbon processing".
- High Re-shapability.




## Stainless Core

- Appropriate stiffness, excellent supportability and torquability.

## Optimal Clearance

• Leads to ideal trackability of micro catheters.




## Lubricity

Full coating of Hydrophilic polymer.

## Visibility

Reduces intraoperative stress.




0.014" 0.016" 0.018"  
(Run & Run)

## Stability

Decrease the need of additional guidewires.


## Product Specification



Catalogue Number	Total Length	Tip Shape	Shaft O.D.	Radiopaque Length	Tapered Length
RUN135-16	135cm	Straight	0.41mm (0.016inch)	3cm	40cm
RUN165-16	165cm				
RUN165-16A		Angle(d)			
RUN180-16	180cm	Straight			
RUN180-16A		Angle(d)			

- Ethylene oxide gas sterilization
- Contents : Guidewire, Shaper, Torque Device, Inserter
- For single use only and do not reuse
- Consult instructions for use

## Tip Image



High Shapability / High Re-shapability

Made in Japan

### PIOLAX MEDICAL DEVICES, INC.

[Overseas Operations Group]  
 2F Pacific Marks Shin-Yokohama, 3-18-20 Shin-Yokohama,  
 Kouhoku-ku, Yokohama-shi, Kanagawa 222-0033 Japan  
 Tel : +81 45 548 9888 / Fax : +81 45 476 5757  
<http://www.piolax-md.co.jp/en/>