

GRASPING FORCEPS




Easy Grasper



Alligator teeth jaws



Pelican



4 prongs

Description

EASY GRASPER: Its rat teeth and its notched jaws provide an excellent gripping force.

ALLIGATOR TEETH JAWS: Excellent gripping force and avoids degrading the grip of the foreign body.

PELICAN: Excellent gripping force and avoids «to drop» the foreign body.

4 PRONGS: Effective grip during food or foreign body withdrawal.

Minimum delivery quantity (MDQ): 1
 Primary packaging: pouch
 Secondary packaging: cardboard box



References	Type	Forceps		Catheter		Minimum diameter of working channel (mm)
		Jaws length (mm)	Maximum opening (mm)	Length (cm)	Diameter (mm)	
CJ-PPDR-23-180	Easy Grasper	5.1	6	180	2.3	≥ 2.8
CJ-2418GH	Alligator teeth jaws	5.8	9.0	230	2.3	
CJ-2423GK	Pelican	11.5	15.5	230	2.3	
CJ-PPF-23-230	4 prongs	45	30	230	2.3	

➤ **INDICATIONS**

This product is used in combination with endoscope to remove foreign corps in human body.

➤ **TECHNICAL INFORMATION**

Materials:

Jaws: Stainless steel
 Handle: ABS
 Catheter: Stainless steel + polyethylene sheath

Contains:

Latex:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Bisphenol A:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No
Product of animal origin:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Nanoparticles:	<input type="checkbox"/> Yes <input type="checkbox"/> No
PVC:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Lead:	<input type="checkbox"/> Yes <input type="checkbox"/> No
Phthalate and DEHP:	<input type="checkbox"/> Yes <input checked="" type="checkbox"/> No		



Sterility:

Sterile: Yes No
 Sterilization mode: Ethylene oxide gaz
 Number of use: Single use



Storage:

Storage conditions: Keep at a temperature of between -20°C to 60°C.
 Keep away from light and moisture.
 Use by date: 2 years

Other data:

GMDN code: 36284
 Barcode: No GS1-128 standard HIBC standard
 With the informations: Reference Lot number Use by date

Use:

Correct and safe use: Cf. Instructions for use
 User advices: Instructions for use, precautions and contraindications - Cf. Instructions for use



Manufacturer:

Manufacturer: Changzhou JIUHONG Medical Instrument
 Country: People's Republic of China

CE Mark:

Notified body: TÜV Rheinland, Certofying body N°: 0197
 Standard: European directive 93/42/EEC - Annex II excluding 4
 Rule: Class IIa medical device

IMPORTER

LIFE PARTNERS EUROPE

161 avenue Gallieni – 93170 Bagnolet
 Phone: + 33 1 49 88 01 01
 Website: www.lifeurope.com

Order:

Fax: + 33 1 48 59 70 72
 E-mail: life@lifeurope.com

Vigilance:

Fax: + 33 1 49 88 83 45
 E-mail: quality@lifeurope.com