

# Renma

0.035" RX PTA Balloon Dilatation Catheter

## Specifications

Catheter design:	Rapid Exchange, coaxial polymer/hypotube shaft	Coating:	SiLX <sup>2</sup> coating
Guide wire compatibility:	0.035" (0.89 mm)	Shaft diameter:	Proximal 2.6Fr Distal 5.0-5.2Fr
Usable shaft length:	135 cm	Introducer compatibility:	6Fr (2.01 mm)
Balloon material:	Polyamide (Nylon)	Nominal Pressure (NBP):	6 atm/ 608 kPa
Balloon folding:	3 folded (Ø 4 – 6 mm); 6 folded (Ø 7 – 10 mm)		



## Ordering information

Renma® Order Codes

Balloon diameter	Balloon length				
	20 mm	40 mm	60 mm	80 mm	100 mm
4.0 mm	BD-C40201DN	BD-C40401DN	BD-C40601DN	BD-C40801DN	BD-C40A01DN
5.0 mm	BD-C50201DN	BD-C50401DN	BD-C50601DN	BD-C50801DN	BD-C50A01DN
6.0 mm	BD-C60201DN	BD-C60401DN	BD-C60601DN	BD-C60801DN	BD-C60A01DN
7.0 mm	BD-C70201DN	BD-C70401DN	BD-C70601DN	BD-C70801DN	BD-C70A01DN
8.0 mm	BD-C80201DN	BD-C80401DN	BD-C80601DN	-	-
9.0 mm	BD-C90201DN	BD-C90401DN	BD-C90601DN	-	-
10.0 mm	BD-CA0201DN	BD-CA0401DN	BD-CA0601DN	-	-

Rated Burst Pressure RBP

Balloon length				
20 mm	40 mm	60 mm	80 mm	100 mm
16 atm/1621 kPa		15 atm/ 1520 kPa		
14 atm/1419 kPa		13 atm/ 1317 kPa		
13 atm/1317 kPa		12 atm/ 1216 kPa		
12 atm/1216 kPa		12 atm/ 1216 kPa		
12 atm/1216 kPa		-		
11 atm/1115 kPa		-		
11 atm/1115 kPa		-		

**TERUMO**

**Terumo Corporation**  
44-1, 2-chome, Hatagaya  
Shibuya-ku, Tokyo 151-0072,  
Japan  
Tel.: +81 3 3374 8111  
Fax: +81 3 3374 8196

**European Headquarters**

**Terumo Europe N.V.**  
Interleuvenlaan 40  
3001 Leuven, Belgium  
Tel.: +32 16 38 12 11  
Fax: +32 16 40 02 49

**Terumo Europe  
Interventional Systems  
Liaison Office**  
Immeuble Krebs  
8, rue Jean-Pierre Timbaud  
CS 10720  
78067 Saint Quentin  
en Yvelines Cedex,  
France  
Tel.: +33 139 30 95 75  
Fax: +33 139 30 95 70

**European Offices**

**Terumo Europe N.V.  
Benelux Sales Division**  
Interleuvenlaan 40  
3001 Leuven,  
Belgium  
Tel.: +32 16 39 25 80  
Fax: +32 16 39 25 99  
The Netherlands:  
Tel.: 0800 0220396  
Fax: 0800 0220414

**Laboratoires Terumo France S.A.**  
Parc Ariane - Bât. Uranus  
1, rue Hélène Boucher  
78284 Guyancourt Cedex,  
France  
Tel.: +33 130 96 13 00  
Fax: +33 130 43 60 85

**Terumo Italia S.r.l.**  
Via Simone Martini, 143 / 145  
00142 Roma RM, Italy  
Tel.: +39 0651 96 14 20  
Fax: +39 0650 30 407

**Terumo UK Ltd**  
Tamesis, The Causeway  
Egham, Surrey  
TW20 9AW, United Kingdom  
Tel.: +44 1784 476 200  
Fax: +44 1784 476 201

**Terumo Deutschland GmbH**  
Ludwig-Erhard-Strasse 6  
65760 Eschborn, Germany  
Tel.: +49 6196 80 230  
Fax: +49 6196 80 23 200

**Terumo Deutschland GmbH**  
Zweigniederlassung Spreitenbach  
Bodenackerstrasse 3  
8957 Spreitenbach, Switzerland  
Tel.: +41 56 419 10 10  
Fax: +41 56 419 10 11

**Terumo Europe España SL**  
Avda. Juan Carlos I, Nº 13 - 7ª Planta  
Edificio Torre La Garena  
28806 Alcalá de Henares (Madrid), Spain  
Tel.: +34 9021 01 298  
Fax: +34 9021 01 358

**Terumo Sweden AB**  
Nya Varvet, Byggnad 90  
Sven Källfelts Gata 18  
426 71 Västra Frölunda, Sweden  
Tel.: +46 3174 85 880  
Fax: +46 3174 85 890

**Terumo Denmark  
Filial of Terumo Sweden AB**  
Tel.: +45 7020 93 80  
Fax: +45 7020 94 80

**Terumo Europe N.V.  
Moscow Representative Office**  
BC Northern Tower  
13Fr., 10 Testovskaya Street  
Moscow 123317, Russia  
Tel.: +7 495 988 4740  
Fax: +7 495 988 4739

**Addresses subject to change,  
please consult our Website.**

© Registered Trademark  
Published by Terumo Europe N.V.