

Misago[®]

RX Self-expanding Peripheral Stent


For a lifetime of flexibility



- + Extended durability
- + Precise deployment
- + Easy to use RX design

New Generation RX Self-expanding stent


- + The special zig-zag cell design provides **optimal flexibility and durability** for the best long term patency for the patient
- + The unique delivery system allows **precise and easy stent deployment**
- + The 0.035"/0,89 mm guide wire compatible Rapid Exchange technology **allows for single operator handling**



✓ *3 radiopaque gold markers
at both ends*


✓ *8 cells for 6-8 mm
9 cells for 9-10 mm*

✓ *Nitinol fabrication*



Misago®

Application:
Occlusive/stenotic
lesions in
iliac/femoro-popliteal
arteries



Ergonomic handle



Thumb wheel



1 2
0.035" (0.89mm)
7X40 mm

TERUMO

Durability compared with compression



In-vitro test at Terumo Corporation, Japan

Ordering information

+
new length
120/150mm
+

| Stent Diameter | Stent length | 40 mm | 60 mm | 80 mm | 100 mm | 120 mm | 150 mm |
|---------------------|--------------|--------------|--------------|--------------|--------------|--------------|--------------|
| 6 mm | | SX-VXA0640LN | SX-VXA0660LN | SX-VXA0680LN | SX-VXA06X0LN | SX-VXA06X2LN | SX-VXA06X5LN |
| 7 mm | | SX-VXA0740LN | SX-VXA0760LN | SX-VXA0780LN | SX-VXA07X0LN | SX-VXA07X2LN | SX-VXA07X5LN |
| 8 mm | | SX-VXA0840LN | SX-VXA0860LN | SX-VXA0880LN | SX-VXA08X0LN | - | - |
| 9 mm | | SX-VXA0940LN | SX-VXA0960LN | - | - | - | - |
| 10 mm | | SX-VXA1040LN | SX-VXA1060LN | - | - | - | - |
| Coaxial part length | | 46 cm | 48 cm | 50 cm | 52 cm | 64 cm | 67 cm |

Effective catheter length: 135 cm for all sizes



Terumo Corporation
44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-0072,
Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8196

European Headquarters

Terumo Europe N.V.
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

**Terumo Europe
Interventional Systems
Liaison Office**
Immeuble Krebs
8, rue Jean-Pierre Timbaud
CS 10720
78067 Saint Quentin
en Yvelines Cedex,
France
Tel.: +33 139 30 95 75
Fax: +33 139 30 95 70

European Offices

**Terumo Europe N.V.
Benelux Sales Division**
Interleuvenlaan 40
3001 Leuven,
Belgium
Tel.: +32 16 39 25 80
Fax: +32 16 39 25 99
The Netherlands:
Tel.: 0800 0220396
Fax: 0800 0220414

Laboratoires Terumo France S.A.
Parc Ariane - Bât. Uranus
1, rue Hélène Boucher
78284 Guyancourt Cedex,
France
Tel.: +33 130 96 13 00
Fax: +33 130 43 60 85

Terumo Italia Srl
Via Simone Martini, 143/145
00142 Roma RM,
Italy
Tel.: +39 0651 96 14 20
Fax: +39 0650 30 407

Terumo UK Ltd
Tamesis, The Causeway
Egham, Surrey
TW20 9AW, United Kingdom
Tel.: +44 1784 476 200
Fax: +44 1784 476 201

Terumo Deutschland GmbH
Ludwig-Erhard-Strasse 6
65760 Eschborn, Germany
Tel.: +49 6196 80 230
Fax: +49 6196 80 23 200

Terumo Deutschland GmbH
Zweigniederlassung Spreitenbach
Bodenackerstrasse 3
8957 Spreitenbach, Switzerland
Tel.: +41 56 419 10 10
Fax: +41 56 419 10 11

Terumo Europe España SL
Avda. Juan Carlos I, Nº 13 - 7ª Planta
Edificio Torre La Garena
28806 Alcalá de Henares (Madrid),
Spain
Tel.: +34 9021 01 298
Fax: +34 9021 01 358

Terumo Sweden AB
Nya Varvet, Byggnad 90
Sven Källfelts Gata 18
426 71 Västra Frölunda, Sweden
Tel.: +46 3174 85 880
Fax: +46 3174 85 890

**Terumo Denmark
Filial of Terumo Sweden AB**
Tel.: +45 7020 93 80
Fax: +45 7020 94 80

**Terumo Europe N.V.
Moscow Representative Office**
JSC Northern Tower
13Fr., 10 Testovskaya Street
Moscow 123317, Russia
Tel.: +7 495 988 4740
Fax: +7 495 988 4739

**Addresses subject to change,
please consult our Website.**

® Registered Trademark
Published by Terumo Europe N.V.



IS76GB0707MDHIS/0611JE