




Angiographic catheter for cardiac application

OUTLOOK[®]

The features you want,
right in your hands




4 Fr.

It's simple: Terumo's **Outlook**[®] set new standards for catheter design and performance.

Outstanding features mean superior torque control, minimized damage to blood vessel walls and less invasive angiography procedures, making **Outlook**[®] ideal for routine use.

A high flow rate, flexibility and kink-resistance, **yet** thin and strong




Stiffer Material


- Different types of Polyamide in tube construction mean a stiffer, more durable shaft.
- Superior backup and shape memory.

Comparison of 4 Fr. Shaft Stiffness

(37°C water temp.)




Comparison of 4 Fr. Torque control



Special Double-Wire Braiding Technology

- Specially designed braiding leads to the proximal shaft.
- Braiding creates gradual flexibility, allowing superb kink-resistance and torque control.




Larger Lumen, Higher Flow

- Ultra-thin wall creates large inner diameter, resulting in higher flow rates.
- 0.038" GW (0.97 mm) compatible.
- Maximum flow rate of 15mL/sec.*
- Pressure limit of 1,200psi (8,274kPa).



*110 cm length catheter, at 37° (10.6 cp)

Features up close



Outlook's unique double-braided structure facilitates high kink-resistance, superior backup and flexibility. Though extremely thin, Outlook® has exceptional torque control and pushability.

Single Catheter (Unit per box 5)

Type	Shape code	Length	Side hole	Outer diameter	Inner diameter	Code number
Amplatz left	AL-1	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4AL1000M
	AL-2	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4AL2000M
	AL-3	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4AL3000M
Amplatz right	AR-1	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4AR1000M
	AR-2	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4AR2000M
	AR-JP	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4ARJP00M
Bypass	Internal mammary short tip	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4BP1R00M
	Judkins left	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4BPJL00M
	Judkins right	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4BPJR00M
	Internal mammary (type JP)	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4BPJP00M
Judkins left	JL-3.5	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4JL3500M
	JL-4.0	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4JL4000M
	JL-5.0	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4JL5000M
	JL-6.0	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4JL6000M
Judkins right	JR-3.5	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4JR3500M
	JR-4.0	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4JR4000M
	JR-5.0	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4JR5000M
Radial	Tiger I	100 cm	1	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4TIG110M
	Tiger II 4.0	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4TR4000M
	Tiger II 4.0	100 cm	1	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4TR4010M
	Tiger II 5.0	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4TR5000M
	BLK 4.0	100 cm	1	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4BLK410M
Straight pigtail	-	110 cm	8	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4SP0081M
	-	90 cm	8	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4SP0089M
	-	65 cm	8	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4SP008GM
	-	65 cm	4	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4SPS04GM
Angled pigtail	145° (5.5 cm)	110 cm	8	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4AP4581M
	155° (5.5 cm)	110 cm	8	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4AP5581M
Multipurpose	2.5 cm	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4MP2500M
	4.0 cm	100 cm	2	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4MP4020M
	El Gamal 2	100 cm	0	4 Fr / 1.40 mm	0.041" / 1.05 mm	RQ-4MP4100M

Tri-Pack Standard Items - References (Unit per box: 5)

Code number	Catheter references	Type	Tip curve	Length	Side hole	Fr size
RQ-MK4226M	RQ-4JR4000M	Judkins right	4.0 cm	100 cm	0	4 Fr / 1.40 mm
	RQ-4JL4000M	Judkins left	4.0 cm	100 cm	0	4 Fr / 1.40 mm
	RQ-4SP0081M	Straight pigtail	-	110 cm	8	4 Fr / 1.40 mm
RQ-MK4227M	RQ-4JR4000M	Judkins right	4.0 cm	100 cm	0	4 Fr / 1.40 mm
	RQ-4JL4000M	Judkins left	4.0 cm	100 cm	0	4 Fr / 1.40 mm
	RQ-4AP4581M	Angled pigtail (145°)	5.5 cm	110 cm	8	4 Fr / 1.40 mm

Penta-pack Standard items - References (Unit per box: 4)

Code number	References catheter items	Introducer (Diameter/Length)	Spring Guide Wire (Diameter/Length)	Spring Mini Guide Wire (Diameter/Length)
RQ-4226B4E42	RQ-4JR4000M / RQ-4JL4000M / RQ-4SP0081M	4 Fr / 100 mm	0.035" (0.89 mm) / 150 cm	0.035" (0.89 mm) / 45 cm
RQ-4227B4E42	RQ-4JR4000M / RQ-4JL4000M / RQ-4AP4581M	4 Fr / 100 mm	0.035" (0.89 mm) / 150 cm	0.035" (0.89 mm) / 45 cm

Dilator length = 155 mm.
Spring Guide Wire: teflon coated, length = 150 cm, J curve = 3 mm, diameter as indicated

Packaging: One delivery box contains 4 Penta-Packs.
Note: For any other possible item combination, please contact your Terumo local sales office.

TERUMO CORPORATION
44-1, 2-chome, Hatakeya,
Shibuya-ku,
Tokyo 151-8572 - Japan
Tel: 81 3 3374 8111
Fax: 81 3 3374 8196

TERUMO EUROPE NV
Researchpark Zone 2 - Hausrode
Interleuvenbaan 40
B-3001 Leuven
Belgium
Tel: 32 16 38 12 11
Fax: 32 16 40 02 49

**TERUMO EUROPE
INTERVENTIONAL SYSTEMS**
Pars Ariane - Bât. Jupiter
1, rue Hélène Boucher
78284 Guyancourt Cedex
France
Tel: 33 1 39 30 95 75
Fax: 33 1 39 30 95 70

TERUMO DEUTSCHLAND GmbH
Hauptstrasse 87
D-65760 Eschborn
Germany
Tel: 49 6196 8023-0
Fax: 49 6196 8023-200

**LABORATOIRES TERUMO
FRANCE S.A.**
Parc Ariane - Bât. Uranus
1, rue Hélène Boucher
78284 Guyancourt Cedex
France
Tel: 33 1 39 30 95 13 00
Fax: 33 1 39 30 43 60 85

TERUMO ITALIA SRL
Via Simone Martini, 143-145,
00142 Roma RM
Italy
Tel: 39 06 51 96 14 20
Fax: 39 06 50 30 40 7

**TERUMO EUROPE
ESPANA SL**
Ay. de Juan Carlos I Plaza 7ª
Edificio Torre La Garena
28806 Alcala de Henares
Madrid Spain
Tel: 34 902 10 12 88
Fax: 34 902 10 13 58

TERUMO UK Ltd
Ternolis
The Crayze Way
Egham
Surrey TW20 6AW
United Kingdom
Tel: 44 178 447 6200
Fax: 44 178 447 6201

**TERUMO EUROPE NV
BENELUX SALES DIVISION**
Researchpark Zone 2 - Hausrode
Interleuvenbaan 40
B-3001 Leuven
Belgium
Tel: 32 16 39 25 80 (Belgium)
Fax: 32 16 39 25 99 (Belgium)
Tel: 0800 0220 036 (NL)
Fax: 0800 0220 414 (NL)

TERUMO SWEDEN AB
Nya Varvet, Byggnad 90
Sven Källviks Gata 18
SE-426 71 Västra Frölunda
Sweden
Tel: 46 31 748 58 90
Fax: 46 31 748 58 90

**MOSCOW REPRESENTATIVE
OFFICE OF TERUMO EUROPE NV**
A. Ivana Franko St,
Moscow, Russia, 121109
Tel: 7 495 234 3303
Fax: 7 495 234 3300

TERUMO DEUTSCHLAND GmbH
Zweig Niederlassung Wetzlingen/ Schweiz
Bahnhofstrasse 88
CH 5430 Wetzlingen/ AG
Tel: +41 56 419 10 10
Fax: +41 56 419 10 11

www.terumo-europe.com
® Registered Trademark
Manufactured by Terumo
Corporation in Japan
© TERUMO CORPORATION