



> All models are pre-loaded onto the Ex-PRESS® Glaucoma Filtration Device Delivery System

Ex-PRESS® Glaucoma Filtration Device Delivery System (EDS)



Ex-PRESS® Glaucoma Filtration Device R-50

Beveled tip for easy insertion

External lumen	400µm
<b>Internal lumen size</b>	<b>50µm</b>
Device length	2.96 mm
Tip shape	Beveled
Back plate shape	Uniform
<b>USA Ordering Info (SKU)</b>	<b>R-50 PL-40053</b>



Ex-PRESS® Glaucoma Filtration Device P Model

Vertical Channel for Posterior Flow

External lumen	400µm
<b>Internal lumen size</b>	<b>50µm / 200µm</b>
Device length	2.64 mm
Tip shape	Decreased bevel angle
Back plate shape	Vertical channel
<b>USA Ordering Info (SKU)</b>	<b>P-50 PL-47053</b> <b>P-200 PL-47203</b>



To order, please contact Alcon customer service at **1-800-TO-ALCON (1-800-862-5266)**.

References  
1. Ex-PRESS® Glaucoma Filtration Device Directions for Use, 2. Maris P.J.G., Ishida K., Netland P.A. Comparison of Trabeculectomy with Ex-PRESS® Miniature Glaucoma Device Implanted Under Scleral Flap. *J Glaucoma* 2007;16:14-19.

**Alcon**

©2010 Alcon, Inc. 6/10

EXP10505VS

