

**MONARCH® IOL Delivery System**

The MONARCH® IOL Delivery System combines a reusable titanium handpiece and a sterile single use cartridge for enhanced implantation of the AcrySof® IOL. Advanced design enables the surgeon to view and verify lens orientation throughout the implantation process. Simplified loading, controlled consistent delivery, and ease of implantation are the benefits that MONARCH® IOL Delivery System can provide.



**MONARCH® II IOL Injector** 8065-977771



**MONARCH® III IOL Injector** 8065-977773

**MONARCH® Cartridges**

10/box



**MONARCH® II Loading Forceps** 560.01

The ALCON®/GRIESHABER® MONARCH® II Loading Forceps is for fully controlled handling of the AcrySof® Single-Piece and AcrySof® Natural IOLs from its packaging into the MONARCH® Delivery System.

**MONARCH® Cartridge/AcrySof® Lens Model/Diopter Range**

8065-977757 A Cartridge <sup>1</sup>	8065-977758 B Cartridge <sup>1</sup>	8065-977759 MONARCH® II C Cartridge <sup>1</sup>	8065-977762 MONARCH® III C Cartridge <sup>1</sup>	8065-977763 D Cartridge <sup>2</sup>	8065-977768 P Cartridge <sup>2</sup>				
MA50BM	SN60AT	SN6AT2	SN60AT	6.0-27.0	SN6AT2	6.0-30.0	SN60WF	6.0-27.0	LI2500
MA60MA	SN60T3	SN6AT3	SA60AT	6.0-27.0	SN6AT3	6.0-30.0	SN6AD1	6.0-27.0	LI3000
MN60MA	SN60T4	SN6AT4	SN60WF	6.0-30.0	SN6AT4	6.0-30.0	SN6AD3	6.0-27.0	LI3500
MA60AC	SN60T5	SN6AT5	SN6AD1	6.0-27.0	SN6AT5	6.0-30.0	SN6AT2	6.0-27.0	LI4000
MN60D3	SN60T6 <sup>+</sup>	SN6AT6	SN6AD3	6.0-27.0	SN6AT6	6.0-27.0	SN6AT3	6.0-25.0	
MN60AC	SN60T7 <sup>+</sup>	SN6AT7	SN60T3	6.0-25.0	SN6AT7	6.0-27.0	SN6AT4	6.0-25.0	
MN6AD1	SN60T8 <sup>+</sup>	SN6AT8	SN60T4	6.0-25.0	SN6AT8	6.0-27.0	SN6AT5	6.0-25.0	
	SN60T9 <sup>+</sup>	SN6AT9	SN60T5	6.0-25.0	SN6AT9	6.0-27.0	SN6AT6	6.0-23.0	
	SN6AD3	SN60WF	SN60T6	6.0-25.0	SNDIT2	6.0-30.0	SN6AT7	6.0-23.0	
	MA30AC	MA60AC	SN60T7	6.0-25.0	SNDIT3	6.0-30.0	SN6AT8	6.0-23.0	
	MA30BA	MN60AC	SN60T8	6.0-25.0	SNDIT4	6.0-30.0	SN6AT9	6.0-23.0	
	SA30AL	SNDIT2	SN60T9	6.0-25.0	SNDIT5	6.0-30.0	SN60T3	6.0-23.0	
	SA60AT	SNDIT3					SN60T4	6.0-23.0	
	SN6AD1	SNDIT4					SN60T5	6.0-23.0	
		SNDIT5					SN60T6	6.0-23.0	
							SN60T7	6.0-23.0	
							SN60T8	6.0-23.0	
							SN60T9	6.0-23.0	
							SNDIT2	6.0-25.0	
							SNDIT3	6.0-25.0	
							SNDIT4	6.0-25.0	
							SNDIT5	6.0-25.0	

<sup>1</sup>Validated with the MONARCH® II IOL Delivery System. <sup>2</sup>Validated with the MONARCH® III IOL Delivery System.