


SUPRAPUBIC CATHETER SETS

Suprapubik Kateter Setleri



spiraliniu galu, drenažo angos vidinėje spiralės pusėje.

✓ Polyurethane

✓ Poliūretan


kateterio dydis

kateterio ilgis

Catheter Fr Kateter Fr	Polyurethane-Pigtail Catheter Set Poliūretan-Pigtail Kateter Set	Length Uzunluk
10	539001	65
12	539002	
14	539003	
16	539004	

Suprapubic Catheter Set:
Catheter + Trocar + Urine Bag + Scalpel + Clamp

su kaniule 12cm, šlapimo maišu 1,5L, skalpeliu ir spaustuku



I.V. Plaster

Pleistras sudarytas iš mikroaktytos neaustinės ar permatomos plėvelės, kuri padengta hipoalergine klijų mase. Absorbuojamasis minkštas tamponas su tinklo dangta neprilimpa prie injekcijos vietos ir nesukelia skausmo nuimant tvarstį.
Dydis - 6 x 8 cm


Prie odos limpanti kateterio fiksacijos sistema.

2 p.d.

SUPRAPUBIC CATHETER SETS

Suprapubik Kateter Setleri

Plasti-med[®]
"supporting healthy life"



fiksuojamasis balionėlis, drenažo angos vidinėje pusėje. ✓ Silicone

✓ Silikon

dydis

ilgis

Catheter Fr Kateter Fr	Silicone-Balloon Catheter Set Silikon-Balonlu Kateter Set	Balloon (ml) Balon (ml)	Length Uzunluk
10	539005	5	40
12	539006	5	
14	539007	5	
16	539008	10	

Suprapubic Catheter Set:
Catheter + Trocar + Urine Bag + Scalpel + Clamp

su kaniule 12cm, šlapimo maišu 1,5L, skalpeliu ir spaustuku



FLOCARE® ACCESSORIES

The Flocare® accessories range includes feeding connectors, adaptors, containers, and pump accessories to aid enteral feeding. All Flocare® Feeding Tubes, Flocare® Feeding Sets and ENFit™ enteral syringes are ENFit™ compatible.

Adaptors	Product description	Code	Units per carton
<p>ENFit™ Transition Adaptor</p> 	<p>Designed to be inserted into a non-ENFit™ feeding tube to enable connection to an ENFit™ feeding set.</p> <p>Recommended for single patient use for the lifespan of the feeding tube.</p>	122073	50
<p>Flocare® Transition Step Connector</p> 	<p>Designed to be attached to an ENFit™ feeding set to enable connection to a non-ENFit™ feeding tube.</p> <p>Recommended for 24-hour, single use, single patient only.</p>	89828	30

Containers	Product description	Code	Units per carton
<p>Flocare® Container 0.5L</p> 	<p>Designed to be used to deliver decanted product, reconstituted powder based formula and/or water via an enteral feeding pump or gravity feeding set.</p> <p>Recommended for 24 hour, single use only.</p>	40441	10
<p>Flocare® Container 1L</p> 	<p>Designed to be used to deliver decanted product, reconstituted powder based formula and/or water via an enteral feeding pump or gravity feeding set.</p> <p>Recommended for 24 hour, single use only.</p>	40455	10
<p>Flocare® Top-fill Container 1.3L</p> 	<p>Designed to be used to deliver decanted product, reconstituted powder based formula and/or water via an enteral feeding pump or gravity feeding set.</p> <p>Top-fill design for ease of use.</p> <p>Recommended for 24 hour, single use only.</p>	70139	10

To order, please call:
Nutricia Customer Care 1800 889 480



Dia.



TUBE FOR TREATMENT OF ESOPHAGEAL VARICES

CODE: 9825 - 9827 (pediatric)

atšaka skrandžio levažui, dekompresijai

manžetės orui pripūsti

galas apvalus , 4 skylutės

du latekso balioniai

su rentgenokontrastiniais žymekliais abiejų balionų distaliniame ir proksimaliniame galuose

Composed of: PVC medžiagos vamzdelis su 3 kanalais

- ♦ 1 PVC tube, 3-way, rounded distal end, with multiple openings to irrigation and suction.
- ♦ Gastric latex balloon in connection valve with "stomach" inscription.
- ♦ Esophageal latex balloon in connection valve with "esophagus" inscription.
- ♦ Radiopaque rings at the proximal and distal end of both balloons.
- ♦ Funnel in proximal end.

Diastalinis galas su ertmėmis praplovimui

36.1.2 Proksimalinis galas piltuvėlio formos

STERILE, SINGLE PACK, SINGLE USE.

AVAILABLE SIZES: Ch 14,0 / 16,0 / 18,0 / 21,0.
Ch 14,0 (pediatric). 36.1 Blekmoro zondas 18CH

PHTHALATE FREE


RECTAL TUBE WITH PVC BALLOON

CODE: 9961

Composed of:

- ♦ 1 PVC tube, rounded distal end, side orifice for irrigation and suction.
- ♦ PVC balloon.
- ♦ Funnel in proximal end.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 6,0 / 8,0 / 10,0.



RECTAL TUBE WITH LATEX BALLOON


CODE: 9962


Composed of:

- ♦ 1 PVC tube, rounded distal end, side orifice for irrigation and suction.
- ♦ Latex balloon.
- ♦ Funnel in proximal end.

STERILE, SINGLE PACK, SINGLE USE.
AVAILABLE SIZES: Ch 8,0 / 10,0.









Rectal Tube


- (1) Made of non-toxic, non-irritant PVC (Medical grade) or DEHP free.
- (2) Color: Transparent or translucent (Surface fog form).
- (3) Closed tip with two lateral eyes.
- (4) Connector (color coded)-Size (Fr/Ch):F18, F20, F22, F24, F26, F28, F30, F32, F34, F36.
- (5) Package: Single in PE or blister package.
- (6) Sterilized by EO GAS.



GT020-100

Rektalinis zondas

- (1) Pagamintas iš netoksiško, nedirginančio PVC (medicinės paskirties) ir be DEHP.
- (2) Spalva: skaidrus arba permatomas.
- (3) Uždaras galiukas su dviem šoninėm angelėms.
- (4) Konektorius (koduotas spalva)- **Dydžiai** (Fr/Ch): F18, F20, F22, F24, F26, F28, F30, F32, F34, F36.
- (5) Pakuotė: Individualus įpakavimas PE pakuotėje arba Blister pakuotėje
- (6) Sterilizuota EO duiomis.



Disposable Specimen Retrieval Bag



Features:

- Ripstop Nylon with polymer coating for leak-proof performance.
- Easy deployment.
- Large opening of the bag ensures the specimen capture easier.
- Proximal and distal tab points make the extraction of bag easier.
- Cost effective choice.

Code No.	Bag Size (mm)	Bag Volume (ml)
QW1-S	110*150	150
QW1-M	120*200	350
QW1-L	250*240	1500

43 p.d. Apsauginis morceliavimo maišas

Disposable Morcellation Bag



Code No.	Bag Material	Bag Size (mm)	Introducer Diameter (mm)	Bag Volume (ml)
QW7-1	Nylon	380*350	/	2000
QW7-2	Nylon	450*360	/	4000
QW4-7	TPU	210*135	12	2100

43 p.d.

Instruction For Use

Disposable Endoscopic Retrieval Bags (Morcellation Bag)



Illustration and Nomenclature

1. Skirtas saugiam morceliavimui ginekologinių operacijų metu.
2. Permatomas.
3. Du įėjimai: vienas-laparoskopui, kitas-morceliatoriui.
4. Tūris 2000 cm³.
5. Vienkartinio naudojimo.
6. Įpakuota po 1-ą vienetą.



Read all information carefully

Failure to properly follow the instructions may lead to serious surgical consequences.


Indications

The Tissue Retrieval Bag (Morcellation Bag) is a disposable device used as a receptacle for the collection of tissue specimens during laparoscopic uterine morcellation procedures.


Device Description

Code No.	Opening (Width*Diameter)	Monitoring Path Diameter	Diameter of Introducer Tube	Volume
43 p.d. QW7-1	380mm*350mm	Φ 20mm	/	2000ml <small>4. Tūris 2000 cm³.</small>
QW7-2	450mm*360mm	Φ 20mm	/	4000ml
QW4-7	210mm*135mm	Φ 6mm	12mm	2100ml


Instructions for use



1. Insert the introducer tube into the abdominal cavity through a trocar, and then use the pushing rod to push the retrieval bag into the abdominal cavity. If without introducer, just rolled the bag and insert the bag into the abdominal cavity.
2. Put the excised uterus into the retrieval bag.
3. Remove the trocar, and then pull the mouth/opening of the retrieval bag out of the body through the remaining incision.
4. Pull the mouth of the monitoring terminal of the bag out of the body through the remaining incision of another trocar.



5. Open the mouth of the monitoring terminal and put the lens into the bag through it.
6. Under the monitor, insert the morcellator into the mouth of the retrieval bag, continuously morcellate the uterus and suck it out of the body;



7. When the uterus in the retrieval bag is crushed and sucked out of the body, it is dragged out of the body together with the bag, and the operation is finished.

Warning and Precautions

1. This product is for Single Use only. Single packed. 5. Vienkartinio naudojimo.
6. Įpakuota po 1-ą vieneta.
2. Do not use if Package has been Damaged, Improperly Stored, or past Expiration Date.
3. The volume of the removed tissue must not exceed the maximum capacity of the retrieval bag.
4. To avoid spillage of tissue fluid or tissue debris, do not damage the retrieval bag while removing excised tissues.
5. Make sure that the excised tissues have been taken out before pulling out the tissue bag. If an intact tissue specimen needs to be taken out, and the bag containing the specimen cannot be easily removed, the incision should be carefully enlarged to facilitate the bag retrieval.

Environmental conditions for transport and storage

Temperature: < 54 °C

Relative Humidity: 0% to 80%

【 Symbols】

	Caution
	Sterilized using ethylene oxide
	Consult instruction for use
	Do not re use
	Do not use if package is damaged
	Catalogue number
	Batch code
	Use by date
	Date of manufacture
	Manufacturer
	Authorized representative in the European Community

**Jiangsu BANA Medical Technology Co., Ltd.**

Address: No.28-3, Guihua Road, Zhonglou Economic Development Zone, Changzhou, 213013 Jiangsu, China.

Tel: 0519-85159001

Fax: 0519-85155036

E-mail: info@jsbana.cn

EC	REP
----	-----

Caretechion GmbH

Address: Niederrheinstr. 71, 40474 Duesseldorf, Germany.

Tel: +49 211 3003 6618

Fax: +49 211 3003 6619

E-mail: info@caretechion.de

Date of Issue: Aug. 01, 2021

Revision No.: C/1

Please see product risk information in the IFU at cookmedical.eu.



51 p.d.

Easytrap®



Sistema per aspirazione polipi

SUCTION POLYP TRAP

Meditalia®

Easytrap®



Meditalia®

SISTEMA ASPIRAZIONE POLIPI SUCTION POLYP TRAP

**Alta performance dei tubi durante l'aspirazione.
Contenitore pre-assemblato in posizione 1**

*High performance of the tube during the aspiration.
Pre-assembled container into position 1*

**Rotazione fluida del contenitore.
I diametri dei fori garantiscono un adeguato
drenaggio dei liquidi aspirati**

*Smooth rotation of the container.
The holes diameter guarantee an adequate
draining of the aspirated liquids*

**Il connettore a forma conica allungata filiforme
facilita l'inserimento dell'aspiratore chirurgico**

*Conical shape of the connector facilitates the insertion
into the surgical aspiration unit*

**Dotato di esclusivo tappo per il trasporto
in laboratorio e completo di un' idonea etichetta
per il riconoscimento del paziente e/o file**

*Provided with an exclusive lid for the transport to
the laboratory labelled with
patients name and/or file number*

Products info

References

Chambers

Material

Accessories

MED/105

4

polycarbonate

secon lid

53. P.D.



OptraGate

Optragate burnos žiodikliai

Moisture control made easy

OptraGate is a latex-free lip and cheek retractor. The aid's unique flexibility and elasticity make it comfortable to wear and assist patients to keep their mouth open.

The main product benefits include:

Easy access to the cavity

- Lips and cheeks are evenly retracted around the mouth
- Easily accessible, manageable treatment field
- Effective and quick relative isolation

Relaxed and more efficient treatment procedure

- Easy placement
- Increased concentration on the actual treatment procedure
- Relieves the work of the dental assistant and dentist.

High patient comfort of wear

- 3-dimensional flexibility for a very good positioning
- Soft, flexible material with recesses to accommodate the upper and lower lip frenulum
- Good compatibility thanks to the latex-free material

minkštos ir lanksčios medžiagos

Efficient Esthetics

OptraGate forms a part of the "Direct Restoratives" product category. The products of this category cover the procedure involved in the direct restoration of teeth – from preparation to restoration care. The products are optimally coordinated with each other and enable successful processing and application. You can find more Information about the efficient esthetics workflow [here](#).

skirti naudoti odontologinių procedūrų, restauracijų metu;

OptraGate can be used for a wide range of indications:

OptraGate is available in the sizes Regular, Small and Junior

The newest colour variations in blue and pink increase the age-appropriate motivation of children and therefore improve patient compliance during the treatment.

OpraGate in blue and pink helps to playfully ease the treatment situation right at the beginning. Young patients get the emotional support they need. Satisfied patients establish a good basis for a good reputation and sustained success.

OpraGate in blue and pink is available in the sizes Junior and Small from now on.

- OpraGate White -
 - Assortment (40 x regular, 40 x small)
 - Refills x 80 (regular, small or junior) **trijų dydžių: vaikiški, maži ir reguliarūs**
 - Refills x 40 (regular, small or junior)
- OpraGate Blue & Pink -
 - Assortment Small x 40
 - Assortment Junior x 20