

# Safety-Lancet Product Range

					
<b>Version</b>	<b>Mini</b>	<b>Normal</b>	<b>Extra</b>	<b>Super</b>	<b>Neonatal</b>
<b>Order no.</b>	85.1015	85.1016	85.1017	85.1018	85.1019
<b>Needle size</b>	28 G	21 G	18 G	Blade 1.5 mm	Blade 1.5 mm
<b>Penetration depth</b>	1.6 mm	1.8 mm	1.8 mm	1.6 mm	1.2 mm
<b>Blood volume*</b>	up to 10 µl	10 - 50 µl	20 - 100 µl	100 - 500 µl	50 - 250 µl
<b>Puncture site</b>	finger, ear, heel	finger, ear, heel	finger, ear, heel	finger, ear, heel	heel, finger
<b>Recommended collection tubes</b>	end-to-end-capillary	Microvette® 100	Microvette® 200, CB 200	Microvette® 200, 300, 500, CB 200, Multivette® 600	Microvette® 200, 300, 500, CB 200
<b>Application</b>	Small children, adults: diabetics, oncology, blood sugar determination	Children, adults: haematocrit, blood count	Children, adults, obese adults: capillary blood gas analytics, full blood count	Children, adults with very large blood volumes, newborn: full blood count, reticulocytes, metabolism test	Children, adults with large collection volumes, newborn: blood count, blood sugar, bilirubin

\* Please be aware that the selection of the appropriate lancet depends upon various factors, such as skin type, blood flow to the puncture site and the correct application of the Safety-Lancet onto the skin. In case of doubt, it is recommended to test a selection of different Safety-Lancets prior to using the products for routine applications. The information given above is a recommendation only and does not, under any circumstances replace medical, scientific or technological advice.

# Saugiu-Lancetu Produktu asortimentas



Versijos	Mini	Normal	Extra	Super	Neonatal
<b>Užsakymo Nr.</b>	85.1015	85.1016	85.1017	85.1018	85.1019
<b>Adatos dydis</b>	28 G	21 G	18 G	Blade 1.5 mm	Blade 1.5 mm
<b>Skverbties gylis depth</b>	1.6 mm	1.8 mm	1.8 mm	1.6 mm	1.2 mm
<b>Kraujo tūris *</b>	up to 10 µl	10 - 50 µl	20 - 100 µl	100 - 500 µl	50 - 250 µl
<b>Punktūros vieta</b>	pirštas, ausis, kulnas nuo galo iki galo kapiliaras	pirštas, ausis, kulnas	pirštas, ausis, kulnas	pirštas, ausis, kulnas	pirštas, kulnas
<b>Rekomendacija</b>	Maži vaikai, suaugusieji: cukriniu diabetu, onkologija, kraujo cukraus nustatymas	Vaikai, suaugusieji: hematokritas, kraujo tyrimas	Vaikai, suaugusieji, nutukę suaugusieji: kapiliarų kraujo dujų analizė, pilnas kraujo kiekis	Vaikai, suaugusieji, kurių kraujo tūris yra labai didelis, naujagimis: pilnas kraujo skaičius, retikulocitai, medžiagų apykaitos bandymas	Vaikai, suaugusieji su dideliu surinkimo kiekiu, naujagimis: kraujo serumas, cukraus kiekis kraujyje, bilirubinas
<b>Paraiška</b>					

\* Atkreipkite dėmesį, kad tinkamo lanksto parinkimas priklauso nuo įvairių veiksnių, tokių kaip odos tipas. kraujas pratekėjo į punkto vietą ir teisingai pritaikė "Safety Lancet" ant odos. Jei kyla abejonų, prieš pradėdami gaminius naudoti įprastoms programoms, rekomenduojama išbandyti įvairius saugos lankstus. Pateikta informacija yra tik rekomendacija ir jokių būdu nekeičia medicinos, mokslinės ar technologinės konsultacijos.