



Instructions for use

MediSafe Solo - 29G needle
1.5mm penetration dept
MediSafe Solo - 23G needle
2.0mm penetration dept

MediSafe® Solo




1



Twist off and remove
the protective cap.

2



Press the MediSafe®
Solo push-button to
activate the device.

Product simplicity – two step operation to obtain required blood

MediSafe® Solo single use, sterile safety lancets are push-button activated medical devices used to perform finger pricks in order to collect micro samples of capillary blood.



MediSafe® Solo

Naudojimo instrukcijos

MediSafe Solo - 29G adata
1.5mm skvarbumo gylis
MediSafe Solo - 23G adata
2.0mm skvarbumo gylis




1



Atsukite ir nuimkite
apsauginį dangtelį.

2



Paspauskite "
MediSafe Solo"
mygtuką įjunkite
prietaisą.

Produkto paprastumas - dviejų etapų operacija norint gauti reikiamą kraują

MediSafe® Solo vienkartiniai sterilūs saugos lantai yra spyruokliniai suaktyvinti medicinos prietaisai, naudojami pirštui atlikti, siekiant surinkti kapiliarinio kraujo mikro mėginius.