

Interventional Cardiology

Product Catalogue | October 2017

**Boston
Scientific**
Advancing science for life™



Table of Contents

1. Stents & Drug Technologies

Abluminal Bioabsorbable Polymer DES	
• SYNERGY™ **	06
Drug-Eluting Stents	
• Promus PREMIER™	08
• PROMUS Element™ Plus	09
Bare Metal Stent	
• REBEL™	11
Drug Coated Balloon	
• AGENT™	13

2. Complex PCI

Balloon Catheters	
Pre-dilatation Balloon Catheters	
• EMERGE™ Monorail™	16
• EMERGE™ OTW	17
• MAVERICK2™ Monorail™	19
Post-dilatation Balloon Catheters	
• NC QUANTUM APEX™	20
• NC EMERGE™	22
Microdilatation Catheters	
• THREADER™ Microdilatation Catheter	25
Guide Catheters	
• Curves Overview	27
• CONVEY™ 5F / 6F / 7F / 8F	28
• MACH 1™ 6F / 7F / 8F	32
• RUNWAY™ 6F	37
Guide Catheter extension	
• GUIDEZILLA™ II Guide Extension Catheter	41
Angiographic Catheters	
• IMPULSE™ 5F / 6F	43
• EXPO™ 5F / 6F / 7F	46
Guidewires	
Standard Lesions / Spring Tip Guidewires	
• CHOICE™ Floppy	50
• LUGE™ Moderate Support Guidewire	50
Complex Lesions / Polymer Tip Guidewires	
• PT2™ LS	51
• CHOICE™ PT Floppy	51
• PT2™ MS	51
Stiff Guidewires	
• CHOICE™ Intermediate	52
• CHOICE™ Standard	52
• PT GRAPHIX™ Intermediate	52
• PT GRAPHIX™ Standard	53
Extra Support Guidewires	
• CHOICE™ Extra Support	54
• CHOICE™ PT Extra Support	54
• MAILMAN™	54
• PT GRAPHIX™ Super Support	55
SENTAI™ Wires	
• MARVEL™ Workhorse Wire	56
• SAMURAI™ Workhorse Wire	56
• SAMURAI™ RC Specialty Crossing Wire	56
• FIGHTER™ Specialty Crossing Wire	56
• HORNET™ Specialty Crossing Wire	57
• HORNET™ 10 Specialty Crossing Wire	57
• HORNET™ 14 Specialty Crossing Wire	57

Extension Guidewires	
• ADDWIRE™ Extension Wire	58
• The 38 MAGNET™ Device	58
• STRETCH™ Extension wire	58

Embolic Protection

• FILTERWIRE EZ™ Embolic Protection System	60
--	----

Accessories

Inflation Devices

• ENCORE™ 26 Inflation Device	62
• ENCORE™ 26 Advantage Kit	62

Y-Adaptors

• GATEWAY™ Plus Y-Adaptor	63
---------------------------	----

Introducer Sheaths

• SUPER SHEATH™ 0.889 mm (0.035") compatible	64
• SUPER SHEATH™ 0.965 mm (0.038") compatible	64

Cutting Balloon™

• FLEXTO™ Cutting Balloon™	66
----------------------------	----

Atherectomy Devices

• ROTALINK™ Advancer	68
• ROTALINK™ Burr	68
• ROTALINK™ Plus (Advancer & Burr)	68
• ROTAWIRE™	69
• ROTABLATOR™ Consoles	69
• Accessories	70

Chronic Total Occlusions

• CROSSBOSS™ & STINGRAY™ LP	72
-----------------------------	----

TRAPPER™ Exchange Device

• TRAPPER™ Exchange Device	74
----------------------------	----

ANGIOJET™ Thrombectomy System

• Spiroflex™	76
• Spiroflex™ VG	76
• Distaflex™	76

3. PCI Guidance

POLARIS™ PCI Guidance System*

• OPTICROSS™ 6 Coronary Imaging Catheter	79
• OPTICROSS™ Coronary Imaging Catheter	79
• POLARIS™ Accessories	79
• COMET™	79
• POLARIS™ Cart System	80
• POLARIS™ Installed System	80

4. Structural Heart

Transcatheter Aortic Valve Implantation (TAVI)

• ACURATE neo™ Aortic Valve	83
• ACURATE neo™ Transfemoral Delivery System	83
• ACURATE neo™ Transapical Delivery System	83
• ACURATE TA™ Transapical Aortic Bioprosthesis	83
• LOTUS™ Introducer Set	84
• SAFARI2™ Pre-Shaped Guidewire	84
• BOLT™ Direct Access Sheath	84

Left Atrial Appendage Closure (LAAC)

• WATCHMAN™	86
-------------	----

1
Stents
and Drug
Technologies

2
Complex
PCI

3
PCI
Guidance

4
Structural
Heart



2 *Complex PCI*

Guidewires

Standard Lesions / Spring Tip Guidewires

Complex Lesions / Polymer Tip Guidewires

Stiff Guidewires

Extra Support Guidewires

SENTAI™ Wires

Extension Guidewires



2. Standard Lesions / Spring Tip Guidewires

CHOICE™ Floppy

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 1213201 M2 08714729877103	182 cm	Straight	0.8	Spring Coil	Light	Hybrid*	3 cm	The MAGNET™	5
H749 1213201 MJ2 108714729877110	182 cm	J	0.8	Spring Coil	Light	Hybrid*	3 cm	The MAGNET™	5
H749 1213201 2 08714729150541	182 cm	Straight	0.8	Spring Coil	Light	Hybrid*	3 cm	ADDWIRE™	5
H749 1213201 J2 08714729155751	182 cm	J	0.8	Spring Coil	Light	Hybrid*	3 cm	ADDWIRE™	5
H749 1211601 2 08714729150589	300 cm	Straight	0.8	Spring Coil	Light	Hybrid*	3 cm	ADDWIRE™	5
H749 1211601 J2 08714729193654	300 cm	J	0.8	Spring Coil	Light	Hybrid*	3 cm	ADDWIRE™	5



LUGE™ Moderate Support Guidewire

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 1213001 M2 108714729877080	182 cm	Straight	0.9	Floppy; Spring Coil w / Polymer Sleeve	Moderate	Hybrid*	3 cm	The MAGNET™	5
H749 1213001 MJ2 108714729877097	182 cm	J	0.9	Floppy; Spring Coil w / Polymer Sleeve	Moderate	Hybrid*	3 cm	The MAGNET™	5
H749 1213001 2 08714729252320	182 cm	Straight	0.9	Floppy; Spring Coil w / Polymer Sleeve	Moderate	Hybrid*	3 cm	ADDWIRE™	5
H749 1213001 J2 08714729252337	182 cm	J	0.9	Floppy; Spring Coil w / Polymer Sleeve	Moderate	Hybrid*	3 cm	ADDWIRE™	5
H749 1210001 2 08714729252207	300 cm	Straight	0.9	Floppy; Spring Coil w / Polymer Sleeve	Moderate	Hybrid*	3 cm	ADDWIRE™	5
H749 1210001 J2 08714729252214	300 cm	J	0.9	Floppy; Spring Coil w / Polymer Sleeve	Moderate	Hybrid*	3 cm	ADDWIRE™	5



Balloon Catheters

Micro-dilatation Catheters

Guide Catheters

Guide Catheter extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

Trapper™ Exchange Device

ANGIOJET™ Thrombectomy System



Complex PCI - 50
Interventional Cardiology

Boston
Scientific

* Hybrid, Hydrophilic coating Distal 3 cm uncoated.

2. Complex Lesions / Polymer Tip Guidewires

PT²™ LS

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 3893101 2 08714729471110	185 cm	Straight	2.5	Polymer Tip	Light	Hydrophilic	2 cm	ADDWIRE™	5
H749 3893101 J2 08714729471134	185 cm	J	2.5	Polymer Tip	Light	Hydrophilic	2 cm	ADDWIRE™	5
H749 3893102 2 08714729471158	300 cm	Straight	2.5	Polymer Tip	Light	Hydrophilic	2 cm	ADDWIRE™	5
H749 3893102 J2 08714729471172	300 cm	J	2.5	Polymer Tip	Light	Hydrophilic	2 cm	ADDWIRE™	5



CHOICE™ PT Floppy

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 1216001 M2 108714729877264	182 cm	Straight	2.1	Polymer Tip	Light	Hydrophilic	35 cm	The MAGNET™	5
H749 1216001 MJ2 108714729877271	182 cm	J	2.1	Polymer Tip	Light	Hydrophilic	35 cm	The MAGNET™	5
H749 1216001 2 08714729150626	182 cm	Straight	2.1	Polymer Tip	Light	Hydrophilic	35 cm	ADDWIRE™	5
H749 1216001 J2 08714729193715	182 cm	J	2.1	Polymer Tip	Light	Hydrophilic	35 cm	ADDWIRE™	5
H749 1215401 2 08714729150633	300 cm	Straight	2.1	Polymer Tip	Light	Hydrophilic	35 cm	ADDWIRE™	5
H749 1215401 J2 08714729193708	300 cm	J	2.1	Polymer Tip	Light	Hydrophilic	35 cm	ADDWIRE™	5



PT²™ MS

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 3893103 2 08714729471196	185 cm	Straight	2.9	Polymer Tip	Moderate	Hydrophilic	2 cm	ADDWIRE™	5
H749 3893103 J2 08714729471219	185 cm	J	2.9	Polymer Tip	Moderate	Hydrophilic	2 cm	ADDWIRE™	5
H749 3893104 2 08714729471233	300 cm	Straight	2.9	Polymer Tip	Moderate	Hydrophilic	2 cm	ADDWIRE™	5
H749 3893104 J2 08714729471257	300 cm	J	2.9	Polymer Tip	Moderate	Hydrophilic	2 cm	ADDWIRE™	5



Balloon Catheters

Micro-dilatation Catheters

Guide Catheters

Guide Catheter extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

Trapper™ Exchange Device

ANGIOJET™ Thrombectomy System



Complex PCI - 51
Interventional Cardiology

Boston
Scientific

5

2. Stiff Guidewires

CHOICE™ Intermediate

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 1213301 M2 108714729877127	182 cm	Straight	1.1	Spring Coil Tip	Light	Hybrid*	3 cm	The MAGNET™	5
H749 1213301 MJ2 108714729877134	182 cm	J	1.1	Spring Coil Tip	Light	Hybrid*	3 cm	The MAGNET™	5
H749 1213301 08714729150558	182 cm	Straight	1.1	Spring Coil Tip	Light	Hybrid*	3 cm	ADDWIRE™	5
H749 1213301 J2 08714729155744	182 cm	J	1.1	Spring Coil Tip	Light	Hybrid*	3 cm	ADDWIRE™	5
H749 1211701 2 08714729150596	300 cm	Straight	1.1	Spring Coil Tip	Light	Hybrid*	3 cm	ADDWIRE™	5
H749 1211701 J2 08714729193661	300 cm	J	1.1	Spring Coil Tip	Light	Hybrid*	3 cm	ADDWIRE™	5



CHOICE™ Standard

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 1213401 M2 08714729877141	182 cm	Straight	2.8	Spring Coil Tip	Light	Hybrid*	3 cm	The MAGNET™	5
H749 1213401 MJ2 108714729877158	182 cm	J	2.8	Spring Coil Tip	Light	Hybrid*	3 cm	The MAGNET™	5
H749 1213401 08714729150565 2	182 cm	Straight	2.8	Spring Coil Tip	Light	Hybrid*	3 cm	ADDWIRE™	5
H749 1213401 J2 08714729155737	182 cm	J	2.8	Spring Coil Tip	Light	Hybrid*	3 cm	ADDWIRE™	5
H749 1211801 2 08714729150602	300 cm	Straight	2.8	Spring Coil Tip	Light	Hybrid*	3 cm	ADDWIRE™	5
H749 1211801 J2 08714729193678	300 cm	J	2.8	Spring Coil Tip	Light	Hybrid*	3 cm	ADDWIRE™	5



PT GRAPHIX™ Intermediate

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 1490201 M2 08714729877202	182 cm	Straight	1.7	Polymer Tip	Moderate	Hydrophilic	2 cm	The MAGNET™	5
H749 1490201 MJ2 08714729877219	182 cm	J	1.7	Polymer Tip	Moderate	Hydrophilic	2 cm	The MAGNET™	5
H749 1490201 2 ?	182 cm	Straight	1.7	Polymer Tip	Moderate	Hydrophilic	2 cm	ADDWIRE™	5
H749 1490201 J2 ?	182 cm	J	1.7	Polymer Tip	Moderate	Hydrophilic	2 cm	ADDWIRE™	5
H749 1491401 2 ?	300 cm	Straight	1.7	Polymer Tip	Moderate	Hydrophilic	2 cm	ADDWIRE™	5
H749 1491401 J2 ?	300 cm	J	1.7	Polymer Tip	Moderate	Hydrophilic	2 cm	ADDWIRE™	5



Complex PCI - 52
Interventional Cardiology

Boston
Scientific

* Hybrid, Hydrophilic coating Distal 3 cm uncoated.

Balloon
Catheters

Micro-
dilatation
Catheters

Guide
Catheters

Guide
Catheter
extension

Angiographic
Catheters

Guidewires

Embolic
Protection

Accessories

Cutting
Balloon™

Atherectomy
Devices

Chronic Total
Occlusions

Trapper™
Exchange
Device

ANGIOJET™
Thrombectomy
System

2. Stiff Guidewires

PT GRAPHIX™ Standard

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 1490301 M2 08714729877226	182 cm	Straight	5.3	Polymer Tip	Moderate	Hydrophilic	2 cm	The MAGNET™	5
H749 1490301 MJ2 08714729877219	182 cm	J	5.3	Polymer Tip	Moderate	Hydrophilic	2 cm	The MAGNET™	5
H749 1490301 2 08714729188469	182 cm	Straight	5.3	Polymer Tip	Moderate	Hydrophilic	2 cm	ADDWIRE™	5
H749 1490301 J2 08714729141754	182 cm	J	5.3	Polymer Tip	Moderate	Hydrophilic	2 cm	ADDWIRE™	5



Balloon Catheters

Micro-dilatation Catheters

Guide Catheters

Guide Catheter extension

Angiographic Catheters

Guidewires

Embolic Protection

Accessories

Cutting Balloon™

Atherectomy Devices

Chronic Total Occlusions

Trapper™ Exchange Device

ANGIOJET™ Thrombectomy System



Complex PCI - 53
Interventional Cardiology

Boston
Scientific

2. Extra Support Guidewires

CHOICE™ Extra Support

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 1213501 M2 108714729877165	182 cm	Straight	0.9	Spring Coil Tip	Extra	Hybrid*	3 cm	The MAGNET™	5
H749 1213501 MJ2 108714729877172	182 cm	J	0.9	Spring Coil Tip	Extra	Hybrid*	3 cm	The MAGNET™	5
H749 1213501 2 08714729150572	182 cm	Straight	0.9	Spring Coil Tip	Extra	Hybrid*	3 cm	ADDWIRE™	5
H749 1213501 J2 08714729193692	182 cm	J	0.9	Spring Coil Tip	Extra	Hybrid*	3 cm	ADDWIRE™	5
H749 1211901 2 08714729150619	300 cm	Straight	0.9	Spring Coil Tip	Extra	Hybrid*	3 cm	ADDWIRE™	5
H749 1211901 J2 08714729193685	300 cm	J	0.9	Spring Coil Tip	Extra	Hybrid*	3 cm	ADDWIRE™	5



CHOICE™ PT Extra Support

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 1216101 M2 08714729877288	182 cm	Straight	2.0	Polymer Tip	Extra	Hydrophilic	35 cm	The MAGNET™	5
H749 1216101 MJ2 108714729877295	182 cm	J	2.0	Polymer Tip	Extra	Hydrophilic	35 cm	The MAGNET™	5
H749 1216101 2 08714729176992	182 cm	Straight	2.0	Polymer Tip	Extra	Hydrophilic	35 cm	ADDWIRE™	5
H749 1216101 J2 08714729177005	182 cm	J	2.0	Polymer Tip	Extra	Hydrophilic	35 cm	ADDWIRE™	5
H749 1215501 2 08714729177012	300 cm	Straight	2.0	Polymer Tip	Extra	Hydrophilic	35 cm	ADDWIRE™	5
H749 1215501 J 08714729176947 2	300 cm	J	2.0	Polymer Tip	Extra	Hydrophilic	35 cm	ADDWIRE™	5



MAILMAN™

Diameter: .014" (0.37 mm)

Catalog Number / Order Number (GTIN)	Length	Tip Configuration	Tip Load	Tip Style	Rail Support	Coating	Tip Radiopacity	Compatible	Quantity
H749 1213601 M2 108714729877189	182 cm	Straight	0.8	Spring Coil Tip	Super	Hybrid*	3 cm	The MAGNET™	5
H749 1213601 MJ2 108714729877196	182 cm	J	0.8	Spring Coil Tip	Super	Hybrid*	3 cm	The MAGNET™	5
H749 1213601 2 08714729165712	182 cm	Straight	0.8	Spring Coil Tip	Super	Hybrid*	3 cm	ADDWIRE™	5
H749 1213601 J2 08714729165729	182 cm	J	0.8	Spring Coil Tip	Super	Hybrid*	3 cm	ADDWIRE™	5
H749 1212001 2 08714729188162	300 cm	Straight	0.8	Spring Coil Tip	Super	Hybrid*	3 cm	ADDWIRE™	5
H749 1212001 J2 08714729165705	300 cm	J	0.8	Spring Coil Tip	Super	Hybrid*	3 cm	ADDWIRE™	5



* Please check availability of this product with your local sales representative.



Complex PCI - 54
Interventional Cardiology

Boston
Scientific

* Hybrid, Hydrophilic coating Distal 3 cm uncoated.

Balloon
Catheters

Micro-
dilatation
Catheters

Guide
Catheters

Guide
Catheter
extension

Angiographic
Catheters

Guidewires

Emboloc
Protection

Accessories

Cutting
Balloon™

Atherectomy
Devices

Chronic Total
Occlusions

Trapper™
Exchange
Device

ANGIOJET™
Thrombectomy
System