

Urine Monovette®

Efficient sample processing in your urine laboratory



Urine Monovette® for hygienic urine collection, transport and analysis

Three functions combined in one product:

1. Hygienic urine collection

Urine can be collected from a sample container by connecting the suction tip and withdrawing the piston

2. Vessel for test strips

The test strip, e.g. Combur9-test, can be immersed directly into the filled Urine Monovette® for urine analysis

3. Centrifuge tube for sedimentation

The Urine Monovette® can be placed into a centrifuge for sedimentation

Volume ml	Order Number	Tube Length/ Ø in mm	Description	Packaging Pcs./Case
8.5	10.258	92/15	-	500
10	10.252	102/15	-	500
10	10.252.001	102/15	paper label	500
10	10.252.030	102/15	white, UV light protection	500
10	10.252.020	102/15	ind. wrapped, sterile	500
10	10.252.015	102/15	50 pieces sub-packaged	500



Urine Monovette® with stabiliser

for bacteriological urine analysis

- The Urine Monovette® prepared with boric acid effectively stabilises micro organisms from sample collection to analysis
- The bacteriostatic effect of boric acid on the most important causes of urinary tract infections has been proven
- Urine collection and sample transport in one container

Volume ml	Order Number	Tube Length/ Ø in mm	Description	Packaging Pcs./Case
10	10.253	102/15	-	500
10	10.253.001	102/15	paper label	500
10	10.253.020	102/15	ind. wrapped, sterile	500