

**AIDIAN**



## QuikRead go<sup>®</sup> Connectivity

Increase patient safety, operational efficiency and traceability

# Increase patient safety, operational efficiency and traceability

Whether your goal is fast, easy and reliable transfer of results to the patient record or full management of near patient testing, the QuikRead go advanced connectivity features are here to help!

Results coupled with patient and operator identification are instantly sent to the laboratory or hospital information system. Bidirectional control of the QuikRead go instruments through a capable management or middleware system ensure your testing environment is compliant to the strictest of rules.

## Enhanced patient safety

- Send results automatically to reduce manual errors
- Get patient lists to the instrument to ensure correct patient identification
- Add comments to results

## Remote user management


- Manage operator lists
- Send information to QuikRead go operators

## Remote instrument management

- Lock/unlock QuikRead go instruments
- Adjust instrument settings and time
- Set QC limits
- Update instrument software
- Follow quality controls remotely
- See instrument messages

## Versatile connectivity

- RJ-45 network cable or wireless connection with a WLAN USB adapter
- Unidirectional connection by LIS01-A2 or bidirectional connection with POCT1-A2



## Contact info

Aidian Oy

aidian@aidian.eu

+358 10 309 3000

For orders: orders@aidian.eu

www.aidian.eu / www.quikread.com

Product	Cat. no.
QuikRead go Instrument	133893
Printer for QuikRead go Instrument	136203
Barcode reader for QuikRead go Instrument	136204
QuikRead go WLAN adapter	148045

# AIDIAN