



[Downloads \(https://www.tscswabs.co.uk/downloads\)](https://www.tscswabs.co.uk/downloads) |  
[Latest News \(https://www.tscswabs.co.uk/news\)](https://www.tscswabs.co.uk/news) |  
[Contact Us \(https://www.tscswabs.co.uk/contact-us\)](https://www.tscswabs.co.uk/contact-us)

Search

## TS/5-42

---

Pack size: 100

Hygiene Swab Sampling - 10ml Neutralising Buffer - 10ml Neutralising Buffer in tube, polystyrene (breakable) shaft with Polyester (dacron) tip swab in peel pouch.



**Description:** 10ml Neutralising Buffer in tube. Includes a polystyrene (breakable) shaft with Polyester (dacron) tip swab in peel pouch.

**Recommended for:** Hygiene and environmental monitoring of equipment or surfaces testing where traces of disinfectants are possibly present

**Ideal for:** maintaining & transporting microbiological samples & eliminating any traces of sanitisers or disinfectants remaining on surfaces.

Excellent recovery of organisms including: TVC, Listeria spp., Staphylococci, Enterobacteriaceae eg Salmonella, E.coli (E coli / ecoli / e. coli), Enterobacter, staphylococcus, klebsiella, E. coli Enterobacter, Enterobacter sakazakii.

**Accurate volumes for serial dilution techniques**

**Shelf Life:** 2 years.

**Sterilisation Method:**

Gamma Irradiation to achieve a SAL of  $10^{-6}$

**Alternatives:** TS/5-40, TS/6-61, TS/6-63BP

Code	Name	Colour	Length	Quantity	Pack size	No. in case	Shelf life
TS/5-42	Swab Sampling Kits - 10ml NB	n/a	n/a	1	100	n/a	2 years

Download product spec sheet ([uploads/images/PDF-Downloads/Current Product Specifications/spec\\_TS5-42.pdf](#)) Download material safety sheet ([uploads/images/PDF-Downloads/Current MSDS Copies/msds\\_5-42.pdf](#)) Download guidance for use ([uploads/images/PDF-Downloads/Current User Guidance Docs/IFU\\_HygieneSwab\\_V3.pdf](#))

1 Add To My Cart

## Full M40 Compliance

### Proback Transport Swab Range

- Complete Sampling Solution
- Safe Collection & Transport
- Rackable tubes for convenience
- Excellent choice of swabs
- 2 Year Shelf Life

Find out more

([news/6/64/Proback-Full-M40-Compliance](https://www.tscswabs.co.uk/news/6/64/Proback-Full-M40-Compliance))



([news/6/64/Proback-Full-M40-Compliance](https://www.tscswabs.co.uk/news/6/64/Proback-Full-M40-Compliance))

## Online Certification

**You can now receive batch specific quality certificates directly to your e-mail inbox.**

Our system will use the product batch number and your email address to send you the batch specific certificate you need within a few minutes. The name of the person requesting the certificate and the email address it was issued to will be endorsed on the front of the certificate for full traceability.

Find out more (<https://www.tscswabs.co.uk/online-certification>)

Copyright © 2013-16 TSC Technical Service Consultants Ltd | All Rights Reserved | Designed by Ignition CBS  
(<http://www.ignitioncbs.co.uk>)

Technical Service Consultants Ltd, Microbiology House, Fir Street, Heywood, Lancashire, OL10 1NW, UK  
T : **+44 (0)1706 620600** F : **+44 (0)1706 620445** E : **sales@tscswabs.co.uk** (<mailto:sales@tscswabs.co.uk>)  
*Our website uses cookies. By using our website and agreeing to this policy, you consent to our use of cookies.*

***Read our cookie policy here. (<https://www.tscswabs.co.uk/cookie-policy>)***