

Katalog-Nr. Kategorie		Bezeichnung	MERLIN Diagnostika
E1-247-100 H		MICRONAUT-S MH Hannover GN 3	
Layout	K		
1 Test			

	1	2	3	4	5	6	7	8	9	10	11	12
A	AMP 0,5	AMP 1	AMP 2	AMP 4	AMP 8	AMP 16	AMP 32	COL 0,5	COL 1	COL 2	COL 4	COL 8
B	AMS 0,5/4	AMS 1/4	AMS 2/4	AMS 4/4	AMS 8/4	AMS 16/4	T/S 0,25/4,75	T/S 0,5/9,5	T/S 1/19	T/S 2/38	T/S 4/76	T/S 8/152
C	PIT 0,25/4	PIT 0,5/4	PIT 1/4	PIT 2/4	PIT 4/4	PIT 8/4	PIT 16/4	PIT 32/4	FOS 16	FOS 32	FOS 64	FOS 128
D	CTX 0,125	CTX 0,25	CTX 0,5	CTX 1	CTX 2	CTX 4	CTX 8	GEN 0,5	GEN 1	GEN 2	GEN 4	GEN 8
E	CAZ 0,5	CAZ 1	CAZ 2	CAZ 4	CAZ 8	CAZ 16	CAZ 32	TOB 0,5	TOB 1	TOB 2	TOB 4	TOB 8
F	CXM 0,5	CXM 1	CXM 2	CXM 4	CXM 8	CXM 16	CIP 0,0625	CIP 0,125	CIP 0,25	CIP 0,5	CIP 1	CIP 2
G	ERT 0,125	ERT 0,25	ERT 0,5	ERT 1	IMP 1	IMP 2	IMP 4	IMP 8	MOX 0,0625	MOX 0,125	MOX 0,25	MOX 0,5
H	GC	MER 0,125	MER 0,25	MER 0,5	MER 1	MER 2	MER 4	MER 8	MER 16	TGC 0,25	TGC 0,5	TGC 1

[Angabe der Konzentrationen in µg/ml]

[Ausgabe von DIN-Kürzeln]

LEGENDE	
AMP	Ampicillin
AMS	Ampicillin/Sulbactam
CTX	Cefotaxim
CAZ	Ceftazidim
CXM	Cefuroxim
CIP	Ciprofloxacin
COL	Colistin
ERT	Ertapenem
FOS	Fosfomycin
GEN	Gentamicin
IMP	Imipenem
MER	Meropenem
MOX	Moxifloxacin
PIT	Piperacillin/Tazobactam
TGC	Tigecyclin
TOB	Tobramycin
T/S	Trimethoprim/Sulfamethoxazol
GC	Wachstumskontrolle