

| Katalog-Nr. Kategorie |   | Bezeichnung                  | MERLIN<br>Diagnostika |
|-----------------------|---|------------------------------|-----------------------|
| E1-096-100 H          |   | MICRONAUT-S MH Hannover GP 3 |                       |
| Layout                | K |                              |                       |
| 1 Test                |   |                              |                       |

|   |               |                            |                |             |             |             |              |              |             |            |               |            |
|---|---------------|----------------------------|----------------|-------------|-------------|-------------|--------------|--------------|-------------|------------|---------------|------------|
|   | 1             | 2                          | 3              | 4           | 5           | 6           | 7            | 8            | 9           | 10         | 11            | 12         |
| A | PEN<br>0,0625 | PEN<br>0,125               | PEN<br>0,25    | PEN<br>0,5  | PEN<br>1    | PEN<br>2    | PEN<br>4     | MUP<br>1     | MUP<br>2    | MUP<br>4   | MUP<br>8      | MPN<br>256 |
| B | AMP<br>0,5    | AMP<br>1                   | AMP<br>2       | AMP<br>4    | AMP<br>8    | AMP<br>16   | OXA<br>0,25  | OXA<br>0,5   | OXA<br>1    | OXA<br>2   | OXA<br>4      | GNH<br>128 |
| C | AMS<br>2/4    | AMS<br>4/4                 | AMS<br>8/4     | CEZ<br>0,5  | CEZ<br>1    | CFL<br>0,25 | CFL<br>0,5   | CFL<br>1     | CBP<br>0,5  | CBP<br>1   | CBP<br>2      | CBP<br>4   |
| D | MER<br>0,25   | MER<br>0,5                 | MER<br>1       | MER<br>2    | MER<br>4    | MER<br>8    | MER<br>16    | DPT<br>0,5   | DPT<br>1    | DPT<br>2   | RAM<br>0,0625 | RAM<br>0,5 |
| E | GEN<br>0,5    | GEN<br>1                   | GEN<br>2       | GEN<br>4    | TET<br>0,5  | TET<br>1    | TET<br>2     | TET<br>4     | VAN<br>1    | VAN<br>2   | VAN<br>4      | VAN<br>8   |
| F | MOX<br>0,25   | MOX<br>0,5                 | MOX<br>1       | MOX<br>2    | LIZ<br>1    | LIZ<br>2    | LIZ<br>4     | LIZ<br>8     | COX<br>2    | COX<br>4   | COX<br>8      | COX<br>16  |
| G | CIP<br>0,125  | CIP<br>0,25                | CIP<br>0,5     | CIP<br>1    | CIP<br>2    | CIP<br>4    | ERY<br>0,125 | ERY<br>0,25  | ERY<br>0,5  | ERY<br>1   | ERY<br>2      | ERY<br>4   |
| H | GC            | T/S<br>0,03125/<br>0,59375 | T/S<br>0,5/9,5 | T/S<br>1/19 | T/S<br>2/38 | T/S<br>4/76 | T/S<br>8/152 | CLI<br>0,125 | CLI<br>0,25 | CLI<br>0,5 | CLI<br>1      | CLI<br>2   |

[Angabe der Konzentrationen in µg/ml]

[Ausgabe von DIN-Kürzeln]

| LEGENDE |                              |
|---------|------------------------------|
| AMP     | Ampicillin                   |
| AMS     | Ampicillin/Sulbactam         |
| CEZ     | Cefazolin                    |
| COX     | Cefoxitin                    |
| CFL     | Ceftarolin                   |
| CBP     | Ceftobiprol                  |
| CIP     | Ciprofloxacin                |
| CLI     | Clindamycin                  |
| DPT     | Daptomycin                   |
| ERY     | Erythromycin                 |
| GEN     | Gentamicin                   |
| GNH     | Gentamicin high level        |
| LIZ     | Linezolid                    |
| MER     | Meropenem                    |
| MOX     | Moxifloxacin                 |
| MUP     | Mupirocin                    |
| MPN     | Mupirocin nasal              |
| OXA     | Oxacillin                    |
| PEN     | Penicillin G                 |
| RAM     | Rifampicin                   |
| TET     | Tetracyclin                  |
| T/S     | Trimethoprim/Sulfamethoxazol |
| VAN     | Vancomycin                   |
| GC      | Wachstumskontrolle           |