

MICRONAUT MIC-Strip

Susceptibility Testing of Bacteria by the Broth Dilution Method to Determine the Minimum Inhibitory Concentration (MIC)

- Color-coded MIC-Strips
- Easy visual reading
- Micro dilution method
- Packing unit
5 plates with 8 strips each, 11 concentrations per strip (40 tests)



MIC-Strip color-coded

NEW

According to the CLSI-EUCAST recommendations of March 2016

MIC-Strip Colistin

Reference method for susceptibility testing of *Enterobacteriaceae*, *Pseudomonas aeruginosa* and *Acinetobacter* spp. related to Colistin according to EN ISO 20776-1.

1	2	3	4	5	6	7	8	9	10	11	12
GC	COL 0.0625	COL 0.125	COL 0.25	COL 0.5	COL 1	COL 2	COL 4	COL 8	COL 16	COL 32	COL 64

Date of revision: 28/04/2016

