

Karl Storz GmbH & Co. · P. O. Box 230 · D-78503 Tuttlingen

ENDOSCOPES FOR
MEDICINE AND TECHNICAL SCIENCE
INSTRUMENTS FOR
OTO-RHINO-LARYNGOLOGY

Ecolab GmbH & Co. OHG
R&D/Health Care
z.H. Herrn Rudolf Glasmacher
Reisholzer Werftstr. 38-42

40589 Düsseldorf

Your ref.

Your letter of

Our ref.
MGr/WHC

Date
29.10.04


Material compatibility Sekusept® Aktiv

Dear Herr Glasmacher,

we write to confirm that if the product Sekusept® aktiv is used properly (in regard to concentration, duration and water quality), instruments made of high-grade steel will not suffer damage.

Chromium-plated brass components may undergo material changes if the coating is damaged due to wear and tear or mechanical action and the underlying metal is able to come into contact with the solution.

Yours sincerely
Karl Storz GmbH & Co. KG



by order
Dipl.-Ing. (FH) Horst C. Weiss
Quality Engineering