


Peripheral Treatment with ePTFE

4 pirkimo punktas,
ePTFE kraujagyslių protezai



FlowLine Bipore

eptfe kraujagysliniai protezai su spirale

- Bipore design with two different fibril lengths for low thrombogenicity
- An additional ePTFE wrap for enhanced suture retention and high burst strength
- Unique guideline indicates diameter and wall thickness of the graft
- Helical reinforcement for increased resistance against kinking and compression
- Excellent and pliable handling and suture behaviour
- Simple and easy removal of the spiral reinforcement

Ordering Information 6 ir 8mm diametras, ilgis nuo 50 iki 80cm.

Length (cm)	Diameter (mm)	Thin Wall	Thin Wall Reinforced
		Catalogue No.	Catalogue No.
10	5	10TW1005N	10TW1005S
20	6	10TW2006N	-
20	7	10TW2007N	-
20	8	10TW2008N	-
40	5	10TW4005N	10TW4005S
40	6	10TW4006N	-
50	6	10TW5006N	10TW5006S
50	7	10TW5007N	10TW5007S
50	8	10TW5008N	10TW5008S
70	5	10TW7005N	10TW7005S
80	6	10TW8006N	10TW8006S
80	7	10TW8007N	10TW8007S
80	8	10TW8008N	10TW8008S

Length (cm)	Diameter (mm)	Standard Wall	Standard Wall Reinforced
		Catalogue No.	Catalogue No.
20	6	10SW2006N	-
20	7	10SW2007N	-
40	5	10SW4005N	-
50	6	10SW5006N	10SW5006S
50	7	10SW5007N	10SW5007S
50	8	10SW5008N	10SW5008S
70	5	10SW7005N	10SW7005S
80	6	10SW8006N	10SW8006S
80	7	10SW8007N	10SW8007S
80	8	10SW8008N	10SW8008S

ARTIVION™

Learn more at artivion.com

All products and indications are not available/approved in all markets. All trademarks are owned by Artivion, Inc. or its subsidiaries. On-X Life Technologies, Inc. Jotec GmbH and Ascyrus Medical GmbH are wholly owned subsidiaries of Artivion, Inc. © 2022 Artivion, Inc. All rights reserved.

Artivion, Inc.

1655 Roberts Blvd., NW
Kennesaw, GA 30144 USA
Phone: +1-888-427-9654
Fax: +770-590-3753

E-mail: inquiries@artivion.com

For contact information by region, please visit. www.artivion.com/contact

 JOTEC GmbH, Lotzenäcker 23, 72379 Hechingen, Germany

JT-BR-0100200-EN V03 12/2022