

RHYTHMIA HDx™ Disposables

MODEL NUMBER	DESCRIPTION
M004 RC64S 0	INTELLAMAP ORION™, 64 electrodes 8.5 F
	<ul style="list-style-type: none"> • 64 low-noise flat printed electrodes • Digitally Printed • 2.5 mm Inter-electrode spacing center-to-center • 0.4 mm² Electrode Area • Usable Length: 115 cm
	EGM Clarity
	<ul style="list-style-type: none"> • 8.5F, bidirectional steerability • Variable diameter (3–22 mm) for use in various anatomical structures
	Mapping Versatility
	<ul style="list-style-type: none"> • 8 smooth, flexible splines • Flushing port designed to prevent clot formation
	Safety in Design
M004 RAPATCH2 0	Location Reference Patch Kit and Location Cable
	<ul style="list-style-type: none"> • Enables Magnetic and Impedance based tracking of catheters with the RHYTHMIA™ system • Contains 10 ECG Electrodes • 1 Location Reference Back Patch • Sold in individual boxes
M004 RAUMBILICAL2 0	RHYTHMIA™ Umbilical Cable ST Female-to-Female BNC Barrel Adaptors



INTELLAMAP ORION Mapping Catheter



Patch Kit



INTELLANAV STABLEPOINT™ Open-Irrigated Ablation Catheter (incl. Local Impedance and Contact Force)

MODEL NUMBER	CURVE STYLE	SHAFT SIZE	TIP SIZE	EFFECTIVE LENGTH
M004 ERFSDS9620 0	Standard	7.5 F (2.5 mm)	7 F (2.33 mm) / 4 mm (straight)	110 cm
M004 ERFSDS9620K2 0	Large	7.5 F (2.5 mm)	7 F (2.33 mm) / 4 mm (straight)	110 cm


Electrode Configuration: Quadripolar;
Electrode Spacing: 4 – 2.5 – 2.5 mm. Compatible Sheath: ≥8.5 F.



RHYTHMIA HDx™ MAPPING SYSTEM COMPONENTS	
MODEL NUMBER	DESCRIPTION
M004 RARC03 0	INTELLANAV STABLEPOINT Catheter Cable 200 cm (6.6. Ft)
M004 RA6301 0	INTELLANAV STABLEPOINT Connection Box – MAESTRO
M004 RA6302 0	INTELLANAV STABLEPOINT Connection Box – IntellaGen
M004 117 0	METRIQ™ Irrigation Tubing Set
M004 GEN100161001	Open-Irrigated Tubing Kit IntellaGen
M004 GEN100199390	INTELLAGEN™ RF Generator Kit
M004 GEN100199400	INTELLAGEN™ Irrigation Pump Kit

BSC CURVES	CURVES FROM OTHER MANUFACTURERS
Standard	Biosense: D Curve Abbott: D Curve, Medium Curl, Medium Sweep
Large	Biosense: F Curve, J Curve Abbott: F Curve, J Curve, Large Curl, Large Sweep, Large

Only compatible with MAESTRO 4000 Controller & METRIQ Irrigation Pump.




Suggested catheter configurations are illustrative representations only and may not reflect actual performance.



INTELLAMAP ORION™
Mapping Catheter

Boston Scientific
Advancing science for life™

Industry-smallest'
0.4 sq mm
flat electrodes



* Calculation utilizes labeled number of 20 electrodes on two commonly used and commercially available electrophysiology mapping catheters (LASSO® and PENTARAY®).

