

Warszawa, dn. 10.12.2024 r.

To Whom It May Concern,

CONFIRMATION LETTER

We, Mindray Medical Poland Sp. z o.o., with the address at ul. Cybernetyki 9, 02-677 Warsaw, Poland, subsidiary of Shenzhen Mindray Bio-Medical Electronics Co., Ltd., ("Mindray") the manufacturer of Mindray Medical supply unit HyPort P30 hereby confirm the following:

- Meets flammability standard UL 94-V1;
- Console system for DPV monitors and oxygen air, vacuum inlets;

1.1.2. • Vertical column with four panels with passing joints arranged in them (for proper positioning of the oxygen flow meter and the suction cup): gas outlets can be installed on 4 panels;

- On the back side of the console can be installed shelf or other equipment; 1.1.3.
- Arrangement of electrical, low current and medical gas connections: 1. The joints are arranged in the vertical columns of the consoles. 2. The location of a particular connection may be optional (compatible with the consumer);
- Inside the console columns there are separate channels: 1. For electrical inputs; 2. For gas inlets;

1.1.14. • Blocking of console movements: pneumatic brakes for the 2-axis of both consoles.

- Single-phase electrical sockets with grounding (220V, 16A) and voltage indicator, as well as covers – 12 pcs. with color coding of sockets. Two types of sockets are installed: powered from the input and powered from the input and backup power source.;

This confirmation letter is (integral) part of the HyPort P30 datasheet. If there are parameters

with different values in the datasheet and confirmation letter, the information in the confirmation letter should be used.

Signature/stamp

Cao Yuxuan

Regional Sales&Marketing Manager
Patient Monitoring and Life Support
Northeast Europe