


Overview Chart

Arthrex Reference Chart for Recommended Shaver Blades per Indication

3.1.4

Knee	Indication	Description	Article	Diameter	Length	Item Number
	ACL/PCL	Torpedo™		5.0 mm	13 cm	AR-8500TD
		Excalibur shaver B.		5.5 mm	13 cm	AR-8550EX
	Healing response / cartilage repair	Power Pick™, 45°, 6 mm drill depth		1.5 mm	13 cm	AR-8150PX-45
	Meniscal tear repair	Torpedo™, curved		4.0 mm	13 cm	AR-8400CTD
	Meniscal tear repair / healing response	Torpedo™		4.0 mm	13 cm	AR-8400TD
	Notchplasty	Bone cutter		4.0 mm	13 cm	AR-8400BC
		Burr, oval, FlushCut™, 12-flute		5.0 mm	13 cm	AR-8500FOT
	Partial meniscectomy (soft tissue)	Dissector		3.8 mm	13 cm	AR-8380DS
		Excalibur		4.0 mm	13 cm	AR-8400EX
Dissector			5.0 mm	13 cm	AR-8500DS	
Synovectomy	Dissector, curved		4.0 mm	13 cm	AR-8400CDS	
	Excalibur		5.0 mm	13 cm	AR-8500EX	
	Dissector		5.5 mm	13 cm	AR-8550DS	




Standartinio seiverio antgalio ilgis yra 13cm, Smulkiuju sanariu 7 cm, o klubo seiverio antgaliu 19 cm

www.Arthrex.com
suvedus koda -
matomi visu antgaliu
ilgiai

3.1.4

Overview Chart

Arthrex Reference Chart for Recommended Shaver Blades per Indication




	Indication	Description	Article	Diameter	Length	Item Number
Shoulder	Acromioplasty	PowerRasp™		5.5mm	13 cm	AR-8550PR
	Debridement	Excalibur, curved		4.0mm	13 cm	AR-8400CEX
		Dissector		4.0mm	13 cm	AR-8400DS
		Excalibur		4.0mm	13 cm	AR-8400EX
		Dissector		5.5mm	13 cm	AR-8550DS
		Excalibur		5.5mm	13 cm	AR-8550EX
		Debridement / rotator cuff repair	Torpedo™		5.0mm	13 cm
	Labral tear	Torpedo™, curved		4.0mm	13 cm	AR-8400CTD
		Torpedo™		4.0mm	13 cm	AR-8400TD
	Rotator cuff repair	Bone cutter		5.0mm	13 cm	AR-8500BC
		Excalibur		5.0mm	13 cm	AR-8500EX
	Rotator cuff repair / acromioplasty	Burr, oval, FlushCut™, 12-flute		5.0mm	13 cm	AR-8500FOT 3.2.4
Small joint for use with AR-8330SJ, small joint Synergy resection handpiece	Ankle SJ /wrist	Dissector, small hub		3.0mm	7 cm	AR-9300DS for use with AR-8330SJ
		Oval burr, 10-flute, small hub		3.0mm	7 cm	AR-9300OBT for use with AR-8330SJ
	Wrist	Control sabre, small hub		2.0mm	7 cm	AR-9200SR for use with AR-8330SJ



Overview Chart

Arthrex Reference Chart for Recommended Shaver Blades per Indication

3.1.4-3.2.4

	Indication	Description	Article	Diameter	Length	Item Number
Hip	Impingement	Round burr, HL, 8-flute		4.0mm	19 cm	AR-6400RBE
		Round burr, concave, HL, 8-flute		5.0mm	19 cm	AR-6500VBE
	Labral tear	Bone cutter, curved, HL		4.2mm	19 cm	AR-6420CBC
		Excalibur, curved, HL		4.2mm	19 cm	AR-6420CEX
		Dissector, HL		4.2mm	19 cm	AR-6420DS
	Labral tear / impingement / synovectomy	Torpedo™, HL		4.2mm	19 cm	AR-6420TD
	Synovectomy	Dissector, curved, HL		4.2mm	19 cm	AR-6420CDS
Elbow	Cubital tunnel repair	Dissector, small hub		3.0mm	7 cm	AR-9300DS for use with AR-8330SJ
	Loose bone extraction	Bone cutter, curved		4.0mm	13 cm	AR-8400CBC
	Osteophyte removal	Oval burr, 8-flute		4.0mm	13 cm	AR-8400OBE
	Synovectomy	Dissector, curved		3.5mm	13 cm	AR-8350CDS
		Torpedo™, curved		4.0mm	13 cm	AR-8400CTD