

Accessories

Wound Retractors



Code	Description	Unit Of Measure
LMQGS	Small Wound Retractor (2.5 - 6cm)	5 Units/Box
LMQGM	Medium Wound Retractor (5 - 9cm)	5 Units/Box
LMQGL	Large Wound Retractor (9 - 14cm)	5 Units/Box
LMQXL	Extra Large Wound Retractor (11 - 17cm)	5 Units/Box

3.pd.

Smoke Evacuation Filter



Code	Description	Unit Of Measure
LM444	Disposable Smoke Evacuation Filter	10 Units/Box



L o c a M e d

PRODUCT INFORMATION DATA SHEET

The following information should be read before using this device.

FOLDED WOUND RETRACTOR

CODE LMQGS (SMALL) LMQGM (MEDIUM) LMQGL (LARGE) LMQGXL (X LARGE)



Description:

The LocaMed Folded Wound Retractor is a sterile and single-use product. It consists of a durable retracting sleeve with adjustable folding ring. The wound retractor is available in three sizes, Small (2.5~6cm incision), Medium (5~9cm incision), Large (9~14cm incision) and X-Large 11~17.5cm incision

Intended use:

The LocaMed Folded Wound Retractor is used as a wound retraction device which can quickly provide the atraumatic retraction, significantly decreases the infection on the site and enables the direct visualization of wound during the open or laparoscopic surgery.

Preparation:

- Select an appropriate size of wound retractor and create a compatible trans-umbilical incision as required according to hospital procedures.
- Remove Folded Wound Retractor from the package.
- Squeeze the blue ring through the incision into the abdomen cavity and adjust the wound



LocaMed Limited. Unit 4 The Enterprise Centre, Alton Road, Farnham GU10 5EH
Tel: +448701 166300 - Fax +448701 164466 - www.locamed.com





L o c a M e d

retractor by making the blue ring evenly attached to the inner abdomen.

- Up-pull the retracting sleeve until the bottom ring tightly touches the inner side of the abdomen, then roll down the upper white ring until it tightly touches the outer abdomen.
- Perform the surgical procedure through the retracted site.
- Remove the in-situ device from the incision by taking out the blue ring from the cavity.

Cautions / Warnings:

- This device is EO sterilised and for single use only. Please do not re-use or re-process to avoid irreversible injury, infection and device malfunction.
- Use of the LocaMed Folded Wound Retractor should be performed only by physicians with adequate training and knowledge of these procedures. In addition, medical literature should be consulted regarding techniques hazards, contraindications, and complications prior to the performance of these procedures.
- The site needs to be inspected for hemostasis after the assembly is removed from cavity.
- Do not use when the package is opened or damaged.
- Store this device in a dry, clean and safe place. Keep the device away from flammable and combustible chemicals or materials.
- Please follow universal guidelines upon medical waste disposal to avoid any potential cross-infection, environmental contamination or injury caused by sharps waste.

This Product Information Data Sheet is subject to change without notice.

Version: 2

Revision Date: 2015/03/30



LocaMed Limited. Unit 4 The Enterprise Centre, Alton Road, Farnham GU10 5EH
Tel: +448701 166300 - Fax +448701 164466 - www.locamed.com



To whom it may concern,

This is to confirm the length of the LMQGM Wound retractor is 20cm

Charles Ford
International Sales Manager





LocaMed
surgery in safe hands

To Whom it may concern,

I can confirm that the LocaMed Polymer Clips comply with the following requirements and specifications:

- Non-absorbable polymer is inert, non-toxic;
- Nonconductive, radiolucent and does not interfere with CT, MRI or X-ray diagnostics;
- sterile.

Charles Ford
International Sales Manager




6 p.d. - Transabdominalinēs adatos genētīnēms amniocentezēms 22Gx15 cm

PAJUNK®



ChibaSono & Chiba

*Fine Needle
Aspiration Biopsy*



ChibaSono & Chiba

Echogenic & Non-echogenic Chiba Needles

Ultrasound guided biopsies have been increasing in the recent years. Together with Dr. Chris Mitchell, PAJUNK® has developed the ChibaSono, a new generation of needles that allow for a precise localisation under ultrasound. Needle and tip design allow the user to focus on the procedure and patient by making the tip visualisation under ultrasound faster, easier and safer. The wide range of ChibaSono and Chiba needles, sizes and lengths offers the optimal tool for percutaneous cholangiography, cytological biopsy and general percutaneous applications. **Needle material – medical stainless steel.**

Nerūdijančio plieno

Clear Plastic Hub

→ Early detection of liquid flow

Su mandrenu

Depth Stop and
Depth Markings

→ Adjustable for easy
determination of puncture depth


→ Markings in 1 cm steps
for easy reading and
identification

Detachable
Plastic Wings

→ Optimum handling for
different gripping techniques

Retaining Clip (optional)

→ Fixation of stylet and needle





Adatos yra echogeniškos, t. y. turi specialią šiurkštumą prie adatos galiuko, kuri yra gerai matoma ultragarsu;

Cornerstone Reflectors

360 degree graduations on the first 20mm of the needle

- Optimised ultrasound visibility of needle shaft¹
- Reliable and optimised needle visualisation at any angle^{2,3}



Echogenic Needle Tip

US tip with two inclination angles

- Improves needle tip visibility under ultrasound



AVAILABLE OPTIONS

Tip options for Chiba:

Beveled Tip

- ▶ Designed for fluoroscopy guided procedures
- ▶ Precision-ground metal stylet
- ▶ Minimized tissue carry-out into the puncture
- ▶ Increased needle stability



US Tip

- ▶ Special grinding with two inclination angles improves needle tip visibility under ultrasound
- ▶ Rounded, blunt stylet – set back from tip



ChibaSono only with US tip

6 p.d. - Transabdominalinēs adatos genetinēms amniocentezēms 22Gx15 cm

ChibaSono & Chiba

ChibaSono

with US tip, Cornerstone Reflectors, depth markings and wings

Size	Item no.	PU
22G x 50 mm (2")	530S070050	25
22G x 90 mm (3 1/2")	530S070090	25
22G x 100 mm (4")	530S070100	25
22G x 125 mm (5")	530S070125	10
22G x 150 mm (6")	530S070150	10
22G x 220 mm (8 2/3")	530S070220	10
20G x 90 mm (3 1/2")	530S095090	25
20G x 120 mm (4 3/4")	530S095120	10
20G x 150 mm (6")	530S095150	10
20G x 220 mm (8 2/3")	530S095220	10
18G x 150 mm (6")	530S120150	10

Chiba

with beveled tip, retaining clip and wings

Size	Item no.	PU
22G x 90 mm (3 1/2")	201S070090	25
22G x 120 mm (4 3/4")	201S070120	10
22G x 150 mm (6")	201S070150	10

Chiba

with beveled tip, depth stop, depth markings, retaining clip and wings

Size	Item no.	PU
22G x 120 mm (4 3/4")	202S070120	10
22G x 150 mm (6")	202S070150	10
22G x 220 mm (8 2/3")	202S070220	10
20G x 120 mm (4 3/4")	202S095120	10
20G x 150 mm (6")	202S095150	10
20G x 220 mm (8 2/3")	202S095220	10
20G x 280 mm (11")	202S095280	10

6 p.d.

Chiba

with US tip, retaining clip and wings

Size	Item no.	PU
22G x 220 mm (8 2/3")	200S070220	10



Scan now!

PAJUNK® GmbH Medizintechnologie
Global Headquarters
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 9291-0
Fax: +49 7704 9291-600
info@pajunk.com
pajunk.com

PAJUNK® Medical Produkte GmbH
Central Europe
Karl-Hall-Straße 1
78187 Geisingen, Germany
Tel.: +49 7704 8008-0
Fax: +49 7704 8008-150
info@pajunk-medipro.com
pajunk.de & pajunk.eu

PAJUNK® Swiss AG
Switzerland
Neuhofstrasse 5A
6340 Baar, Switzerland
Tel.: +41 417 609 627
info@pajunk.ch
pajunk.ch

PAJUNK® Medical System L.P.
United States of America
4575 Marconi Drive
Alpharetta, GA 30005, USA
Tel.: +1 888 972 5865
Fax: +1 678 514 3388
info@pajunk-usa.com
pajunkusa.com

PAJUNK® UK Medical Products Ltd
United Kingdom
Unit D1 The Waterfront
Goldcrest Way, Newburn Riverside
Newcastle upon Tyne, NE15 8NY, UK
Tel.: +44 191 264 7333
info@pajunk.co.uk
pajunk.co.uk

PAJUNK®
Real People. Real Care.



Ecological Action



Social Commitment



Responsible Management



Cardboard min. 80% Recyclable Material

pajunk.com/sustainability

pajunk.com

HS130000 a

PAJUNK®



PAJUNK® GmbH
Medizintechnologie
Karl-Hall-Strasse 1
78187 Geisingen
Germany

Made in Germany

PAJUNK®

REF 530S070150 ChibaSono 22G x 150mm



(01)04048223128480(11)251224(17)301223(10)1615

2025-12-24
2030-12-23

LOT 1615

STERILE EO

MD QTY 10



PAJUNK®



PAJUNK® GmbH
Medizintechnologie
Karl-Hall-Strasse 1
78187 Geisingen
Germany

Made in Germany

