

## Superior nystagmus observation

VF405 provides ideal conditions for fixation-free observation and recording of eye movements during

1.1 Punktas. Paskirtis

1.2 Punktas. Aparatūrinė dalis

### The Combi Goggle

The VF405 system incorporates the industry's leading goggle design, also known from the clinical VN415/VO425 FireWire® VNG systems.

- Light weight with easy fitting magnetic cover
- Built in fixation light
- Disposable foam cushions provide a sanitary solution
- Comfortable
- Completely light tight fit
- One or two cameras – upgrade as you like

The FireWire® cameras connect directly to a laptop or a desktop computer. Infrared illumination and high resolution cameras with focus adjustment allow for superior image quality.

- Clear detection of torsional movement
- IEEE 1394a FireWire® standard





**Video Recording**

The VF405 displays eye images in partial or full screen mode; simply double-click on the image to toggle between image sizes.  
 Recording indicator with elapsed time.  
 Pre-labelled tests or create your own test names  
 Select High resolution recording for extremely detailed images  
 Remote controlled video recording and fixation light with footswitch  
 Video files are stored in AVI file format

**Video Playback**

Advanced Video Player with playback at variable speed down to one picture at a time  
 Highlight videos of interest  
 Notes can be created to indicate events of interest within the video

**CUSTOMER SERVICE**

**Phone:** +45 6371 3555  
**Fax:** +45 6371 3522  
**E-mail:** [info@interacoustics.com](mailto:info@interacoustics.com)

**REMOTE SUPPORT**

Do you need a quick solution to a technical problem?  
[Download support client](#)

**CONNECT WITH US!**

**YouTube:** [Go to our channel](#)  
**LinkedIn:** [Connect with us!](#)  
**Facebook:** [Follow us](#)

*\* Please note that laptops and computer monitors are only used to showcase our software. These are not part of the Interacoustics brand. The binocular version comes with two infrared cameras and the monocular comes with one infrared camera. FireWire and the FireWire symbol are trademarks of Apple Computer, Inc., registered in the U.S. and other countries. The FireWire logo is a trademark of Apple Computer, Inc.*

Goto: [Frontpage](#) ▶ [VF405](#)

**Eye Movement Recording**

[VN415 / VO425](#)  
[EN415 / EO425](#)  
[VF405](#)  
[EyeSeeCam](#)

**Caloric Irrigators**

[Air Fx](#)  
[Aqua Stim](#)

**Rotary Chair**

[Nydiag 200](#)

**Evoked Potentials**

[VEMP module](#)

**Visuel Stimuli**

[Digital Lightbar](#)

**Interacoustics integration software**

[Interacoustics Suites](#)  
[OtoAccess™ Database](#)  
[Electronic Medical Record](#)

Our Interacoustics Software Suites combines with Electronic Medical Records (EMR), OtoAccess™ Interacoustics database, NOAH or XML and even features combined reporting.

**Ear Tips and Disposables**

[Download General Catalogue](#)  
[Download Disposables Brochure](#)  
[Download Ear Tip Flyer](#)

**INTERACOUSTICS SUPPORT**  
 Are you a professional?  
 Go to our Extranet website for downloads and more

**EVENTS**  
 Find local exhibitions and events

Join our youtube channel for exciting videos and tutorials

**1.4 Punktas. Kamera. Sujungimas su kompiuteriu tiesioginis per IEEE1394 FireWire; rezoliucija 640 x 240, 105 Hz; 640 x 480, 60 Hz.**

**Technical Specifications**

**Camera**

IEEE1394 FireWire®

**Resolution**

640 x 480 Pixels @ 60 Hz  
640 x 240 Pixels @ 105 Hz

**Images**

60 or 105 images per second

**Safety standards**

IEC 60601-1 (EMC)  
IEC 60601-1-2 (EMC)

**Dimensions (L x W x H)**

302 x 216 x 131mm

**Physical Specifications**

Goggle Weight with one camera

240g (non-occluded view)  
320g (occluded view)

with two cameras

305g (non-occluded view)  
385g (occluded view)

Dispensing box with 24 pcs. of disposable goggle foam pads

**Hardware requirements**

Intel Core 2 Duo processor 1.8 GHz or better  
Minimum 2 GB of DDR3 RAM  
160GB harddisk available  
Two USB ports available  
Monitor resolution: 1024 x 768  
Laptop PC: One 34 mm PCIeExpressCard slot available  
Desktop PC: One free slot for IEEE1394a FireWire® PCI board

**Software requirements**

Windows® 7 (32-bit or 64-bit)  
Windows® Vista Business SP  
Windows® Vista Ultimate operating system  
Windows® XP SP2  
Microsoft .Net 2.0 Framework  
Windows Installer 3.x

**1.3 Punktas. Akiniai. Keičiamos akinių pagalvėlės dėžutė su 24 vnt.; Akinių su dviem video kamerom svoris 240 g, kai akiniai atidengti, ir 320 g, kai akiniai uždengti dangteliu/**

**1.6 Punktas. Reikalavimai kompiuterinei ir programinei įrangai. procesorius Intel Core 2 Duo 1,8 GHz arba geresnis, Minimaliai 2 GB DDR3 RAM, 2 USB portai, Monitoriaus rezoliucija 1024x768, nešiojamas kompiuteris su vienu 34 mm PCIeExpressCard, Stacionarus kompiuteris turi IEEE 1394a FireWire, Operacinė sistema Windows 7 (32 bitų arba 64 bitų) su papildymais, leidžiantis dirbti su įrangos programinę įranga ir klinikoje naudojamomis programomis.**

**1.10 Punktas. Standartai. Atitinka standartus IEC 60601-1, IEC 60601-1-1; IEC 60601-1-2. Turi CE ženklą.**

	VF405m Basic	VF405b Basic	VF405m Extended	VF405b Extended
	✓	✓	✓	✓
	✗	✗	✓	✓
Video storage	✗	✗	✓	✓
Video play back	✗	✗	✓	✓
Elapsed test time	✗	✗	✓	✓
Report editor	✓	✓	✓	✓
Sound recording	✗	✗	✓	✓
Test design template	✗	✗	✓	✓
Footswitch operated video capture	✗	✗	✓	✓
Footswitch operated fixation light	✗	✗	✓	✓
Toggle one camera/two cameras	✗	✓	✗	✓
SW controlled Fixation light	✓	✓	✓	✓
External camera enabled	✗	✗	✓	✓
Toggle picture quality high/normal	✓	✓	✓	✓

\* Licensed features with each configuration

**Included parts**

	VF405m Basic	VF405b Basic	VF405m Extended	VF405b Extended
Camera FireWire® complete	1	2	1	2
Combi Mask VOG	1	1	1	1
FireWire® SDM express card 2P	1	1	1	1
Power supply 10W	1	1	1	1
FireWire® PCI-Card IEEE1394		1	1	1

EX-10000 Express card kit	1	1	1	1
Cable fasteners	3	3	3	3
Disposable goggle foam pads	24	24	24	24
Screwdriver SE 1.5 x 60mm 105	1	1	1	1
Web camera w. microphone	0	0	1	1
Installation sheet	1	1	1	1
Combined CE & Operation Manual	1	1	1	1
Installation CD	1	1	1	1

1.8 Punktas. Apmokymas ir instaliavimas. -instaliavimo kompaktinis diskas; - akiniai su viena video kamera ir laidais; - keičiamos akinių paraloninės pagalvėlės - 1 dėžutė; - FireWire PCI Express adapteris plokštė; - pacientų archyvavimo duomenų bazės programa; - kojinis pedalas vaizdo įrašymui paleisti; - instrukcijos anglų ir lietuvių kalbomis.

CUS

Pho

Fax: +45 6371 3522

E-mail: [info@interacoustics.com](mailto:info@interacoustics.com)

to a technical problem?  
[Download support client](#)

[LinkedIn: Connect with us!](#)  
[Facebook: Follow us](#)

\* Please note that laptops and computer monitors are only used to showcase our software. These are not part of the Interacoustics brand.  
 The binocular version comes with two infrared cameras and the monocular comes with one infrared camera.  
 FireWire and the FireWire symbol are trademarks of Apple Computer, Inc., registered in the U.S. and other countries. The FireWire logo is a trademark of Apple Computer, Inc.

Goto: [Frontpage](#) ▶ [VF405](#) ▶ [Specifications](#)

## Eye Movement Recording

[VN415 / VO425](#)  
[EN415 / EO425](#)  
[VF405](#)  
[EyeSeeCam](#)

## Evoked Potentials

[VEMP module](#)

## Caloric Irrigators

[Air Fx](#)  
[Aqua Stim](#)

## Rotary Chair

[Nydiag 200](#)

## Visuel Stimuli

[Digital Lightbar](#)

## Interacoustics integration software

[Interacoustics Suites](#)  
[OtoAccess™ Database](#)  
[Electronic Medical Record](#)

Our Interacoustics Software Suites combines with Electronic Medical Records (EMR), OtoAccess™ Interacoustics database, NOAH or XML and even features combined reporting.

## Ear Tips and Disposables

[Download General Catalogue](#)  
[Download Disposables Brochure](#)  
[Download Ear Tip Flyer](#)

## INTERACOUSTICS SUPPORT

Are you a professional?  
 Go to our Extranet website for downloads and more

## EVENTS

Find local exhibitions and events

Join our youtube channel for exciting videos and tutorials

---

# VF405 Video Frenzel

For superior  
observation and  
recording of  
eye movements



# Superior nystagmus observation

VF405 provides ideal conditions for fixation-free observation and recording of eye movements during

- Spontaneous testing
- Positional and Dix-Hallpike testing
- Head Impulse testing
- User defined testing

## The Combi Goggle

The VF405 system incorporates the industry's leading goggle design, also known from the clinical VN415/VO425 FireWire® VNG systems.

- Light weight with easy fitting magnetic cover
- Built in fixation light
- Disposable foam cushions provide a sanitary solution
- Comfortable
- Completely light tight fit
- One or two cameras - upgrade as you like


The FireWire® cameras connect directly to a laptop or a desktop computer. Infrared illumination and high resolution cameras with focus adjustment allow for superior image quality.

- Clear detection of torsional movement
- IEEE 1394a FireWire® standard




Light tight goggle with sanitary, disposable face cushions

1.5 Punktas. Testai. Įrašytų tyrimo vaizdų peržiūra, atsakymo pacientui redagavimas, garso įrašymas, testų dizaino sukūrimas.



Playback mode



Automatic or manual camera brightness adjustment



Synchronized videos of eye movements and patient movements

## Video Recording

The VF405 displays eye images in partial or full screen mode; simply double-click on the image to toggle between image sizes.

- Recording indicator with elapsed time.
- Pre-labelled tests or create your own test names
- Select High resolution recording for extremely detailed images
- Remote controlled video recording and fixation light with footswitch
- Video files are stored in AVI file format

## Video Playback

- Advanced Video Player with playback at variable speed down to one picture at a time
- Highlight videos of interest
- Notes can be created to indicate events of interest within the video



Playback controls

Synchronized replay  
of situation video, eye  
image(s), and sound

Integrated patient  
database

### External Situation Video


The VF405 includes an external camera for recording of sound and video images from the examination.

- Synchronized sound, eye image video and external situation video.

### The Software

The VF405 is an independent software module operated from the OtoAccess™ patient database. The user interface is designed for maximum user satisfaction with intuitive procedures and a supportive design.

- View all sessions for every patient in the database
- Multi-language interface
- Report editor with word processing function



Available tests panel



1.5 Testai.  
Spintaninis  
nystagmus testas;  
Pozicinis  
testas;Dix-Hallpike  
testas;Vartotojo  
apibrėžtas testas;



## Interacoustics a/s

Interacoustics is a world leading diagnostic solutions provider in the field of hearing and balance assessment. We help the professional audiology world reach new milestones through continuous developments and a constant focus on integration and direct customer value.

## With you at all times

Whenever you purchase audiological equipment with the Interacoustics brand on it, you are not just guaranteed a quality product but also highly trained support personnel with specialized knowledge. We operate in over 100 countries worldwide through a carefully selected network of distributors and service centres to ensure that you receive total support and service.

## Product specifications

All technical and hardware specifications concerning the VF405 can be downloaded from our website.

Read more at  
[interacoustics.com](http://interacoustics.com)

Go online to  
explore our  
full product  
range

## Interacoustics a/s

Drejervaenget 8  
DK-5610 Assens

T +45 6371 3555  
F +45 6371 3522

[info@interacoustics.com](mailto:info@interacoustics.com)  
[interacoustics.com](http://interacoustics.com)

FireWire and the FireWire symbol are trademarks of Apple Computer, Inc., registered in the U.S. and other countries. The FireWire logo is a trademark of Apple Computer, Inc.



# Interacoustics