



TECHNICAL DATA SHEET

Product Name: Container Label Steam

Product Description:

Our 'MEDSTER SafeCheck Container Label Steam' is designed to ensure that sterilized materials are accurately recorded. There are four lines on the product to record sterilization date, sterilization number, operator name information. When exposed to steam sterilization conditions, the indicator ink will change from the initial color to the reference color. LOT and expiry date information can be found on the product label. The shelf life of our product is 3 years.

Complying Harmonized Standards

EN ISO 13485: 2016

EN ISO 14971: 2019

ISO TR 24971:2020

ISO 20417:2021

ISO 15223-1:2021

ISO 11140-1:2014





Other Regulations, Directive and Guidelines

Medical Devices Regulation (EU) 2017/745

Product Properties:

Product Size:	3,5cm x 7,5cm
Paper Type:	American Bristol
Sticky Feature:	N/A
Cellphone:	Not Included
Environmental impact:	Lead and metal free
Sterilization Type:	Steam

Indicator Properties:

Product Code	Product Description	Printed Indicator Type	Indicator Initial Color	Indicator Conversion Color
80.100	Container Label Steam	Steam Indicator	 	 

Chemical Components:

Chemical Ingredients Name	% of ingredients	CAS No
Water	%35	7732-18-5
Diacetone Alcohol	%3	123-42-2
Proprietary Combination	%62	N/A

Packaging Information:

Product Code	Product Description	Package Contents	Inner Box Quantity
80.100	Container Label Steam	100 pcs	50 packs

Storage, Storing, and Disposal Information:

- ☐ Do not expose to direct sunlight and keep it dry.
- ☐ Between 0°C and 40°C and Max. Store at 70% RH humidity.
- ☐ Store in its original packaging.
- ☐ Dispose of in accordance with the place of use and local regulations.



Medster Tıbbi Cihaz ve Sağlık Hizmetleri LTD. ŞTİ.

Başkent OSB Malıköy 19. Cad No.51 Sincan Ankara TURKEY

Tel: +90 312 387 05 50 Fax: +90 312 387 05 51 www.medster.com.tr

TDS.80.100/REV.01/02.01.2023