



Q Path® Coverquick mounting agents for microscopy

A range of mounting media for manual and automatic mounting of slides. The Q Path® Coverquick range consists of four different mounting media with different viscosities to meet all the need in microscopy. The number after the name indicates the viscosity: Q Path® Coverquick 1000 has the lower viscosity (around 260 mPaS) and Q Path® Coverquick 4000 has the higher viscosity (around 920 mPaS).

- A complete range of mounting agents for histology and cytology
- The entire range is suitable for manual mounting
- For automatic mounting the Q Path® Coverquick 1000 and 2000 are recommended
- Delivered in pet amber bottle, only Q Path® Coverquick 3000 is delivered in a aluminium bottle
- Dries quickly in 15 to 20 minutes

Description	Pk	Cat. No.
Q Path® Coverquick 1000, 260 mPaS	500 ml	05547528.
Q Path® Coverquick 2000, 470 mPaS	500 ml	05547530.
Q Path® Coverquick 2000, 470 mPaS	1 l	05547531.
Q Path® Coverquick 3000, 610 mPaS	500 ml	05547537.
Q Path® Coverquick 4000, 920 mPaS	500 ml	05547539.
Q Path® Coverquick 4000, 920 mPaS	1 l	05547540.




Histology mounting media



Description	Pk	Cat. No.
Q Path® Coverquick 1000, 260 mPaS	500 ml	05547528.
Q Path® Coverquick 2000, 470 mPaS	500 ml	05547530.
Q Path® Coverquick 2000, 470 mPaS	1 l	05547531.
Q Path® Coverquick 3000, 610 mPaS	500 ml	05547537.
Q Path® Coverquick 4000, 920 mPaS	500 ml	05547539.
Q Path® Coverquick 4000, 920 mPaS	1 l	05547540.
Q Path® Isomount 2000, 470 mPaS	500 ml	05547535.
Q Path® Isomount 20000, 470 mPaS	1 l	05547536.
Safemount mounting media Q Path®	500 ml	00647520.




Mounting media, HistoChoice®



HistoChoice® Mounting Media was designed for use with HistoChoice® Clearing Agent, but can also be used as a separate, stand-alone medium for standard slide preparation. It is suitable for both thick and thin specimen mounting. It can also be used with immersion oil.

- Low viscosity means easier application resulting in bubble-free preparations which saves valuable laboratory time

Description	Pk	Cat. No.
Mounting media, HistoChoice® 120 ml	120 ml	H157-120ML
Mounting media, HistoChoice® 475 ml	475 ml	H157-475ML



E-BUSINESS

More than just a webshop for all of your laboratory needs

Please visit us at: <http://eu.vwr.com/e-business>