

## Transferpette® S, adjustable volume

Volume range µl	Volume step µl	A* ≤ ± %	CV* ≤ %	Increment µl	Recommended type of tip, µl
0.1 - 1	1	2	1.2	0.001	0.1 - 20
	0.5	4	2.4		
	0.1	20	12		
0.1 - 2.5	2.5	1.4	0.7	0.002	0.5 - 20
	1.25	2.5	1.5		
	0.25	12	6		
0.5 - 10	10	1	0.5	0.01	0.5 - 20
	5	1.6	1		
	1	7	4		
2 - 20	20	0.8	0.4	0.02	2 - 200
	10	1.2	0.7		
	2	5	2		
5 - 50	50	0.8	0.3	0.05	2 - 200
	25	1.2	0.5		
	5	4	2		
10 - 100	100	0.6	0.2	0.1	2 - 200
	50	0.8	0.4		
	10	3	1		
20 - 200	200	0.6	0.2	0.2	2 - 200
	100	0.8	0.3		
	20	3	0.6		
100 - 1000	1000	0.6	0.2	1	50 - 1000
	500	0.8	0.3		
	100	3	0.6		
500 - 5000	5000	0.6	0.2	5	500 - 5000
	2500	0.8	0.3		
	500	3	0.6		
1000 - 10000	10000	0.6	0.2	10	1000 - 10000
	5000	0.8	0.3		
	1000	3	0.6		

\* A = Accuracy, CV = Coefficient of Variation



Final test values related to the nominal capacity (maximum volume) indicated on the instrument, and the indicated volume steps obtained when instrument and distilled water are equilibrated at ambient temperature (20 °C/68 °F) and with smooth operation. According to DIN EN ISO 8655.