

EF-0912 – CE-Declaration of Conformity

Konformitätserklärung/ Déclaration de Conformité
Version 5

Adolis

Trusted Solutions,
Passionate People



2.

We declare under our sole responsibility that the product listed below

Wir erklären in alleiniger Verantwortung, dass das unten aufgeführte Produkt
Nous déclarons sous notre propre responsabilité que le dispositif médical

Product description, Code, Type, Model

Produktbezeichnung; Code, Typ, Modell
Nom, code, type, modèle

Action 3 NG

manufactured by:

hergestellt durch:
fabriqué par:

Invacare France Operations

Route de Saint Roch

37230 FONDETTES

FRANCE

meets all the provisions of the of the Medical Device Directive 93/42/EEC, Annex VII, which apply to it.

allen Anforderungen der Medizinprodukte-Richtlinie 93/42/EWG, Anhang VII, entspricht, die anwendbar sind.
remplit toutes les exigences de la directive sur les dispositifs médicaux 93/42/CEE, annexe VII, qui le concernent.

Device Classification according to Annex IX

Class I

Produktklassifizierung nach Anhang IX
Classification selon l'annexe IX

Klasse I

Classe I

Applied harmonised standards,
national standards or other
normative documents

EN12183 :2014

Angewandte harmonisierte Normen, natio- nale
Normen oder andere normative Doku- mente

Normes harmonisées, normes nationales et
autres documents normatifs appliqués

Initial date of first Declaration of Conformity 2009

Erstmalige Ausstellung der Konformitätserklärung
Année de première Déclaration de Conformité

Issue date: 15 SEP.2017

Operations Manager (Name printed and sign.)

Erstellt am:

Geschäftsführer
Le Directeur des Opérations

Date de délivrance:

ERIC MICHEL 

Quality Manager (Name printed and sign.)

Leiter Qualitätsmanagement
Le Directeur de Qualité

YVAN CHEVREZ 

EF-0912- CE Atitikties Deklaracija

/Invacare logotipas/

Mes su visa atsakomybe pareiškiame, kad gaminys aprašytas žemiau

Gaminio aprašymas, kodas, tipas, modelis

Action 3 NG

Pagamintas

Invacare France Operations
Rout de Saint Roch
37230 Fondettes
Prancūzija

atitinka visas medicinos prietaisų direktyvos 93/42 / EEB VII priedo nuostatas

Prietaiso klasifikacija pagal priedą IX

I klasė

Taikomi suderinti standartai, nacionaliniai standartai ar kiti norminiai dokumentai

EN12183:2014

Pirmos atitikties deklaracijos išdavimo data

2009

Išdavimo data 2017 m. rugsėjo 15 d.

Procesų vadovas Eric Michel /parašas/

Kokybės vadovas Yvan Chevrer / parašas /



Vertimo tikslas:

UAB „TEIDA“
viešoji įmonė, vadybininkė
Miglė Jonikė

2018-05-16

Features and options

stūmimo rankenos 4.




Bi-material cover*

The new black mesh cover has been oversized to cover the back canes better and allows for easy repositioning of the lumbar supporting strap.



Transit chair

The Transit version offers seat heights of 460, 485 and 510 mm with 12" plastic rear wheels.



Tubular armrest

New, slimline, tubular armrests can be selected to give the **Action³ NG** a more active look.




One arm lever drive

This control system allows users with low muscular tone to more easily propel the **Action³ NG** using one arm.



Tie-down hook

A brightly coloured hook allows simple and secure fixation in a vehicle.



Front castor

Front castors available in 5", 6" and 8" pneumatic or solid tyre, with various options of tyre width also available.




Headrests

A neckrest, headrest or cheek support are all available, depending on user needs. These elements can all be adjusted on the backbrace, changing the height or depth to ensure the most optimum position.



12.

pripiūami priekiniai ratai



knija i vestim tber

UAB 'TEIDA' viešųjų pirkimų vadovė Miglė Jonikė

*Action³ NG shown

2019-05-20

Features and options



Backrest types

On top of the fixed backrest (0° or 7°), **Action^{NG}** offers 4 different types of backrests to suit user's needs :

- ① folding,
- ② angle adjustable (-15°/+15°),
- ③ ratchet reclining (0°/30°) and
- ④ gas spring reclining (0°/30°)



sidyuli 6.



Dual handrim propulsion*

Specifically for hemiplegic use, this handrim with QR mechanism makes the chair easily foldable.


ranky abomas 7.



Armrest assortment options


palilianos ranky abomas

① Swing away detachable, ② swing away detachable height-adjustable, ③, ④ and ⑤ removable height-adjustable, ⑥ comfort detachable height-adjustable and ⑦ hemi.



Anti-tip device

Provides extra safety and comes either fixed (top) or removable (bottom) versions.



Transparent tray

Useful for mealtimes and leisure activities such as reading.



Matrix® backrest

Matrix Elite and **Elite Deep** can be added as an option to provide greater comfort.



Buckle belt

This optional **BodyPoint®** belt brings increased user safety whilst seated as well as increased comfort and postural support.

*Action^{NG} shown



UAB 'TEIDA'
Viešykų pirėmių vadybininkė
Miglė Jonikė
2019-05-20


Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	<i>18.</i> Seat width	Seat depth	Seat angle	Seat height	Backrest angle	<i>Support above</i> Backrest height	Front frame angle
Self Propulsion	380 / 405 / 430 / 455 / 480 / 505 mm	400 / 450 / 500 mm	0° / 3° / 6°	360 / 385 / 410 / 435 / 460 / 485 / 510 mm	0° / 7° / 0° - 30° / -15° - +15°	325 - 510 mm	80° / 0 - 80° / 90°
Transit	380 / 405 / 430 / 455 / 480 / 505 mm	400 / 450 mm	0° / 3° / 6°	460 / 485 / 510 mm	0° / 7° / 0° - 30° / -15° - +15°	380 - 460 mm	80° / 0 - 80° / 90°


	Length of footrests	<i>parentis</i> Armrest height	Driving unit width	Total width, folded	Total length inc. footrests	Total length without footrests	Total height
Self Propulsion	200 - 400 / 330 - 460 / 440 - 510 mm	200 - 280 / 240 / 170 - 300 / 210 - 350 mm	Seat width + 190 mm	295 mm	980 - 1105 mm	835 - 960 mm	820 - 1050 mm
Transit	200 - 400 / 330 - 460 / 440 - 510 mm	200 - 280 / 240 / 170 - 300 / 210 - 350 mm	Seat width + 160 mm	260 mm	950 - 1050 mm	710 - 910 mm	920 - 1050 mm

	Max. safe slope	Radius	<i>Naustotaps</i> Max. user weight	Total weight from	Mass of the heaviest part, from
Self Propulsion	9°	870 mm	125 kg	13.0 kg	7.1 kg
Transit	9°	870 mm	125 kg	13.5 kg	10.3 kg



16. *svēta vārdiņa*

Frame colours



Upholstery

Black Nylon



Scan here to access the website.

Please notice colors may vary slightly from those displayed above.

Tel: +41 61 487 70 70
hqeurope@invacare.com
www.invacare.eu.com

Invacare International GmbH
Benkenstrasse 260
4108 Witterswil
Switzerland



Yes, you can.



© 2018 Invacare International GmbH. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. A10022NG - EU - 10/2018

UAB 'EIDA' - Viesuļi pirkumu vadybininke Miglė Jonikė
2019-05-20