



Si - tec Occlusions - spray

occlusionsspray is a universal color indicator to test the occlusal contacts and accuracy of fit of crowns and bridges.



- * Easy to administer and leaves a thin colored film with can easily be removed with water
- * Si-tec occlusions spray can be used for approximal contacts when trial seating crowns and bridges
- * Contains physiologically safe ingredients and is filled with environmentally neutral propellants.
- * Contents 75 ml, green

Si - tec GmbH - Dental attachments

Leharweg 2 - 58313 Herdecke (Germany)
Phone: 0049 (0) 2330/806940
FAX : 0049 (0) 2330/8069420
www.info@si-tec.de e-mail: info@si-tec.de