

cobas 6500 installation kit for LAS **cobas**[®]

cobas 6500 installation kit for LAS

| REF | SYSTEM |
|--------------|---|
| 07438451 001 | cobas u 601 cobas u 701 cobas 6500 |

English

Intended use

The **cobas 6500** installation kit for LAS is intended to make **cobas u 601** urine analyzer or **cobas u 701** microscopy analyzer or the combination of both analyzers ready for uni-directional connection to the laboratory automation systems.

For professional use only.

Precautions and warnings

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory equipments.

Disposal of all waste material should be in accordance with local guidelines.

Handling

The installation, upgrade/connection has to be performed by a Roche field service engineer.

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Storage and stability

Store at 5-40 °C.

Relative humidity: 75 % at 30 °C (non-condensing).

Altitude and pressure: maximum 2000 m above sea level, 80 to 106 kPa

Materials provided

- 1 x input connection unit
- 1 x Water Container (External Water Supply)
- 1 x insulating plate
- 1 x connecting plate
- 1 x inside plate
- 1 x LAS cover plate
- 1 x waste water tube (5 m)
- 1 x gravitational waste water tube (5 m)
- 1 x washer water tube (5 m)
- 1 x LAS cable 8/9 L1000
- 5 x hex socket head cap screw M6 x 16
- 3 x hex socket button head cap screw M4 x 6
- 3 x insulating washer
- 2 x cable tie
- 1 x tube fitting, diameter 6-M5
- 2 x tube fitting, diameter 8-1/8 inch
- 1 x LAS Back Cover Plate
- 1 x connecting plate c8100
- 1 x inside plate c8100
- 1 x gender changer

Materials required (but not provided)

- [REF] 06390498001, **cobas u 601** urine analyzer
- [REF] 06334601001, **cobas u** pack
- [REF] 06390579001, **cobas u** calibration strip
- [REF] 06390501001, **cobas u 701** microscopy analyzer
- [REF] 06390552001, **cobas u** cuvette
- [REF] 06390510001, **cobas 6500** installation kit
- [REF] 07524811001^{a)}, **cobas 6500** table

- [REF] 04590457001^{b)}, Computer Table (PC Table D)

- General laboratory equipment

a) Separately required for **cobas u 601** urine analyzer and **cobas u 701** microscopy analyzer.

b) For optional use.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

For further information, please refer to the appropriate operator's manual for the relevant analyzer and Method Sheets of all necessary components.

Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard:

| | |
|------------|---|
| CONTENT | Contents of kit |
| SYSTEM | Analyzers/Instruments on which reagents can be used |
| REAGENT | Reagent |
| CALIBRATOR | Calibrator |
| → | Volume after reconstitution or mixing |
| GTIN | Global Trade Item Number |

COBAS and COBAS U are trademarks of Roche.

All other product names and trademarks are the property of their respective owners.

Additions, deletions or changes are indicated by a change bar in the margin.

© 2017, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim
www.roche.com

