



Inoveca Medtechna <inoveca@gmail.com>

Dimensions between eyepieces in head

Eugene E.C. Dams <eugene@labomedeuropa.com>
Kam: Inoveca <inoveca@gmail.com>

2020 m. sausio 16 d. 11:31

Dear Kestutis,

Please find attached two pictures of the inside of the Prima DNT ergo head.
If I calculate, the heart distance between the two circles is 22 mm.
Dimension A is 3 mm
Dimension B is 19 mm
Two times half the circle (B) is $2 * 9,5 \text{ mm} = 19 \text{ mm}$ plus A (3 mm) makes 22 mm.
Hopefully this makes sense.

Kind regards / Met vriendelijke groet / Mit freundlichen Grüßen,

Eugene Dams

Sales manager



-----Oorspronkelijk bericht-----

Afzender: Inoveca <inoveca@gmail.com>

Verstuurd: Woensdag 15 Januari 2020 5:30

Aan: Eugene E.C. Dams <eugene@labomedeuropa.com>; Labomed Europe B.V.
<info@labomedeuropa.com>

Onderwerp: Fwd: Dimensions between eyepieces in head

Sorry for my mistake with 23mm.

I ask documentation about STEREO BASE of Prima DNT. Please see Leica M320 - Stereo base 24 mm.

Regards
Kęstutis Sirgedas
Direktorius / CEO



UAB "Inoveca medtechna"

Verkių g.29

LT-09108 Vilnius, Lithuania

Tel. +370 60010001

www.inoveca.lt

facebook.com/InovecaMedtechna



----- Forwarded message -----

Nuo: **Inoveca** <inoveca@gmail.com>


Date: 2020-01-14, an, 17:13

Subject: Dimensions between eyepieces in head

To: Eugene E.C. Dams <eugene@labomedeuropa.com>, Labomed Europe B.V. <info@labomedeuropa.com>

Dear Eugene,

I ask manufacturer technical documents for tender, that distance between eyepieces is 23 mm. see photo.




[Cituojamas tekstas paslėptas]

2 priedai (-ų)



Prima DNT ergo head distance 160120.jpg
90K



Dimensions ergo head 160120.jpg
59K