

## Curettes | Küretten | Curetas | Curettes |

### 12.1. Gracey



Gracey 1/2 • BEUL-05601 • BE-1005001 • LT-1004810 • PremierAir-1004660  
Short shank for incisors and cuspids



Gracey 3/4 • BEUL-05603 • LT-1004811 • PremierAir-1004661  
Short contra-angle shank for incisors and cuspids



Gracey 5/6 • BEUL-05605 • BE-1005010 • LT-1004812 • PremierAir-1004662  
Medium shank for cuspids and molars



Gracey 7/8 • BEUL-05607 • BE-1005015 • LT-1004813 • PremierAir-1004663  
Medium contra-angle shank for bicuspid and molars, buccal and lingual surfaces



Gracey 9/10 • BEUL-05609 • LT-1004814  
Long contra-angle shank for molars with more acute blade angle



Gracey 11/12 • BEUL-05611 • BE-1005020 • LT-1004815 • PremierAir-1004665  
Triple angulation for easy access to buccal/mesial and lingual/mesial surfaces of premolars and molars



Gracey 13/14 • BEUL-05613 • BE-1005025 • LT-1004816 • PremierAir-1004666  
Triple angulation for easy access to distal/posterior surfaces of premolars and molars



Gracey 11/14 • LT-1004817  
For easy access to a facial or lingual sextant with a single instrument.





Gracey 15/16 • BEUL-05653 • BE-1005026 • LT-1004819  
Angled shank, designed to reach mesial surfaces of posterior teeth



Gracey 17/18 • BEUL-05655 • LT-1004820  
Highly angled shank for reaching distal surfaces of posterior teeth

### Gracey-Rigid



After probing, InSite™ scalers provide added visual confirmation that you have reached the base of the pocket.

### Gracey Variations

Pocket Scalers (PS) and Mini Pocket Scalers (MPS) have a slender and longer terminal shank to provide access into deeper pockets. MPS has the same slender tip design as PS, but with half the blade length. InSite PS and MPS have a depth gauge at 5mm through 10mm offering an excellent visual guide.

**Big Easy Ultralite** Pocket Scalers (PS) and Mini Pocket Scalers (MPS) do NOT feature the depth gauge found in InSite Scalers

**Rigid** Gracey designs have a thicker shank than the original Graceys, but are overall slightly shorter for increased strength when removing heavy calculus.



# SCALING AND PERIO→

## WITH CONTROLLED INSTRUMENT MOVEMENTS GENTLE AND PAINLESS

### 12.2.



**A INSTRUMENT**  
For removing supragingival deposits in all quadrants



**DS-001A**



**PS INSTRUMENT**  
For removing calculus in deep periodontal pockets in all quadrants  
“EMS PERIO SLIM TIP HAD BEST INTERPROXIMAL AND SUBGINGIVAL ACCESS” – CRA (Clinical Research Associates) Newsletter, USA, June 1998



**DS-016A**



**B INSTRUMENT**  
For cleaning lingual surfaces with heavy stain and calculus



**DS-003A**



**C INSTRUMENT**  
For removing deposits and stain – especially suitable for removing bonding residue, e.g. in orthodontics



**DS-002A**

## IMPLANT CLEANING WITH PIEZON®





**PI INSTRUMENT**  
For the effective and gentle cleaning of implants and corresponding restorations



**FS-295**




**DT-065A**


**STARTER KIT**  
4 PI instruments  
1 file holder (120°)  
1 flat wrench

**REFILL**  
4 PI instruments






**PL1 INSTRUMENT**  
Left-curved instrument for periodontal debridement – especially suitable for difficult-to-reach interproximal surfaces



**DS-030A**



**PL4 INSTRUMENT**  
Left-curved instrument with small ball (> 0.6-mm diameter) on instrument tip for the periodontal debridement of furcations and concavities



**DS-033A**



**PL2 INSTRUMENT**  
Right-curved. Otherwise same as PL1 instrument



**DS-031A**



**PL5 INSTRUMENT**  
Right-curved. Otherwise same as PL4 instrument




**DS-034A**




**PL3 INSTRUMENT**  
Instrument for cleaning and disinfecting periodontal pockets – especially suitable for follow-up and recall visits




**DS-032A**



**HPL3/DPL3 INSTRUMENTS**  
Diamond-tipped perio instruments for in-depth root debridement using direct vision (> “Flap”) – for smoothing restoration overhangs as well as expanding root furcations



**DS-075A**



**DS-065A**

HPL3 70-µm GRAIN→

DPL3 15-µm GRAIN→

# BacterX®→

READY-TO-USE CHLORHEXIDINE SOLUTION (0,2 %)  
WITH AND WITHOUT ALCOHOL  
TO SUPPORT PERIODONTAL TREATMENT





## SONICflex, prophylaxis tips



**Important information regarding material numbers:**  
**XX:** Tip for SONICflex 2003/2000 (long thread)  
**XXA:** Tip for SONICflex 2008 (short thread)

### SONICflex for prophylaxis

Thanks to its numerous tips, SONICflex is a perfect partner for prophylaxis purposes. All the tips available, which can be used for a wide range of indications, can be found on pages 38–41.



**Scaler tip no. 5 Universal**  
 5: 0.571.5171 / 5A: 1.005.8949



**Scaler tip no. 6 sickle**  
 6: 0.571.5181 / 6A: 1.005.8950



**Para tip no. 60 | long, straight tip, rounded surface**  
 60: 0.571.7402 / 60A: 1.006.1934



**Para tip no. 61 | angled tip, left, rounded surface**  
 61: 0.571.7412 / 61A: 1.006.1935



**Scaler tip no. 7 perio**  
 7: 0.571.5191 / 7A: 1.005.8951



**Scaler tip no. 81 perio, extra long**  
 8: 0.571.5371 / 8A: 1.006.1953



**Para tip no. 62 | angled tip, right, rounded surface**  
 62: 0.571.7422 / 62A: 1.006.1936



**Implant Set | wedge-shaped tips, single-use plastic head**  
 Set: 1.003.8167 / Set A: 1.006.2027  
 Implant Refill (30 pcs): 1.003.8168

## PiezoLED, prophylaxis tips

### PiezoLED for prophylaxis

PiezoLED is fullyadjustable to suit individual requirements, which also makes it ideal for prophylaxis purposes. An overview of the various tips associated with its broad range of applications can be found on page 42–43.




**PIEZO Scaler Tip 201**  
 For removal of supragingival calculus in all quadrants  
 Mat. no. 1.007.4024



**PIEZO Scaler Tip 202**  
 For removal of supra- and subgingival calculus in all quadrants up to a pocket depth of 4 mm  
 Mat. no. 1.007.4026



**PIEZO Scaler Tip 203**  
 For removal of subgingival deposit on root surfaces. Also suitable for deep periodontal pockets and periodontal treatments/chest-ups  
 Mat. no. 1.007.4028



## Cleaning powder for subgingival and supragingival treatment

PROPHYpearls, PROPHYflex Powder and PROPHYflex Perio Powder are perfectly suited to the KaVo PROPHYflex 4 and allow for tooth cleaning that is both thorough and gentle.



### 12.3.

#### KaVo PROPHYpearls

For an all-round clean and gentle prophylaxis

- Indication: effective elimination of discolouration, plaque and biofilm
- Active ingredient: calcium carbonate with a particle size of 60–70 µm, biodegradable, no salt, sugar or gluten
- Highlights: round shape and a uniform size for optimum cleaning with minimal abrasion, 4 flavours and 2 neutral variants

**Sticks (80 pcs per 15 g)**  
 Neutral: Mat. no. 1.070.1826  
 Mint: Mat. no. 1.070.1828  
 Peach: Mat. no. 1.070.1829  
 Orange: Mat. no. 1.070.1830  
 Blackcurrant: Mat. no. 1.070.1831  
 Bottles (4 pcs per 250 g)  
 Neutral: Mat. no. 1.070.1798

#### KaVo

#### PROPHYflex Powder

Polishing can be so quick.

- Indication: highly effective elimination of discolouration, plaque and biofilm
- Active ingredient: sodium bicarbonate with a particle size of 50–70 µm, water-soluble
- Highlights: simple, fast and highly effective treatment with 4 flavours

**Sticks (80 pcs per 15 g)**  
 Orange: Mat. no. 1.007.0014  
 Berry: Mat. no. 1.007.0015  
 Cherry: Mat. no. 1.007.0016  
 Mint: Mat. no. 1.007.0017

#### KaVo

#### PROPHYflex Perio Powder

For subgingival and supragingival treatment

- Indication: efficient biofilm management in periodontal pockets and clearing of crowns, bridges and implants
- Active ingredient: fine glycine powder (18–22 µm), protects gums
- Highlights: very gentle and not too abrasive, suitable for efficient cleaning of pockets as well as supragingival applications on implants


**Bottles (4 pcs per 100 g)**  
 Mat. no. 1.009.3732

#### Perio Kit for PROPHYflex 4


10 perio tips for safe and targeted pocket cleaning thanks to the stable, thintip\* with depth markings and Perio Adapters as well as PROPHYflex Perio Powder and accessories  
 Perio Kit: Mat. no. 1.071.9403\*\*  
 Refill Set: Mat. no. 1.070.0287\*\*




#### Supragingival



#### Subgingival



#### Subgingival



\* Each tip can be reprocessed up to 10 times.  
 \*\* Use other mat. no. for PROPHYflex 3.

# AIR N GO

EASY



## 12.4.

### AIR-N-GO® easy

The new generation of dual purpose Air Polisher: SUPRA and PERIO

- AIR-N-GO® easy with KaVo® adapter **F10100**
- AIR-N-GO® easy with Sirona® adapter **F10101**
- AIR-N-GO® easy with W&H® adapter **F10102**
- AIR-N-GO® easy with Bien-Air® adapter **F10103**
- AIR-N-GO® easy with Midwest® adapter **F10105**

Delivered with a convertible AIR-N-GO® easy handpiece, a SUPRA nozzle, 1 intro kit of 20g bags of «CLASSIC» powders, 1 AIR-N-GO® easy maintenance kit (replacement o'rings, syringe of silicone grease with brushes, cannula, cleaning drill).

\*KaVo, Sirona, W&H, Bien-Air and Midwest are registered trademarks.

Class IIa medical devices (GMED) - CE 0459



AIR-N-GO® easy brochure

## AIR-N-GO® EASY ACCESSORIES



### AIR-N-GO® easy PERIO KIT

**F10120**

Turns your handpiece into a Perio system to treat periodontal disease (tooth & implants) on pockets from 3 to 8mm.

Delivered with a AIR-N-GO easy nozzle and a bottle of AIR-N-GO PERIO glycine powder.

Class IIa medical devices (GMED) - CE 0459



## 12.5.

### AIR-N-GO® SUPRA nozzle

**F10125**

For all prophylaxis supra-gingival treatments.



### AIR-N-GO® PERIO easy nozzle

**F10126**

Turns your handpiece into a Perio system to treat non surgical periodontal disease (tooth & implants) in pockets from 3 to 8mm.



### AIR-N-GO® PERIO nozzle

**F10127**

Turns your handpiece into a Perio system to treat non surgical and surgical periodontal disease (tooth & implants) in pockets from 8 to 10mm.



### AIR-N-GO® PERIO Maintenance nozzle

**F10128**

Turns your handpiece into a system adapted for periodontal maintenance in pockets until 4mm.



### AIR-N-GO® easy maintenance kit

**F10121**