


3 pirkimo dalis- Nosies švirkštas

LMA MAD Nasal™ Needle-free intranasal drug delivery device.

The unique design of LMA MAD Nasal™ enables a seal to be formed with the nostril and atomized particles of medication to be delivered to the nasal mucosa. The particles are then absorbed through the mucosal tissue and into the bloodstream, avoiding first pass metabolism.



Time from presentation to the paediatric ER department to analgesia was significantly reduced with IN vs IV drug delivery ($p < 0.000$).¹

IN vs IV drug delivery

Within a few minutes of delivery, serum levels of some medications delivered intranasally are comparable to injectable levels.²

IN vs PR drug delivery

With LMA MAD Nasal™, seizure medication can be administered intranasally:

- ▶ Significantly reduces the likelihood of recurrent seizures in the ED³
- ▶ Significantly reduces the need for additional medication to treat ongoing seizures³
- ▶ Avoids the embarrassment associated with rectal drug delivery²

IN vs IM & Oral drug delivery

Utilising low molecular weight lipophilic drugs intranasally ensures that molecules are easily absorbed across the nasal mucosa for comparable efficacy vs IM and oral drug delivery.^{4,5,6}

“The intranasal route of administration presents several advantages in small children. It is easily administered to a reluctant child, cannot be spit out and allows a more rapid absorption than when given orally.”

Talon M. et al., 2009

Vaistai įpurškiami dulksnos pavidalu

minkštas kūgio formos kamštis

Soft conical plug

The plug forms a seal with the nostril preventing expulsion of fluid.

Malleable stylet

The malleable stylet allows 180° positioning of the nasal plug.

Atomization spray

The spray atomizes drugs into a fine mist of particles 30-100 microns in size.⁴

Pressure

High applied pressure ensures that drugs are atomized into a fine mist of particles through the tip of the plug.

lankstus stiletas,
keičia padėtį
180 laipsnių kampu

Accurate dosing

The syringe enables the accurate measurement of drugs to be delivered.

LMA MAD Nasal™: Product specification

PO number	Quantity / case	Quantity / box	Description
MAD100	100	25	Intranasal mucosal atomization device with 3 ml syringe
MAD100OS	100	25	Intranasal mucosal atomization device with 3 ml oral syringe
MAD110	100	25	Intranasal mucosal atomization device with 1 ml syringe
MAD110OS	100	25	Intranasal mucosal atomization device with 1 ml oral syringe
MAD130	100	25	Intranasal mucosal atomization device with 1 ml syringe and vial adaptor
MAD130OS	100	25	Intranasal mucosal atomization device with 1 ml oral syringe and vial adaptor
MAD140	100	25	Intranasal mucosal atomization device with 3 ml syringe and vial adaptor
MAD140OS	100	25	Intranasal mucosal atomization device with 3 ml oral syringe and vial adaptor
MAD300	100	25	Intranasal mucosal atomization device without syringe

Distributed by:

TELEFLEX HEADQUARTER INTERNATIONAL, IRELAND
Teleflex Medical Europe Ltd., IDA Business and Technology Park,
Dublin Road, Athlone, Co Westmeath, Ireland
Phone +353 (0)9 06 46 08 00 Fax +353 (0)14 37 07 73
orders.intl@teleflex.com



Manufactured by:
Wolfe Tory Medical, Inc.
79 West 4500 South
Salt Lake City, Utah 84107



Consult IFU on this website:
www.LMACO.com



Printed on paper from
sustainable sources

Teleflex www.teleflex.com LMA International www.lmaco.com Germany www.lma.de USA www.lmana.com
Italy www.lmaco.it China www.lmaco.cn Canada www.lmavitaid.com Australia & New Zealand www.lmapacmed.com

Copyright © 2013 Teleflex Incorporated. All rights reserved. LMA, LMA MAD, LMA MAD Nasal, MAD, MAD Nasal and LMA Better by Design are trademarks or registered trademarks of Teleflex Incorporated or its affiliates.

94 06 95 - 00 00 01 - REV 05/13 - LMA / PDF - 11 13 PDF

References:

1. Borland M. et al. Emergency Medicine Australasia 2008; **20**: 515-520.
2. Wolfe T.R., Macfarlane T.C. American Journal of Emergency Medicine 2006; **24**: 343-346.
3. Holsti M. et al. Pediatric Emergency Care 2007; **23** (3): 148-153.
4. Talon M. et al., J Burn Care Research 2009; **30**: 599-605.
5. Pires A. et al. J Pharm Pharmaceut Sci 2009; **12** (3): 288-311.
6. Ahmad S. et al. Lancet 2006; **367**: 1591-1597.



For the latest
clinical evidence on
LMA MAD Nasal™

www.doctorevidence.com/lma



For the latest
educational videos
and clinician
testimonials

[www.youtube.com/
LaryngealMaskAirway](http://www.youtube.com/LaryngealMaskAirway)



For the latest
news from LMA

[www.facebook.com/
LMAInternational](http://www.facebook.com/LMAInternational)



For product
information and
access to product
instructions for use

www.lmaco.com



Mediware Disposable syringes with Luer-Lock

- The high quality 3-pieced special Luer-Lock syringe
- Without cannula
- Crystal-clear
- With extended, non-erasable scale
- Sterile and singly packed

Description for Mediware Disposable syringes with Luer-Lock

- The high quality 3-pieced special Luer-Lock syringe
- Without cannula
- Crystal-clear
- With extended, non-erasable scale
- Sterile and singly packed

3 ml švirkštas

Volume	PZN	SU	SU	REF
2,5/3 ml		3000 pcs.	100 pcs.	I3 0501
5 ml		2000 pcs.	100 pcs.	I3 0502
10 ml		1200 pcs.	100 pcs.	I3 0503
20 ml		900 pcs.	100 pcs.	I3 0504
30 ml		450 pcs.	50 pcs.	I3 0505
50/60 ml		400 pcs.	50 pcs.	I3 0506

Nosies švirkštas, skirtas leisti vaistus

Unikalus LMA MAD Nasal™ dizainas leidžia užkimšti šnervę ir į nosies gleivinę įpurškšti vaistų daleles, kurios yra sugeriamos per gleivinę, išvengiant pirmojo prasiskverbimo efekto.

Suskaldytų dalelių įpurškimas

Vaistai yra įpurškiami 30-100 mikronų dydžio dalelėmis dulksnos pavidalu.¹

Minkštas kūgio formos kamštis

Kamščiu užkemšama šnervė, taip išvengiant skysčių išteklėjimo.

Lankstus stiletas

Lankstusis stiletas leidžia nosies kamščio padėtį keisti 180° kampu.

Tikslus dozavimas


Švirkštu galima tiksliai išmatuoti įpurškiamų vaistų kiekį.

"Kraujagyslėmis išvagota nosies gleivinė ir uoslės audiniai, turintys tiesioginį ryšį su centrine nervų sistema, leidžia per nosį įleidžiamus vaistus greitai pernešti į kraujotaką ir smegenis, šiuo atžvilgiu beveik prilygstantis intraveninei terapijai."

Wolfe T. & Braude D., 2010

Slėgis

Taikomas aukštas slėgis užtikrina, kad vaistai per kamščio galiuką būtų suskaldyti į smulkias daleles dulksnos pavidalu.



Bone and Bone Marrow Biopsy Device 8 pd- Kaulų čiulpų infuzinė biopsinė adata 18 G



Adata smailiu galu

MYELO SLIM Qty/Box 10pcs

Add **E** for Explant version (VVI1E)

CODE	SIZE
VVI11410/50	14Gx10/50mm
VVI11430/70	14Gx30/70mm
VVI11510/50	15Gx10/50mm
VVI11530/70	15Gx30/70mm
VVI11610/50	16Gx10/50mm
VVI11630/70	16Gx30/70mm
VVI11810/50	18Gx10/50mm

On request Vigeo is available to provide needles with different diameter and length


Adatos ilgis yra reguliuojamas nuo 10 iki 50mm

- Cannula and stylet designed to allow an easier penetration of the bone cortex
- Adjustable depth stop for sternal or paediatric application, which can be removed for iliac crest use
- Luer slip and luer-lock syringes fit securely

Galimybė prijungti švirkštą su Luer Lock ir lygiavertę jungtimi

MYELO

CLINICAL ONCOLOGY



Ergonomic handle

MYELO Qty/Box 10pcs

Add **E** for Explant version (VVIE)

CODE	SIZE
VVI1410/50	14Gx10/50mm
VVI1430/70	14Gx30/70mm
VVI1510/50	15Gx10/50mm
VVI1530/70	15Gx30/70mm
VVI1610/50	16Gx10/50mm
VVI1630/70	16Gx30/70mm
VVI1810/50	18Gx10/50mm

On request Vigeo is available to provide needles with different diameter and length

Adatos kaulų čiulpų biopsijai

MYELO

Skurtos paprastesnei kaulų penetracijai

Nuimamas ilgių reguliatorius

Kodas VVI

Nuimamas reguliatorius be sriegio

Gylio ribotuvas

CM ženkliniai




Atraumatinis briaunotas adatos galiukas

Kodas	Dydis	Pak.
VVI1435/100	14Gx35-100 mm	20
VVI1425/90	14Gx25-90 mm	20
VVI1410/75	14Gx10-75 mm	20
VVI1535/100	15Gx35-100 mm	20
VVI1525/90	15Gx25-90 mm	20
VVI1510/75	15Gx10-75 mm	20
VVI1625/90	16Gx25-90 mm	20
VVI1610/75	16Gx10-75 mm	20
VVI1810/75	18Gx10-75 mm	20

12pd- Tracheostominių vamzdelių fiksatorius: mažas, vidutinis, didelis

Product features














- For securely holding a tracheostomy tube, a laryngectomy tube or a Stoma button with eyelets in place
- Provide an optimum wearing comfort, as they are made of soft, skin-compatible material (soft fabric, foam, Polyurethan (PU))
- Available in different widths and with different fasteners (hook and slide (narrow and wide), velcro fastener, velcro/ hook)
- Individually adjustable
- Variants for adults and children
- PRIMEDIFIX® with hook enables fast fixation to the eyelets of a tube flange or Stoma Button (A-E)
- PRIMEDIFIX® with velcro fastener is especially suitable for sensitive skin (F-I)



- EASYFIX® with Velcro fastener on one side and a hook on the other enables even motor impaired patients (e.g. after a stroke) to easily fix the neck strap to a tube or stoma button, even with just one hand (J-M)


Indication

- Tracheostomy/ Laryngectomy
- Radiation/ MRI

		PU	Article number
A	PRIMEDIFIX® neck strap, 1-piece, white, with narrow hook, width = 12 mm, length = 26-47 cm, non sterile		10 pieces 200229
B	PRIMEDIFIX® neck strap, 1-piece, skin-colored, with narrow hook, width = 12 mm, length= 26-47 cm, non sterile		10 pieces 200227
C	PRIMEDIFIX® neck strap, 1-piece, transparent, with wide hook, Polyurethan (PU), width = 10 mm, length = 27-50 cm, non sterile		10 pieces 200219
D	PRIMEDIFIX® neck strap, 1-piece, skin-colored, with wide hook, width = 12 mm, length= 26-47 cm, non sterile		10 pieces 200226
E	PRIMEDIFIX® neck strap, 1-piece, white, with wide hook, width = 12 mm, length= 26-47 cm, non sterile		10 pieces 200228
F	PRIMEDIFIX® neck strap for children, 2-part, white, with velcro fastener width = 20 mm, length = 24 cm, non sterile		10 pieces 200203
G	PRIMEDIFIX® neck strap, 2-part, white, with velcro fastener width = 24 mm, length = 47 cm, non sterile 2 dalių, su Velcro užsegimais ilgis, 47cm		10 pieces 200211
H	PRIMEDIFIX® neck strap, 2-part, white, with velcro fastener width = 16 mm, length = 45,5 cm, sterile		10 pieces 21451
I	PRIMEDIFIX® neck strap, 2-part, white, with velcro fastener, width = 26 mm, length = 45,5 cm, sterile		10 pieces 21450
J	EASYFIX® neck strap, 1-piece, skin-colored, with wide hook and velcro fastener, width = 12 mm, length = 26-47 cm, non sterile		10 pieces 200216
K	EASYFIX® neck strap, 1-piece, white, with wide hook and velcro fastener, width = 12 mm, length = 26-47 cm, non sterile		10 pieces 200217
L	EASYFIX® neck strap, 1-piece, skin-colored, with narrow hook and velcro fastener, width = 12 mm, length = 26-47 cm, non sterile		10 pieces 200214
M	EASYFIX® neck strap, 1-piece, white, with narrow hook and velcro fastener, width = 12 mm, length = 26-47 cm, non sterile		10 pieces 200215

1. Product Name and Product Code:

AL-1700.V005- Neonatal Anesthesia Circuit Corrugated PE Tube
Kvėpavimo sistema naujagimiams

**2. Manufacturer:**

Altera Tibbi Malzeme Sanayi ve Ticaret A.Ş.

3. Sterilization Method :

Ethylene oxide sterilization. This product line is provided both as sterile and nonsterile.

4. Packaging:

Nylon packaging. ISO Class 8 Clean Room Manufactured



PRODUCT DESCRIPTION

PD08-03

5. Intended Use and Functional Description :

Disposable circuits used to conduct medical gases from the anaesthesia system to the patient.

Anaesthetic gas exits the anaesthesia machine (via the common gas outlet) and then enters a breathing circuit. The function of the circuit is to deliver oxygen and anaesthetic gases to the patient and to eliminate CO₂. The CO₂ may be eliminated by gas inflow or by soda lime absorption.

6. Shelf life :

5 years from the date of production .

7. Product Description:

7.1. Classification: Class IIa Rule 2

7.2. This product does not contain any metallic parts

7.3. Spesification of the product :

Component PN	Component Description	Qty.	UoM
150.11.036	ELBOW CONNECTOR LUER CAPalkūninė jungtis su dangteliu	1	Pcs.
150.11.226	NEONATAL RED CAP	1	Pcs.
151.03.128	10M/22 F STRAIGHT CONNECTOR tiesi jungtis	2	Pcs.
151.03.140	10M-10F STRAIGHT CONNECTOR tiesi jungtis	2	Pcs.
151.03.177	NEONATAL Y CONNECTOR ANGLED WITH CO ₂ PORT PP	1	Pcs.
151.05.01.300	10 MM CORRUGATED TUBING PE 160 CM CUT	2	Pcs.

Jungtus su CO₂ anga

Sistemos ilgis 160cm , 10mm diametras (2 gofruopti vamzdžiai)

8. Standarts:

EN ISO 13485:2003, 93/42/EEC Medical Device Directive, EN ISO 5367, EN ISO 5356-1, EN ISO 11607-1:2009, EN ISO 11607-2:2009, EN ISO11135-1:2007, EN ISO 23328-1, EN ISO 23328-2, EN ISO 594-1

9. Waste Method:

After use product will be “contaminated medical waste” and package will be “packaging waste” so they should be handled according to relevant national and international standards and regulations.

Date: 19.03.2020

STATEMENT

To whom it may concern,

We are producing different modifications of breathing circuits in accordance with international standards in the direction of customer demands.

MEDİTERA TİBBİ MALZEME SAN. VE TİC. A.Ş.
İbni Melek OSB Mah. TOSBİ YOL 4 Sok. No:29 35900 Tire-İZMİR
Tire V.D: 060 029 8703 Tic. Sicil No: 2346
Mersis No: 0060 0298 7030 0016
Tel: 0 232 513 51 10 Fax: 0 232 513 51 14

KONFIDENCIALU

/Altera logotipas/

UAB Skirgesai, Lietuva
Deklaraciją pagal užklausimą dėl produktų modifikacijų

2020 03 19

Visiems suinteresuotiems asmenims

Mūsų kompanija gamina įvairių modifikacijų kvėpavimo sistemas, kurios atitinka tarptautinius standartus ir atskirų klientų poreikius.

Linkėjimai
Guliz SARAC

/parašas/

Produktų atitikties vadovas
Altera Tibbi Malz. San. Ve Tic. A.S.

/Įmonės rekvizitai/



REF 450-P



TRACOE vario Tracheostomy Tube

with atraumatic insertion system, spiral-reinforced, adjustable neck flange, low pressure cuff, scale and 15 mm connector

- Including:
- ◆ lubricating gel
 - ◆ obturator and neck strap
 - ◆ in sterile package

As percutaneous set in REF 420

22pd- Tracheostominiai vamzdeliai su manžete Trachoe Vario



REF 460-P



TRACOE vario Tracheostomy Tube

with atraumatic insertion system, adjustable neck flange, low pressure cuff, X-ray contrast line, scale and 15 mm connector

- Including:
- ◆ lubricating gel
 - ◆ obturator and neck strap
 - ◆ in sterile package

As percutaneous set in REF 423

su
atraumatinė
įvedimo
sistema,
žemo slėgio
manžetė,
rentgeno
kontrastinis

obturatorius, kaklo raištis



REF 470-P



TRACOE vario extract Tracheostomy Tube

with atraumatic insertion system, adjustable neck flange, low pressure cuff, subglottic suction line, X-ray contrast line, scale and 15 mm connector

- Including:
- ◆ lubricating gel
 - ◆ 2 connectors for suction pump
 - ◆ obturator and neck strap
 - ◆ in sterile package

As percutaneous set in REF 422



REF 464




TRACOE vario Tracheostomy Tube

with adjustable neck flange, fenestrated, X-ray contrast line, scale and 15 mm connector

- Including:
- ◆ obturator and neck strap
 - ◆ in sterile package

Technical Data


TRACOE vario Techniniai duomenys, išmatavimai



REF 450, REF 455, **REF 460** *in sizes 06 – 11*
 REF 464 *in sizes 07 – 09*
 REF 470 *in sizes 07 – 10*
 REF 450-P, REF 460-P *in sizes 07 – 09*
 REF 470-P *in sizes 08 – 09*

Size	ID mm	OD mm	A mm	C mm	B mm	TL mm	OB mm	BA °	CD mm
06	6.0	8.2	30	28	12	70	76	100	23
07	7.0	9.7	33	32	15	80	86	100	26
08	8.0	11.2	35	35	20	90	97	100	28
09	9.0	12.3	38	35	27	100	106	100	30
10	10.0	13.7	38	37	35	110	120	100	32
11	11.0	14.8	38	37	35	110	121	100	34

TRACOE vario XL



REF 451, REF 461, REF 471 *in sizes 07 – 10*
 REF 451-P, REF 461-P *in sizes 07 – 09*
 REF 471-P *in sizes 08 – 09*

Size	ID mm	OD mm	A mm	C mm	B mm	TL mm	OB mm	BA °	CD mm
07	7.0	9.7	35	36	33	104	110	100	23
08	8.0	11.2	37	43	40	123	128	100	28
09	9.0	12.3	38	43	45	126	134	100	30
10	10.0	13.7	38	52	45	135	146	100	32

ID: inside diameter (clear width) at bottom of cannula; **OD:** outside diameter at bottom of cannula; **A:** length; **C:** length;
B: length; **TL:** length along centre line; **OB:** length along outer bend; **BA:** bending angle; **CD:** cuff diameter

To order, specify **REF + size**, e. g. **REF 450-08**

25 pirkimo dalis: Neinvazinis endolaringinis švirkštas ir

26 pirkimo dalis: Neinvazinis pediatriškas nosies/burnos gleivinių
aerozolinis švirkštas

MADgic

Laryngo-Tracheal Mucosal Atomization Device
To support smooth emergence from anesthesia,
difficult and awake intubations.

A versatile choice to administer topical anesthesia* to the upper airway.

The MADgic Laryngo-Tracheal Mucosal Atomization Device from Teleflex provides clinicians with the ability to deliver atomized topical anesthesia directly to the airway mucosa. This helps to suppress the airway responses to emergence phenomena such as coughing or laryngospasm^{1,2} and to reduce the cardiovascular responses to mechanical stimulation of the airway³. Topical anesthesia administered by the MADgic Device supports increased patient compliance during difficult and awake intubations⁴.

Coverage, Comfort, Convenience

The MADgic Device allows for optimal coverage by enabling a gentle dispersion of atomized particles across a broad area of the airway mucosa. The malleable stylet allows for targeted

delivery of medication to the deeper pharyngeal structures and to the glottis⁵.

The MADgic Device provides enhanced patient comfort by allowing less stressful intubation and extubation⁶. It also helps to reduce the risk of patient injury by featuring a convenient narrow profile designed to comfortably fit through the vocal cords, into the nasal cavity, or down a supraglottic airway device.

Safety

The single-use MADgic Device enhances patient safety by helping to reduce the risk of cross-contamination while also helping to reduce the risk of toxicity by facilitating exact dosing.

Benefits



The Clinician

Help in managing emergence issues and difficult airway intubation.



Your Institution

Single-use reduces risk of cross-contamination.




The Patient

Increase patient comfort and reduced risk of patient injury.

* For use with drugs approved for intranasal and oropharyngeal delivery.

Teleflex



To learn more visit Imaco.com/products

MADgic Device Specifications

Typical Particle Size	30–100µm	vaistų įpurškimas mikronais
System dead space	0.12 ml MAD720/0.27 ml MAD600/0.18 ml MAD700	
Tip Diameter	0.19" (4.8 mm)	
Applicator Length	4.9" (12.4 cm) MAD720/8.9" (22.6 cm) MAD600 and MAD700	
Stileto ilgis		

MADgic Laryngo-Tracheal Mucosal Atomization Device

ITEM NUMBER	PRODUCT DESCRIPTION	QTY/BOX
MAD600	MADgic Device with 3 ml Syringe	25
MAD700	MADgic Device without ml Syringe	25
MAD720	MADgic Device without 3 ml Syringe	25

References

1. Mihara T, Uchimoto S, Morita S, Goto T. The efficacy of lidocaine to prevent laryngospasm in children: a systematic review and meta-analysis. *Anaesthesia* 2014; 69, 1388–1396.
2. Diachun CA, Tunink BP, Brock-Utne JG. Suppression of cough during emergence from general anesthesia: laryngotracheal lidocaine through a modified endotracheal tube. *J Clin Anesth.* 2001 Sep;13(6):447-51.
3. Hamaya Y., Dohi S. Differences in Cardiovascular Response to Airway Stimulation at Different Sites and Blockade of the Responses by Lidocaine. *Anesthesiology* 2000; 93:95-103.
4. Lohse J.A., Piepho T., Noppens R.R. Awake video laryngoscopy – A revolution in the management of the anticipated difficult airway? *Trends in Anesthesia and Critical Care.* Feb 2016 1-6.
5. Doyle J. Airway Anesthesia, Theory and Practice. *Anesthesiology Clinics* 33 (2015) 291-304.
6. Leung Y., Vacanti FX. Awake without complaints: maximizing comfort during awake fiberoptic intubation. *J Clin Anesth.* 2015 Sep;27(6):517-9

Teleflex, the Teleflex logo, MAD and MADgic are trademarks or registered trademarks of Teleflex Incorporated or its affiliates, in the U.S. and/or other countries. Information in this document is not a substitute for the product Instructions for Use. The products in this catalogue may not be available in all countries. Please contact your local representative. All data current at time of printing (07/2016). Subject to technical changes without further notice.

© 2016 Teleflex Incorporated. All rights reserved. 94 06 96 - 00 00 01 - REV A - 07 16 PDF

Teleflex

Distributed by:

Teleflex Headquarters International, Ireland · Teleflex Medical Europe Ltd. · IDA Business & Technology Park
Dublin Road · Athlone · Co Westmeath · Tel. +353 (0)9 06 46 08 00 · Fax +353 (0)14 37 07 73 · orders.intl@teleflex.com
United Kingdom Tel. +44 (0)14 94 53 27 61 · info.uk@teleflex.com
South Africa Tel. +27 (0)11 807 4887 · assist.africa@teleflex.com



Mediware Disposable syringes with Luer-Lock

- The high quality 3-pieced special Luer-Lock syringe
- Without cannula
- Crystal-clear
- With extended, non-erasable scale
- Sterile and singly packed

Description for Mediware Disposable syringes with Luer-Lock

- The high quality 3-pieced special Luer-Lock syringe
- Without cannula
- Crystal-clear
- With extended, non-erasable scale
- Sterile and singly packed

3 ml švirkštas

Volume	PZN	SU	SU	REF
2,5/3 ml		3000 pcs.	100 pcs.	I3 0501
5 ml		2000 pcs.	100 pcs.	I3 0502
10 ml		1200 pcs.	100 pcs.	I3 0503
20 ml		900 pcs.	100 pcs.	I3 0504
30 ml		450 pcs.	50 pcs.	I3 0505
50/60 ml		400 pcs.	50 pcs.	I3 0506

Strength in Details

29 pd- Sistema invazinio kraujo spaudimo matavimui

Clinically, blood pressure is of fundamental importance in the physiological parameter for humans. For critical patients, a subtle change in blood pressure may be the harbinger of a developmental trend. Therefore, Antmed uses the world’s leading pressure chip, combined with our high-end manufacturing processes, because, we understand, the slightest change in blood pressure may be an indication of something critical.



Vienkartinė sistema invazinio kraujo spaudimo matavimui

: Disposable Pressure Transducers



• Single Channel DPT



• Double Channel DPT



• Triple Channel DPT



• Single Channel DPT (Precision Flush Device)



• Single Channel DPT (Arterial and Venous pipeline)



• Single Channel DPT (Self-flush valve)

: Blood Sampling System


Praplovimo vožtuvas

Self-Flush Blood Sampling Interface

Self-flush after blood sampling can enhance the efficiency of clinical work. Needle-free connector can avoid cross infection

Negative Pressure Suction Device



Negative pressure suction device is used for rapid blood sampling and propelling



Precision Flush Device

High precision flow control and good product durability assure long-term usage effectiveness

Single Channel Transducer

P/N	Flow rate	Description	Specification	Diagram
PT111101	3ml/hr	1-infusion set 1-core part	30pcs/case	
PT111201	30ml/hr	1-flush device 1-stopcock		
PT111103	3ml/hr	1-infusion set 1-core part 1-flush device		
PT111203	30ml/hr	2-stopcock 1-120cm Pressure connecting tube 1-30cm Pressure connecting extension tube		
PT111105	3ml/hr	1-infusion set 1-core part 1-flush device		
PT111205	30ml/hr	2-stopcock 1-150cm Pressure connecting tube 1-30cm Pressure connecting extension tube		
PT121101	3ml/hr	1-infusion set 1-core part 1-flush device	30pcs/case	
PT121201	30ml/hr	1-stopcock		
PT121103	3ml/hr	1-infusion set 1-core part 1-flush device		
PT121203	30ml/hr	2-stopcock 1-120cm Pressure connecting tube 1-30cm Pressure connecting extension tube		
PT121105	3ml/hr	1-infusion set 1-core part 1-flush device		
PT121205	30ml/hr	2-stopcock 1-150cm Pressure connecting tube 1-30cm Pressure connecting extension tube		


Į rinkinį įeina:

- Infuzinė sistema
- Praplovimo vožtuvas (3 ml/h greičiu)
- Daviklis
- 3 krypčių kranelis
- Aukšto slėgio linija
- Aukšto slėgio prailginimo linija

Single/Double/Triple available

技术参数:

- 压力电阻硅芯片, 符合美国AAMI BP22标准
- +/-20mmHg零压输出
- -50~300mmHg使用压力范围
- 4.95~5.05uV/V/mmHg灵敏度
- 168小时使用寿命保证
- 输入电压: 2~10伏直流
- 保存温度: -18~60℃
- 操作温度: +10~+40℃
- 过压保护: -400~6000mmHg
- 3毫升、30毫升/小时或无冲洗阀, 三种选择
- 环氧乙烷灭菌



SPECIFICATIONS:

- Fully piezoresistive silicon pressure sensor, AAMI BP22 compliant
- +/-20mmHg Offset
- -50~300mmHg operating pressure range, Customized pressure range
- 4.95~5.05uV/V/mmHg sensitivity
- 168 hours operating life
- Supply Voltage: 2~10 VDC
- Storage Temperature: -18~60℃
- Operating Temperature: +10~40℃
- Overpressure: -400~6000mmHg
- 3cc, 30cc/hr or without flush device
- ETO sterilization

连续 精确 安全 可靠 简便
Consistent Accurate Safe Reliable Easy-to-Use



antmed
Partner For Life


Shenzhen Antmed Co., Ltd.

18 Jinhui Ave., Pingshan New District, Shenzhen 518122, China
Tel: +86 755 8606 0992 Fax: +86 755 8606 0991
E-mail: info@antmed.com http://www.antmed.com

30pd- Epidurinė adata Tuohy tipo 18 G 120 mm

EpiLong Tuohy

The Tuohy cannula in this kit is equipped with a plastic stylet, making it correspondingly less expensive. The standard epidural catheter is optionally available with closed tip and three lateral openings or with a central opening.



Epidural anaesthesia

Kit components

Tuohy cannulas

Product	Size	Item no.	PU
TuohySono			
with metal stylet and wings	17 G x 90mm	1085-4M090	25
	18 G x 50mm*	1085-4K050	25
	18 G x 90mm	1085-4K090	25
	18 G x 120mm	1085-4K120	10
Tuohy			
with metal stylet and wings	21 G x 90mm	1150-4F090	25
	19 G x 50mm	1150-4I050	25
	19 G x 80mm	1150-4I080	25
	19 G x 90mm	1150-4I090	25
	18 G x 50mm*	1150-4K050	25
	18 G x 80mm	1150-4K080	25
	18 G x 90mm	1150-4K090	25
	18 G x 120mm	1150-4K120	10
	18 G x 150mm	1150-4K150	10
	17 G x 80mm	1150-4M080	25
	17 G x 90mm	1150-4M090	25
	17 G x 120mm	1150-4M120	10
Dydis	17 G x 150mm	1150-4M150	10
	16 G x 80mm	1150-4O080	25
	16 G x 90mm	1150-4O090	25
	16 G x 120mm	1150-4O120	10
	19 G x 90mm	1150-6I091	25
	18 G x 80mm	1150-6K081	25
Tuohy	18 G x 90mm	1150-6K091	25
	17 G x 80mm	1150-6M081	25
	17 G x 90mm	1150-6M091	25
	16 G x 90mm	1150-6O091	25
	with plastic stylet and wings		


* with 0.5 cm graduation



TuohySono cannula with metal stylet




Tuohy cannula with metal stylet



Tuohy cannula with plastic stylet

Crawford cannulas

Product	Size	Item no.	PU
CauDal			
with 45° grinding, wings and graduated depth stop	23 G x 25mm	1149-12D025	25
	22 G x 40mm	1149-12E040	25



Crawford cannula

spinal anesthesia needles

31 pd- Spinalinè adata Quincke tipo 22 G 120 mm



Atrau-com®

- Multi-angle grinded tip with special bending
- Extra thin wall
- Atraumatic

pencil point


- Perfect conic elliptical bevel enables easy access to desired location, without cutting or harming tissues.


quincke tip


- Quincke bevel is formed using state of the art technology. Smooth needle point ensures easy insertion into tissues.
- Perfect matching of stylet and needle bevel prevent coring. Angular tip prevents tissue damage.
- Crystal clear needle hub enables monitoring fluid flow and helps to keep CSF loss at a minimum level.
- Thin wall needle lets faster fluid flow.

Tulip Design

kodas Aprašymas dydis

SIZES	reference	product name	gauge	color	length mm	length inch
<p>quincke</p>  <p>Quincke</p>	TSPQ1875	18 G – 75 mm Spinal Needle Quincke Bevel	18G		75mm	3"
	TSPQ1890	18 G – 90 mm Spinal Needle Quincke Bevel	18G		90mm	3½"
	TSPQ1990	19 G – 90 mm Spinal Needle Quincke Bevel	19G		90mm	3½"
	TSPQ2075	20 G – 75 mm Spinal Needle Quincke Bevel	20G		75mm	3"
	TSPQ2090	20 G – 90 mm Spinal Needle Quincke Bevel	20G		90mm	3½"
	TSPQ2190	21 G – 90 mm Spinal Needle Quincke Bevel	21G		90mm	3½"
	TSPQ2240	22 G – 40 mm Spinal Needle Quincke Bevel	22G		40mm	1½"
	TSPQ2275	22 G – 75 mm Spinal Needle Quincke Bevel	22G		75mm	3"
	TSPQ2290	22 G – 90 mm Spinal Needle Quincke Bevel	22G		90mm	3½"
	TSPQ22120	22 G – 120 mm Spinal Needle Quincke Bevel	22G		120mm	4¾"
	TSPQ2390	23 G – 90 mm Spinal Needle Quincke Bevel	23G		90mm	3½"
	TSPQ2490	24 G – 90 mm Spinal Needle Quincke Bevel	24G		90mm	3½"
	TSPQ2575	25 G – 75 mm Spinal Needle Quincke Bevel	25G		75mm	3"
	TSPQ2590	25 G – 90 mm Spinal Needle Quincke Bevel	25G		90mm	3½"
	TSPQK2590	25 G – 90 mm Spinal Needle Quincke Bevel with Guide	25G		90mm	3½"
	TSPQK25103	25 G – 103 mm Spinal Needle Quincke Bevel with Guide	25G		103mm	4"
	TSPQK25120	25 G – 120 mm Spinal Needle Quincke Bevel with Guide	25G		120mm	4¾"
	TSPQ2690	26 G – 90 mm Spinal Needle Quincke Bevel	26G		90mm	3½"
	TSPQK2690	26 G – 90 mm Spinal Needle Quincke Bevel with Guide	26G		90mm	3½"
	TSPQK26103	26 G – 103 mm Spinal Needle Quincke Bevel with Guide	26G		103mm	4"
	TSPQK26120	26 G – 120 mm Spinal Needle Quincke Bevel with Guide	26G		120mm	4¾"
	TSPQ2790	27 G – 90 mm Spinal Needle Quincke Bevel	27G		90mm	3½"
	TSPQK2790	27 G – 90 mm Spinal Needle Quincke Bevel with Guide	27G		90mm	3½"
	TSPQK27103	27 G – 103 mm Spinal Needle Quincke Bevel with Guide	27G		103mm	4"
	TSPQK27120	27 G – 120 mm Spinal Needle Quincke Bevel with Guide	27G		120mm	4¾"
	TSPQ2990	29 G – 90 mm Spinal Needle Quincke Bevel	29G		90mm	3½"
	TSPQK2990	29 G – 90 mm Spinal Needle Quincke Bevel with Guide	29G		90mm	3½"
	TSPQK29120	29 G – 120 mm Spinal Needle Quincke Bevel with Guide	29G		120mm	4¾"

SIZES	reference	product name	gauge	color	length mm	length inch
<p>pencil point</p> 	TSP2290	22 G – 90 mm Spinal Needle Pencil Point	22G		90mm	3½"
	TSPPK2290	22 G – 90 mm Spinal Needle Pencil Point with Guide	22G		90mm	3½"
	TSP2390	23 G – 90 mm Spinal Needle Pencil Point	23G		90mm	3½"
	TSP2490	24 G – 90 mm Spinal Needle Pencil Point	24G		90mm	3½"
	TSP2590	25 G – 90 mm Spinal Needle Pencil Point	25G		90mm	3½"
	TSPPK2590	25 G – 90 mm Spinal Needle Pencil Point with Guide	25G		90mm	3½"
	TSPPK25103	25 G – 103 mm Spinal Needle Pencil Point with Guide	25G		103mm	4"
	TSPPK25120	25 G – 120 mm Spinal Needle Pencil Point with Guide	25G		120mm	4¾"
	TSP2690	26 G – 90 mm Spinal Needle Pencil Point	26G		90mm	3½"
	TSPPK2690	26 G – 90 mm Spinal Needle Pencil Point with Guide	26G		90mm	3½"
	TSPPK26103	26 G – 103 mm Spinal Needle Pencil Point with Guide	26G		103mm	4"
	TSP2790	27 G – 90 mm Spinal Needle Pencil Point	27G		90mm	3½"
	TSPPK2790	27 G – 90 mm Spinal Needle Pencil Point with Guide	27G		90mm	3½"
	TSPPK27103	27 G – 103 mm Spinal Needle Pencil Point with Guide	27G		103mm	4"
	TSPPK27120	27 G – 120 mm Spinal Needle Pencil Point with Guide	27G		120mm	4¾"
	TSP2990	29 G – 90 mm Spinal Needle Pencil Point	29G		90mm	3½"
	TSPPK2990	29 G – 90 mm Spinal Needle Pencil Point with Guide	29G		90mm	3½"

SIZES	reference	product name	gauge	color	length mm	length inch
<p>atrau-com®</p> 	TSPAK2590	25 G – 90 mm Atraucom Spinal Needle with Guide	25G		90mm	3½"
	TSPAK25100	25 G – 100 mm Atraucom Spinal Needle with Guide	25G		100mm	4"
	TSPA2650	26 G – 50 mm Atraucom Spinal Needle	26G		50mm	2"
	TSPAK2688	26 G – 88 mm Atraucom Spinal Needle with Guide	26G		88mm	3½"
	TSPAK26100	26 G – 100 mm Atraucom Spinal Needle with Guide	26G		100mm	4"

Spinalinės anestezijos adatos

Atraucom tipo

- Daugiakampių galiuku, specialiai išlenkta
- Ypač plona sienele
- Atraumatinė

Pieštukinės

- Ideali konusinė eliptinė nuopjova leidžia lengvai pasiekti lokaciją be pjovimo ir audinių pažeidimų

Quincke tipo galiuku

- Quincke nuopjova suformuota specialia technologija
- Tolygus adatos galiukas užtikrina lengvą įvedimą
- Idealiai suderintas stiletas su adatos nuopjova. Kampuotas galiukas apsaugo audinius nuo pažeidimų
- Itin skaidri adatos jungtis leidžia stebėti skysčio tekėjimą ir padeda palaikyti jo netekimą minimaliame lygyje.
- Plona sienelė leidžia skysčiams tekėti greičiau.

BAKTERİ FİLTRELERİ Bacterial Filter

Steril	Sterile
70 mm	70 mm
BFE %99,997	BFE %99,997
(bakteriyal filtre etkinliği)	(bacterial filtration efficiency)
VFE % 99,9998	BFE %99,9998
(viral filtre etkinliği)	(viral filtration efficiency)



Tip/Type	Ref No
Bakteri/Bacterial	140 101
HME	140 120
HMEF	140 110

SPIROMETRE FİLTRELERİ

Spirometer Filter

Her türlü cihaza uygun
BFE %99,990
(bakteriyal filtre etkinliği)
VFE % 99,9979
(viral filtre etkinliği)

Suitable for all kind of
ventilation machines
BFE %99,990
(bacterial filtration efficiency)
VFE % 99,9979
(viral filtration efficiency)



Bağlantı Çapları	İç çap / Inn dia	Dış çap / Out dia
	21, 26, 30	24, 26, 30, 34

	Ref No
Tip/Type	Sterile
Kademeli / Graded	141 101
Elips / Eliptic	141 102
Dairesel / Circular	141 103

SPIROMETRE AĞIZLIĞI Spirometer Mouthpiece

65 mm boy	65 mm lenght
28, 30, 33 mm dış çap	28, 30, 33 mm external size

Tip/Type	Ref No
28 mm	141 228
30 mm	141 230
33 mm	141 233



BURUN MANDALI

Nose Clip



Nosies spaustukai



CODE / KOD
141 110



Tissue and Organ Retrieval Bags




42 poz
X SERIES

W SERIES

Q SERIES

*Safe, Fast and Easy to Use for laparoscopic
operations with three different models
and large dimension alternatives*



Endobags of can be closed off more than one in abdomen. They are produced two different design as structure of asymmetric mouth and circular mouth. Endobags are manufactured from polyuretan that streches 80 and 100 micron.

- ✓ **Extra wide and asymmetric mouth structure facilitates insertion of specimen**
Intuitional design facilitates the extraction of specimen
- ✓ **Extraordinary durability and closure**
Maximize the specimen's safety
Minimize the risk of contamination
Provides complete isolation for cancerous tissues
- ✓ **Antibacterial, Non-toxic and Non-allergic**
Provides safe us efor all patient groups
- ✓ **Flexible and Transparent**
Visible specimen



introducer with double handle



opening direction indicator (intra - abdominal)



asymmetric mouth design



circular mouth design



opening direction indicator (push handle)

REFERENCE	dydis SIZE	išmatavimai DIMENSIONS	IN A BOX	tūris VOLUME
DLR-XS05H	X-SMALL	50mm x 100mm	12 PCS	100ml
DLR-S10H	SMALL	65mm x 150mm	12 PCS	200ml
DLR-M10H	MEDIUM	75mm x 200mm	12 PCS	400ml
DLR-L10H	LARGE 10	100mm x 200mm	12 PCS	800ml
DLR-L12H	LARGE12	125mm x 230mm	12 PCS	1500ml
DLR-XL15H	X-LARGE	130mm x 250mm	12 PCS	2000ml
DLR-XXL15H	XX-LARGE	150mm x 250mm	12 PCS	2500ml

42 poz

Sterilus, supakuotas po 1 vnt.

tiemann catheter 43 P.D.

sterile

features:

- made of soft and flexible PVC, resistant to bending and twisting
- atraumatic, slightly rounded special tip facilitating insertion of the catheter
- tip with two side holes with soft rounded edges
- cone-shaped curved tip, bent at an angle of 45°
- transparent drain enabling visual inspection
- colour-coded connectors for size recognition, numerical size identification on the packaging
- single-use
- DEHP-free
- EO sterilized
- individual packaging: paper-foil

indication:

- catheterization of the urinary bladder in the case of a narrowed urethra

ref.	length	size	colour	intermediate packaging	bulk packaging
CT-08-40	400 mm	CH 8	blue	60 pcs	10 x 60 pcs
CT-10-40	400 mm	CH 10	black	60 pcs	10 x 60 pcs
CT-12-40	400 mm	CH 12	white	50 pcs	10 x 50 pcs
CT-14-40	400 mm	CH 14	green	50 pcs	10 x 50 pcs
CT-16-40	400 mm	CH 16	orange	50 pcs	10 x 50 pcs
CT-18-40	400 mm	CH 18	red	40 pcs	10 x 40 pcs
CT-20-40	400 mm	CH 20	yellow	40 pcs	10 x 40 pcs
CT-22-40	400 mm	CH 22	purple	40 pcs	10 x 40 pcs
CT-24-40	400 mm	CH 24	dark blue	40 pcs	10 x 40 pcs



urine bag hanger

non-sterile

features:

- made of strong and durable plastic
- special clamp prevents tube bending
- suitable for round and square bed frames
- DEHP-free

indication:

- hanging standard urine collection bags

ref.	packaging unit	bulk packaging
WNM01	50 pcs	4 x 50 pcs



urine bag

sterile

features:

- latex-free bag
- capacity 2000 ml
- made of 0.22 mm thick, latex-free PVC (double wall)
- double sealed PVC to prevent spilling of fluids
- effective inlet valve
- tight T-valve, can be easily operated with one hand
- connecting drain with universal tip, standard length 90 cm (other lengths available: 110, 130, 150 cm)
- drain made of material preventing its bending or twisting and allowing easy flow of urine
- precise, easy-to-read scale (every 100 ml)
- strengthened holes for hanging, matching standard hangers
- white back wall making urine easy to notice
- single-use
- EO sterilized
- individual packaging: foil

indication:

- 24-hour collection of urine and other body fluids

ref.	capacity	length of the drain	intermediate packaging	bulk packaging
WMD2000-90	2000 ml	90 cm	10 pcs	25 x 10 pcs
WMD2000-150	2000 ml	150 cm	10 pcs	25 x 10 pcs





Tiemann Kateteriai

- Pagaminti iš minkšto lankstaus PVC, atsparaus lankstymui ir sukimui
- Atraumatinis, šiek tiek užapvalintas specialus galiukas leidžia įvesti kateterį
- Galiukas su dvejomis šoninėmis angomis švelniais užapvalintais kraštais
- Kūgio formos užlenktas galiukas, lenktas 45° kampu
- Permatomas – suteikia galimybę vizualiniam patikrinimui
- Spalva koduoti antgaliai – dydžio atpažinimui, dydžių identifikacija skaičiais ant pakuotės
- Vienkartinio naudojimo
- Be DEHP
- EO sterilizuota
- Individuali pakuotė

Indikacijos:

- Šlapimo pūslės kateterizavimui susiaurėjusios šlaplės atveju

Ref. Kodas	Ilgis	Dydis	Spalva	Pakuotė	Pakuotė
CT-08-40	400 mm	CH 8	Mėlyna	60 vnt	10 x 60 vnt
CT-10-40	400 mm	CH 10	Juoda	60 vnt	10 x 60 vnt
CT-12-40	400 mm	CH 12	Balta	50 vnt	10 x 50 vnt
CT-14-40	400 mm	CH 14	Žalia	50 vnt	10 x 50 vnt
CT-16-40	400 mm	CH 16	Oranžinė	50 vnt	10 x 50 vnt
CT-18-40	400 mm	CH 18	Raudona	40 vnt	10 x 40 vnt
CT-20-40	400 mm	CH 20	Geltona	40 vnt	10 x 40 vnt
CT-22-40	400 mm	CH 22	Violetinė	40 vnt	10 x 40 vnt
CT-24-40	400 mm	CH 24	Tamsiai Mėlyna	40 vnt	10 x 40 vnt

Zabrze, 29.09.2020

Statement

We, Zarys International Group Sp. z o.o. sp.k. hereby declare that our product, Tiemann catheter, is non-pyrogenic.

PRODUCT MANAGER
ZARYS International Group sp. z o.o. sp.k.
Bożena Smolnik

Zabrze, 29.09.2020

PATVIRTINIMAS

Mes, Zarys International Group Sp. z o.o.sp.k. pareiškiame, kad produktas Tiemann kateteris yra nepirogeniškas.

Produkto vadybininkas
Zarys International Group Sp. z o.o.sp.k.
Bożena Smolnik

URETERAL STENTS & SETS


ÜRETERAL STENT ve SETLERİ



Trusted, high performance design. Proven TPU material. Firm during placement for easier stent passage, then softenes at body temperature for increased patient comfort. Used to provide temporary drainage from the ureteropelvic junction to the bladder. Supplied sterile in peel-open packages. Intended for single use. **Caution:** These stents are not intended as permanent indwelling devices. Periodic evaluation via cystoscopic or radio-graphic means is suggested. Must not remain indwelling for longer than 3 months.

Güvenilir, yüksek performans dizayn ile TPU konstrüksiyonu uygulamada kolaylık ve malzemenin vücut içinde yumuşamasına olanak tanımaktadır. UPJ ile mesane arasında geçici drenaj sağlamak amacı ile kullanılmaktadır. Steril paketlerde sunulmaktadır. Tek kullanımlık olarak dizayn ve imal edilmiştir.

Dikkat: Bu stentler kalıcı tıbbi cihaz değildir ve vücutta kalıcı kullanıma uygun değildir. Sistoskop veya radio-grafik kontroller önerilmektedir. Yerleştirilmiş olarak yerinde 3 aydan fazla kalmamalıdır.



Poz. 44

Ø FR	Guide Wire Kilavuz Tel Ø inch	Length Uzunluk cm	Distal End Closed Bir Ucu Kapalı 1 5	Both Ends Open İki Ucu Açık 2 5	Distal End Closed SET Bir Ucu Kapalı SET 1 3 4 5	Both Ends Open SET İki Ucu Açık SET 2 3 4 5
3	0,018	8	500001	500085	510001	510085
	0,018	10	500002	500086	510002	510086
	0,018	12	500003	500087	510003	510087
	0,018	14	500004	500088	510004	510088
	0,018	16	500005	500089	510005	510089
	0,018	18	500006	500090	510006	510090
	0,018	20	500007	500091	510007	510091
	0,018	22	500008	500092	510008	510092
	0,018	24	500009	500093	510009	510093
	0,018	26	500010	500094	510010	510094
4	0,028	10	500014	500098	510014	510098
	0,028	12	500015	500099	510015	510099
	0,028	14	500016	500100	510016	510100
	0,028	16	500017	500101	510017	510101
	0,028	18	500018	500102	510018	510102
	0,028	20	500019	500103	510019	510103
	0,028	22	500020	500104	510020	510104
	0,028	24	500021	500105	510021	510105
	0,028	26	500022	500106	510022	510106
	0,028	28	500023	500107	510023	510107
	0,028	30	500024	500108	510024	510108

Ø FR	Guide Wire Kılavuz Tel Ø inch	Lenght Uzunluk cm	Distal End Closed Bir Ucu Kapalı 1 5	Both Ends Open İki Ucu Açık 2 5	Distal End Closed SET Bir Ucu Kapalı SET 1 3 4 5	Both Ends Open SET İki Ucu Açık SET 2 3 4 5
4,7	0,032	8	500025	500109	510025	510109
	0,032	10	500026	500110	510026	510110
	0,032	12	500027	500111	510027	510111
	0,032	14	500028	500112	510028	510112
	0,032	16	500029	500113	510029	510113
	0,032	18	500030	500114	510030	510114
	0,032	20	500031	500115	510031	510115
	0,032	22	500032	500116	510032	510116
	0,032	24	500033	500117	510033	510117
	0,032	26	500034	500118	510034	510118
	0,032	28	500035	500119	510035	510119
	0,032	30	500036	500120	510036	510120
4,8	0,035	8	500037	500121	510037	510121
	0,035	10	500038	500122	510038	510122
	0,035	12	500039	500123	510039	510123
	0,035	14	500040	500124	510040	510124
	0,035	16	500041	500125	510041	510125
	0,035	18	500042	500126	510042	510126
	0,035	20	500043	500127	510043	510127
	0,035	22	500044	500128	510044	510128
	0,035	24	500045	500129	510045	510129
	0,035	26	500046	500130	510046	510130
	0,035	28	500047	500131	510047	510131
	0,035	30	500048	500132	510048	510132
5	0,035	8	500049	500133	510049	510133
	0,035	10	500050	500134	510050	510134
	0,035	12	500051	500135	510051	510135
	0,035	14	500052	500136	510052	510136
	0,035	16	500053	500137	510053	510137
	0,035	18	500054	500138	510054	510138
	0,035	20	500055	500139	510055	510139
	0,035	22	500056	500140	510056	510140
	0,035	24	500057	500141	510057	510141
	0,035	26	500058	500142	510058	510142
	0,035	28	500059	500143	510059	510143
	0,035	30	500060	500144	510060	510144
6	0,038	8	500061	500145	510061	510145
	0,038	10	500062	500146	510062	510146
	0,038	12	500063	500147	510063	510147
	0,038	14	500064	500148	510064	510148
	0,038	16	500065	500149	510065	510149
	0,038	18	500066	500150	510066	510150
	0,038	20	500067	500151	510067	510151
	0,038	22	500068	500152	510068	510152
	0,038	24	500069	500153	510069	510153
	0,038	26	500070	500154	510070	510154
	0,038	28	500071	500155	510071	510155
	0,038	30	500072	500156	510072	510156
7	0,038	24	500081	500165	510085	510065
	0,038	26	500082	500166	510082	510166
	0,038	28	500083	500167	510083	510167
	0,038	30	500084	500168	510084	510168

HIGH PRESSURE SYRINGE

63 pd- Accutron MR injektoriaus vienkartinų švirkštų rinkinys



P/N: 400305

Compatible with Medtron Accutron MRI Injector

Suderinama su Medron Accutron MRI injektoriumi

Contents: 1-65ml Syringe

1vnt - 65ml švirkštas

1-200ml Syringe

1vnt- 200ml švirkštas

1-MRI Connecting Tube

1vnt- MRI prijungimo vamzdelis

1-Long Spike

ilgas smaigas

1-Short Spike

trumpas smaigas

Packing: 20PCS/CASE

Pressure: 350PSI

Slėgis: 350PSI



antmed
Partner For Life



- Bar-bell shaped pyrolytic carbon coated ceramic with exceptionally biocompatibility with clear imaging characteristics
- Unique shapes designed to prevent migration **negalintis migruotis**
- No allergy response or hypersensitivity reported **nealergizuojantis**
- Designed for safe, reliable, and precise placement during core biopsies
- Permanent and clear visibility with mammography, ultrasound and MRI **naudojamas mamografijoje, su UG ir MRI**
- No interference with MRI visualization, lower level of artifact than metal markers
- No interference with Tomosynthesis visualization
- Available in three distinct size configurations
- Available for vacuum assisted and core biopsy needle devices

WHEN TO USE

- Marking after removal of tumor
- Marking microcalcification area after breast vacuum assisted biopsy (VAB)
- Marking of suspicious breast lesions after critical core biopsy
- Orientation for planning radiation therapy



Ultrasound



Mammography

BIOMARC Qty/Box 10pcs			
CODE	SIZE	MARKER	DESCRIPTION
040220	8Gx165,3mm	2x4mm	Flexible Applicator, 8G Mammotome®
040203	9Gx165,3mm	2x4mm	Flexible Applicator, 11G Mammotome®
040204	9 G x 1 2 0 m m	2x4mm	Flexible Applicator, 9G Eviva® / 10G Spirotome®
040228	9Gx363,7mm	2x4mm	Flexible Applicator, 9G Atec®
040239	9Gx363,7mm	2x4mm	Flexible Applicator, 10G SenoRX® / 10G Bexcore®
040306	13Gx117,7mm	2x4mm	Rigid Applicator, 13G
040320	14Gx117,7mm	1x4mm	Rigid Applicator, 14G
040300	14Gx117,7mm	1x3mm	Rigid Applicator, 14G
040310	17Gx117,7mm	0,9x5mm	Rigid Applicator, 17G

Kodas

Dydis

Markeris

On request, Vigeo is available to provide needles with different diameter and length

PAJUNK®

ChibaSono & Chiba
Fine Needle Biopsy Products

69 pd- Transabdominalinēs adatos amniocentezēms



MADE IN GERMANY

ChibaSono


The echogenic biopsy needle

Echogeniška adata

Ultrasound guided biopsies have been increasing in the recent years. Together with Dr. Chris Mitchell, Pajunk has developed the ChibaSono, a new generation of needles that allow for a precise localisation under ultrasound. Needle and tip design allow the user to focus on the procedure and patient by making the tip visualisation under ultrasound faster, easier and safer.



Approved for
prenatal diagnostics



ChibaSono with beveled tip

Stylet ground exactly flush

- ➔ Increased needle stability
- ➔ Avoids tissue carry-over into the puncture area



ChibaSono with US tip

Special needle cut and round blunt stylet, set back from tip


- ➔ Double reflex for improved needle tip visibility under ultrasound




The essentials:

- ➔ Optimum sonographic visibility
- ➔ High precision needle tip
- ➔ Excellent gliding properties
- ➔ Depth control ensures patient and user safety

Depth markings
(selected products only)



OptiView = Optimum echogenicity from needle shaft to tip



Due to their geometry, arrangement and placement, the three-dimensional Cornerstone impressions ensure outstanding sonographic visibility.¹


- ➔ The ultrasonic waves reflect without obstruction and are easier to identify irrespective of insertion angle²
- ➔ Outstanding needle shaft and tip visibility¹
- ➔ 360° clear identification over a total length of 20 mm^{1,2}
- ➔ High echogenicity at steep insertion angle³

¹ Edgcombe, Hocking, Sonographic identification of needle tip ..., 2010; 35(2): 207–211

² Hebard, Hocking, Echogenic technology can improve needle visibility ..., 2011; 36(2): 185–189


³ Uppal, Sondekoppam, Ganapathy, Effect of beam steering on ..., 2014; 61(10): 909–915

ChibaSono with beveled tip




ChibaSono needle with **beveled tip**,
Cornerstone Reflectors, depth stop,
graduation, retaining clip and wings

Size	Item no.	PU
18 G x 90 mm	532S120090	25
18 G x 150 mm	532S120150	10
20 G x 90 mm	532S095090	25
20 G x 150 mm	532S095150	10



ChibaSono with US tip



ChibaSono needle with **US tip**,
Cornerstone Reflectors, graduation
and wings

Dydis

Kodas

Size	Item no.	PU
22 G x 50 mm	530S070050	25
22 G x 90 mm	530S070090	25
22 G x 100 mm	530S070100	25
22 G x 125 mm	530S070125	10
22 G x 150 mm	530S070150	10
22 G x 220 mm	530S070220	10
20 G x 90 mm	530S095090	25
20 G x 120 mm	530S095120	10
20 G x 150 mm	530S095150	10
20 G x 220 mm	530S095220	10
18 G x 150 mm	530S120150	10






SWAB, SUCTION SWAB & TOOTBRUSH FOR ORAL CARE

ORAL CARE PROGRAM: Coremec disposable accessories for a appropriate hygiene procedure for mouth cavity.

Coremec has a complete line of swabs which will facilitate a proper appliance of oral hygiene and is focused on reducing hospital acquired infections for patients. The swabs and Brushes for oral hygiene "oral care" are packaged individually and are available in neutral version or treated with a refreshing mint substance. All items have a long plastic handle in order to prevent the possibility of contamination by contact. The foam tip has a curved shape suitable to easily remove the residue adhering to the teeth, helping to prevent dental plaque formation and the risk of VAP (ventilator associated pneumonia). Our Coremec Sponge Swabs are also available with the suction channels with a valve for regulating flow.



76 poz.

Art. A9000 Disposable oral SWAB Untreated (rose)
Individually wrapped - non sterile Medical device class I



250 un./Box

Art. A9100 Disposable oral SWAB mint treated (green)
Individually wrapped - non sterile Medical device class I



250 un. /Box

Vienkartinis šepetėlis su siurbimo galimybe
Art. C9010 Disposable oral SUCTION SWAB Untreated (green)
Individually wrapped - non sterile Medical device class IIa
Individuali pakuotė



50 un. /Box

Art. C9110 Disposable oral SUCTION SWAB mint treated (green)
Individually wrapped - non sterile Medical device class IIa



50 un. /Box

Art. D9011 Disposable oral SUCTION TOOTHBRUSH Untreated gr.
Individ. wrapped - non sterile Medical device class IIa



50 un. /Box

Art. D9111 Disposable oral SUCTION TOOTHBRUSH mint treated (gr.)
Individually wrapped - non sterile Medical device class IIa




50 un. /Box

Ingredients treated sponge:
Agua, sodium lauryl sulfate, mint flavour, ethyl alcohol, sodium saccharine.
CEE 93/42 - 2007/47

LATEX FREE

Used for stone manipulation and removal in the urinary tract. The nitinol memory wire allows the basket to retain its shape following extreme torsion. The unique tipless design aids in calyceal stone capture by allowing the basket to open at the base of the calyx or in the ureter and encircle the stone. The tipless

Poz. 82



LUMINA® NITINOL STONE BASKETS		
CODE	DESCRIPTION	QTY/ BOX
7743090	LUMINA STONE BASKET 3.0 FR. 4 WIRE 90 CM (H)	5/BOX
7744090	LUMINA STONE BASKET 4.0 FR. 4 WIRE 90 CM (H)	5/BOX
7753090	LUMINA STONE BASKET 3.0 FR. 5 WIRE 90 CM (H)	5/BOX
7754090	LUMINA STONE BASKET 4.0 FR. 5 WIRE 90 CM (H)	5/BOX
7743011	LUMINA STONE BASKET 3.0 FR. 4 WIRE 110 CM (H)3 Fr 4 vielučių 110 cm	5/BOX
7744011	LUMINA STONE BASKET 4.0 FR. 4 WIRE 110 CM (H)helical pynimo	5/BOX
7753011	LUMINA STONE BASKET 3.0 FR. 5 WIRE 110 CM (H)	5/BOX
7754011	LUMINA STONE BASKET 4.0 FR. 5 WIRE 110 CM (H)	5/BOX
7743013	LUMINA STONE BASKET 3.0 FR. 4 WIRE 130 CM (H)	5/BOX
7744013	LUMINA STONE BASKET 4.0 FR. 4 WIRE 130 CM (H)	5/BOX
7753013	LUMINA STONE BASKET 3.0 FR. 5 WIRE 130 CM (H)	5/BOX
7754013	LUMINA STONE BASKET 4.0 FR. 5 WIRE 130 CM (H)	5/BOX
7843090	LUMINA STONE BASKET 3.0 FR. 4 WIRE 90 CM (S)	5/BOX
7844090	LUMINA STONE BASKET 4.0 FR. 4 WIRE 90 CM (S)	5/BOX
7853090	LUMINA STONE BASKET 3.0 FR. 5 WIRE 90 CM (S)	5/BOX
7854090	LUMINA STONE BASKET 4.0 FR. 5 WIRE 90 CM (S)	5/BOX
7843011	LUMINA STONE BASKET 3.0 FR. 4 WIRE 110 CM (S)	5/BOX
7844011	LUMINA STONE BASKET 4.0 FR. 4 WIRE 110 CM (S)	5/BOX
7853011	LUMINA STONE BASKET 3.0 FR. 5 WIRE 110 CM (S)	5/BOX
7854011	LUMINA STONE BASKET 4.0 FR. 5 WIRE 110 CM (S)	5/BOX
7843013	LUMINA STONE BASKET 3.0 FR. 4 WIRE 130 CM (S)	5/BOX
7844013	LUMINA STONE BASKET 4.0 FR. 4 WIRE 130 CM (S)	5/BOX
7853013	LUMINA STONE BASKET 3.0 FR. 5 WIRE 130 CM (S)	5/BOX
7854013	LUMINA STONE BASKET 4.0 FR. 5 WIRE 130 CM (S)	5/BOX
7943090	LUMINA STONE BASKET 3.0 FR. 4 WIRE 90 CM (N)	5/BOX
7944090	LUMINA STONE BASKET 4.0 FR. 4 WIRE 90 CM (N)	5/BOX
7953090	LUMINA STONE BASKET 3.0 FR. 5 WIRE 90 CM (N)	5/BOX
7954090	LUMINA STONE BASKET 4.0 FR. 5 WIRE 90 CM (N)	5/BOX
7943011	LUMINA STONE BASKET 3.0 FR. 4 WIRE 110 CM (N)	5/BOX
7944011	LUMINA STONE BASKET 4.0 FR. 4 WIRE 110 CM (N)	5/BOX
7953011	LUMINA STONE BASKET 3.0 FR. 5 WIRE 110 CM (N)	5/BOX
7954011	LUMINA STONE BASKET 4.0 FR. 5 WIRE 110 CM (N)	5/BOX
7943013	LUMINA STONE BASKET 3.0 FR. 4 WIRE 130 CM (N)	5/BOX
7944013	LUMINA STONE BASKET 4.0 FR. 4 WIRE 130 CM (N)	5/BOX
7953013	LUMINA STONE BASKET 3.0 FR. 5 WIRE 130 CM (N)	5/BOX
7954013	LUMINA STONE BASKET 4.0 FR. 5 WIRE 130 CM (N)	5/BOX
7043090	LUMINA STONE BASKET 3.0 FR. 4 WIRE 90 CM (N-ZERO TIP)	5/BOX
7053090	LUMINA STONE BASKET 3.0 FR. 5 WIRE 90 CM (N-ZERO TIP)	5/BOX
7043011	LUMINA STONE BASKET 3.0 FR. 4 WIRE 110 CM (N-ZERO TIP)	5/BOX
7053011	LUMINA STONE BASKET 3.0 FR. 5 WIRE 110 CM (N-ZERO TIP)	5/BOX
7043013	LUMINA STONE BASKET 3.0 FR. 4 WIRE 130 CM (N-ZERO TIP)	5/BOX
7053013	LUMINA STONE BASKET 3.0 FR. 5 WIRE 130 CM (N-ZERO TIP)	5/BOX
7133090	LUMINA STONE GRASPER 3.0 FR. 90 CM (TRIPRONG)	5/BOX
7133012	LUMINA STONE GRASPER 3.0 FR. 120 CM (TRIPRONG)	5/BOX
7131238	LUMINA PCN BASKET (PERK) 12.0 FR 38CM	5/BOX
7131037	LUMINA PERK BASKET(PCNL) 10 FR 37 CM	5/BOX
7131038	LUMINA PERK BASKET(PCNL) 10 FR 38 CM	5/BOX
7041512	LUMINA STONE BASKET 1.5FR 4WIRE 120CM NITINOL-(Z)	5/BOX
7041712	LUMINA STONE BASKET 1.7FR 4WIRE 120CM NITINOL-(Z)	5/BOX
7041912	LUMINA STONE BASKET 1.9FR 4WIRE 120CM NITINOL-(Z)	5/BOX
7042212	LUMINA STONE BASKET 2.2FR 4WIRE 120CM NITINOL-(Z)	5/BOX
7042412	LUMINA STONE BASKET 2.4FR 4WIRE 120CM NITINOL-(Z)	5/BOX



Dokumentą elektroniniu parašu pasirašė DOVILĖ ANDRIJAUSKAITĖ
Data: 2020-10-02 13:55:41
Paskirtis: Pirkimo Nr. 504288
Vieta: Energetikų g. 8, LT-52461, Kaunas
Kontaktinė informacija: UAB "Skirgesa", Viešųjų pirkimų specialistė