

SelectBag® Citrate

SelectBag® Citrate is a unique acetate-free concentrate with all the beneficial properties of citrate. With the BiCart Select Citrate system it is easy to optimize the dialysis fluid composition. Sodium and bicarbonate levels are accurately controlled to your settings while the SelectBag Citrate container allows you to easily vary potassium, calcium and glucose.

- Improved acid-base status
- Easier electrolyte individualization
- Less thrombogenic dialysis
- Reduced inflammation
- Improved tolerance to treatment
- Improved efficiency



BiCart Select® Citrate System

SelectBag® Citrate container



BiCart Select® combi-pak (1 BiCart® cartridge + 1 SelectCart® cartridge)



Product description:	Electrolytes (1:200) for preparation of dialysis fluid	Sodium bicarbonate and sodium chloride cartridges for preparation of dialysis fluid
Product code:	See table below	EU 107617 / WW 109574
Units per box:	6	4
Specifications		
Contents:	1000 ml concentrated solution	BiCart cartridge: 720 g NaHCO ₃ SelectCart cartridge: 1200g NaCl
pH of solution:	~ 1.5	
Dilution ratio:	1:200	
Capacity:	200 liters of dialysis fluid can be produced	At Na ⁺ /HCO ₃ ⁻ settings of 140/34, 190 liters of dialysis fluid can be produced
Duration:	6 hrs 40 min Based on a Q _D of 500 ml/min	
Storage temperature:	+4°C to +30°C	Store below +40°C
Shelf life:	1 year from date of manufacture. Refer to the expiry date stated on the product label	2 years from date of manufacture. Refer to the expiry date stated on the product label
Packaging:	The box is designed with a top opening for easy access	Both cartridges are packed together in a low-density polyethylene (LDPE) bag with a convenient handle
Cartridge material:		Polypropylene (PP)
Container material:	Low density ethylene (LDPE)	
Box dimensions:	Approx. 318 x 256 x 178 mm	Approx. 298 x 199 x 332 mm
Box weight:	Approx. 7 kg	Approx. 9 kg

SelectBag® Citrate			Contribution to the dialysis fluid				
Product code	Product name	Dilution ratio	K ⁺ mmol/l	Ca ²⁺ mmol/l	Mg ²⁺ mmol/l	Citrate mmol/l	Gluc. g/l
114683	CX150G	1:200	1.00	1.50	0.50	1	1
114845	CX165G	1:200	1.00	1.65	0.50	1	1
114684	CX175G	1:200	1.00	1.75	0.50	1	1
114685	CX250G	1:200	2.00	1.50	0.50	1	1
114687	CX260G	1:200	2.00	1.60	0.50	1	1
114846	CX265G	1:200	2.00	1.65	0.50	1	1
114688	CX275G	1:200	2.00	1.75	0.50	1	1
114686	CX2550G	1:200	2.50	1.50	0.50	1	1
114694	CX2565G	1:200	2.50	1.65	0.50	1	1
114849	CX2575G	1:200	2.50	1.75	0.50	1	1
114689	CX350G	1:200	3.00	1.50	0.50	1	1
114693	CX365G	1:200	3.00	1.65	0.50	1	1
114690	CX375G	1:200	3.00	1.75	0.50	1	1
114691	CX450G	1:200	4.00	1.50	0.50	1	1
114848	CX465G	1:200	4.00	1.65	0.50	1	1
114692	CX475G	1:200	4.00	1.75	0.50	1	1

SelectBag® Citrate			Contribution to the dialysis fluid				
Product code	Product name	Dilution ratio	K ⁺ mmol/l	Ca ²⁺ mmol/l	Mg ²⁺ mmol/l	Citrate mmol/l	Gluc. g/l
114695	CX150	1:200	1.00	1.50	0.50	1	0
114850	CX165	1:200	1.00	1.65	0.50	1	0
114696	CX250	1:200	2.00	1.50	0.50	1	0
114851	CX265	1:200	2.00	1.65	0.50	1	0
114697	CX275	1:200	2.00	1.75	0.50	1	0
114698	CX350	1:200	3.00	1.50	0.50	1	0
114852	CX365	1:200	3.00	1.65	0.50	1	0
114699	CX375	1:200	3.00	1.75	0.50	1	0

These specifications are subject to change without notice. For further information and operating instructions please refer to the operator manual

The BiCart Select system is compatible with Gambro dialysis systems.

CE 0086 The products meet the applicable provisions of Annex I (Essential Requirements) and Annex II (Full quality assurance system of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.

Gambro, BiCart Select, SelectBag, BiCart, SelectCart are trademarks belonging to the Gambro Group.

Gambro Lundia AB
PO Box 10101, SE-220 10 Lund, Sweden
Phone + 46 46 16 90 00, partner@gambro.com, www.gambro.com

