

Certofix[®] Central Venous Catheter Sets

Centrinės venos kateterizavimo rinkiniai su V tipo adata ir Nitinolo styga (atspari persilenkimui)

Uncompromising Quality for Safe Infusion Therapy




Central Venous Catheters

Certofix® Paed

... for children's safety

Customized for the little ones

B. Braun's long experience and global expertise in infusion therapy are incorporated in its paediatric range of central venous catheters. Uncompromising quality and state-of-the-art features combine to meet the special requirements of infants and children.




User benefits

- Safe and effective positioning
 - kink-proof guide wire
 - real-time ECG-lead position verification
- Gentle for sensitive patients
 - Teddy Bear catheter fixation makes stitching unnecessary
 - kink-proof guide wire avoids any damage of the intima


Product specifications

- Radio-opaque polyurethane catheter with
 - soft tip
 - length markings
 - junction hub for stitch fixation if required
 - transparent extension lines
 - color-coded hub
 - white slide clamps
- Main components
 - Certofix® central venous catheter with side clamps
 - kink-proof guide wire with soft "J" tip and straight tip on dispenser
 - Seldinger introducer needle
 - ECG lead cable
 - dilator
 - scalpel
 - soft, self-adhesive catheter fixation
- Accessories depending on set configuration
 - Safsite® valve or Discofix® three-way stopcock
 - 3 or 5 ml Omnifix® luer lock syringe
 - infusion extension line with side clamp (mono sets only)
 - Safsite® valve
 - movable fixation wing with clamp




Kink-proof guide wire

- Highly flexible and resilient guide wire made of Nitinol
- Guide wire with the just the right stiffness for easy advancing



Teddy Bear fixation

- Safe and easy catheter fixation without stitching
- Soft on the skin of infants and children
- Gentle removal possible



Introducer guide for guide wire

- Improves threading of the "J" tip of the guide wire into the Seldinger needle
- Allows dispenser-free advancing of the guide wire through the introducer

REF

Spindis, G

llgis, cm

Product Type	Code No. (REF)	Cath. Lumen ø G	Flow rate (ml/min) *	Length (cm)	Guide wire ø / Length (cm)	Valve / Plug
Mono Paed S 110*	4160177	22	D 15	10	0,5/25	–
Duo Paed S 408	4166906	22/22	D 18; P 18	8	0,46/50	Safsite® Valve
Duo Paed S 413	4166922	22/22	D 10; P 10	13	0,46/50	Safsite® Valve
Duo Paed S 420	p.d.17 4166949	22/22	D 9; P 9	20	0,46/50	Safsite® Valve
Duo Paed S 508	4167112	18/20	D 30; P 15	8	0,46/50	Safsite® Valve
Duo Paed S 513	4167139	18/20	D 22; P 11	13	0,46/50	Safsite® Valve
Duo Paed S 520	p.d.18 4167155	18/20	D 18; P 9	20	0,46/50	Safsite® Valve
Trio Paed S 508	4167228	20/22/22	D 30; M1 15; P 15	8	0,46/50	Safsite® Valve
Trio Paed S 513	4167244	20/22/22	D 29; M1 10; P 10	13	0,46/50	Safsite® Valve
Trio Paed S 520	4167260	20/22/22	D 18; M1 6; P 6	20	0,46/50	Safsite® Valve

* D (distal); M1 (middle1); M2 (middle2); M3 (middle3); P (proximal)

Additional information:

- The sales unit of Certofix® Paed Sets is 10 pieces
- All catheters are made of PUR