

## > oxygen therapy

## > flowmeters with pre-adjusted flows



065

**DEBSON TM O<sub>2</sub> flowmeter** with a CCO humidifier (250 cc) and a direct probe.



1.5.3.

066

**DEBSON TM Medical Air flowmeter** with an outlet tubing nipple and a complete 1.5m rail mounting system.



067

**Direct probes** for O<sub>2</sub> (left) and Medical Air (right).



068

**Outlet tubing nipples**, 1/4 G (left) and 9/16" (middle). 1/4 G to 9/16" **intermediate adaptor** for connecting a 1/4G outlet flowmeter to a 9/16" inlet humidifier or nebulizer (right).



069

**Complete rail mounting systems** for O<sub>2</sub> and Medical Air including:

- a crimped angled remote probe,
- 1.5m of antistatic hose (other lengths available upon request),
- a metal sliding piece,
- a EURO TM universal rail adaptor.