

# HANDLING PIPETTE

## PIPETTE FOR HANDLING EMBRYOS AND OOCYTES

The handling pipette is a single-use IVF laboratory micropipette used for transferring oocytes and embryos from one dish to another.

### FEATURES


- IVF laboratory micropipette for dish-to-dish transfer of oocytes and embryos
- will not scratch dishes
- available in three lengths: 45, 55 and 94 mm
- available in two packaging variations: single packed and packed by ten pieces



ORIGINAL PICTURE

### GYNO INFO HANDLING PIPETTE

Length	45, 55 or 94 mm ( $\pm 3$ mm)
Tube OD	1.25 mm ( $\pm 0.05$ mm)
Tube ID	1.00 mm ( $\pm 0.05$ mm)
Sterile	sterilized by gamma irradiation
Packaging	250/shipping box (1/pouch) 1250/shipping box (10/pouch)
Quality control	MEA and LAL
Order number	#3145 (45 mm, 1/pouch) #3100 (55 mm, 1/pouch) #3125 (94 mm, 1/pouch) #3145 TP (45 mm, 10/pouch) #3100 TP (55 mm, 10/pouch) #3125 TP (94 mm, 10/pouch)



GMP Belgium Handling Pipette 09/2016

