

# Eppendorf Safe-Lock Tubes



Catalog No. 0030121597

Inquire

Product Information

- 2,0 mL
- Biopur®, colorless, 100 tubes, individually packed

Show more Product Information



- ⚙ Request lot-specific certificates (not applicable for "Eppendorf Quality")
- ⚙ Learn more about purity grades available from Eppendorf

You will find additional download material at the bottom of this page

## Features

- Hinged Eppendorf Safe-Lock lid prevents unintentional lid opening during incubation and storage
- g-Safe®: Depending on Safe-Lock Tube variant (see manual) exceptional centrifugation stability up to 30,000 × g prevents sample loss due to tube breakage and provides extra safety when working with hazardous samples
- Manufactured without mold release agents, softeners or biocides – Materials which could have a negative effect on bioassays
- g-Safe®: Exceptional centrifugation stability up to 25,000 × g prevents sample loss due to tube breakage and provides extra safety when working with hazardous samples
- Precise lid sealing for minimal evaporation rates during long-term storage
- Large frosted lid and surface on the side for easy labeling
- Ensured functionality from -86 °C to 100 °C
- Autoclavable when open (121 °C, 20 min)
- Also available in Eppendorf Quality and in purity grades of lot-certified PCR clean and Forensic DNA Grade
- Also available in variants Eppendorf Protein LoBind or Eppendorf DNA LoBind for maximum recovery of valuable samples and in amber for light-sensitive samples

## Downloads: Eppendorf Safe-Lock Tubes

### Content

Brochures (12)	English ×
Certificates (6)	
Manuals (2)	Flyer Eppendorf Purity Grades Selection Guide
	PDF .3MB