


Zibo Eastmed Healthcare Products Co., Ltd

Add: No.118 Huaguang Road, Zhangdian District, Zibo City 255000, Shandong China

Infusion Sets with Plastic Spike

<ul style="list-style-type: none">- PVC tubing L = 150cm, ϕ 2.8*3.8mm- Vented drip chamber 48mm with medical filter & vented spike- Fluid filter 15μ- PE roller clamp with spike housing- 1ml = 20drops- Latex-free ball injection port- Male luer lock or luer slip with injection Needle 21G x 1 1/2''		<ul style="list-style-type: none">- Non-pyrogenic- For single use- ETO sterilized <div>CE 0197</div>	
	No.	Designation	Materials
	1	Plastic spike with cap	ABS + PE
	2	Air vent	ABS
	3	Drip chamber 48mm	PVC
	4	Filter	ABS, nylon
	5	PE roller clamp with spike housing	PE
	6	PVC tubing ϕ 2.8*3.8mm 180cm	PVC
	7	Latex-free ball injection port	Rubber
	8	Male luer lock or luer slip	ABS
	9	Injection needle 21G x 1 1/2''	PE + ABS + stainless steel

PACKING:

1. Individual Poly bag packed
2. Packing is adapted to the ethylene oxide sterilization. The cover mentions the necessary indications for content identification: designation, item code, manufacture date and expire date and lot number.
3. Packing details: 1pc/PE bag, 25pcs/middle PE bag, 500pcs/ctn, carton size is 61*38*32.5cm, gross weight is 13kgs, net weight is 12kgs.

QUALITY:

1. The product quality is guaranteed by inspections at all levels of manufacture, all our products are manufactured in 100,000 grade clean room, the final control is made at the level of packaging.
2. The packaged infusion set is sterilized by ethylene oxide according to a validated method in the physical and microbiological field. The shelf life is 5 years from the manufacture date.
3. The product conforms to MDD 93/42/EEC amended by 2007/47/EC and matins a quality assurance system according to ISO 13485.