

Aesculap Power Systems Skin Graft Mesher

The ultimate tool for skingraft expansion
and preparation



- Unique design eliminates need to change cutters mid-procedure
- Thin and thick graft cutter adjustments
- Skin carriers in three expansion ratios




AESCULAP®

Skin Graft Mesher


- Ideal addition to ACCULAN® Dermatome
- Single cutter cuts skin graft into geometrically defined pattern to enlarge the skin surface
- Graft expansion predetermined by skin carrier, eliminating the need to change cutter mid-procedure
- Continuous feed carrier prevents skin carrier from being loaded incorrectly and allows two carriers to be used without interruption for more uniformed grafts
- Unique cutter adjustment accommodates thin and thick grafts
- Easy cutting roller exchange
- Keyless disassembly




Skin Graft Mesher



Cutter Adjustment



Continuous Feed Carrier Guide



Keyless Disassembly

Skin Carriers

- Available in three expansion ratios for flexibility: 1:1.5, 1:3 & 1:6
- Packaged single, sterile, 10 per box

Ordering Information:

Item #	Description
BA720R	Accu-Mesh™ Skin Graft Mesher includes BA726 ratchet & BA725 cutter
BA721	Prepackaged Sterile Skin Graft Carrier expansion ratio 1:1.5, 10 per box
BA722	Prepackaged Sterile Skin Graft Carrier expansion ratio 1:3, 10 per box
BA723	Prepackaged Sterile Skin Graft Carrier expansion ratio 1:6, 10 per box

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.

B | BRAUN
SHARING EXPERTISE

Aesculap, Inc.

3773 Corporate Parkway
Center Valley, PA 18034
Phone: 800-282-9000

<http://www.aesculapusa.com>