


# AESCALAP® MECHANICAL PROCESSING

## ECCOS® HOLDER

Acculan 4

<p><b>GB485R</b></p> <p>ECCOS® HOLDER for GA875</p>  56 mm	<p><b>GB486R</b></p> <p>ECCOS® HOLDER for GA878</p>  56 mm	<p><b>GB688R 1.6</b></p> <p>ECCOS® HOLDER <b>Laikiklis</b> for BA720R</p>  106 mm
 	 	 
<p><b>GB481R</b></p> <p>ECCOS® HOLDER for saw blades and other accessories</p>  56 mm	<p><b>GB690R</b></p> <p>ECCOS® HOLDER 5 mm for drilling wire guide GA672207, GA344244</p>	<p><b>GB459R</b></p> <p>ECCOS® MOUNTING KIT</p>
 	 	
<p><b>JF148R</b></p> <p>MINI TRAY WITH LID</p> <p>outer size: 145 x 44 x 25 mm inner size: 140 x 39 x 18 mm</p> <p>for storage of Acculan 4 Dermatome rod and flaps during automatic cleaning process</p> 		