

Eleganza 4			
picture	description	code	EUR
Weighing scales			
	With scales and iBoard Basic	SC:1	
X ray modification			
	X ray design of the backrest; X ray cassette holder	XR:03	
Castors type			
	4 double castors Tente Integral light 150 mm one antistatic castor	CA:77	
	4 castors Tente Integral, 150 mm + 5th castor one antistatic castor	CA:73	
	4 double castors Tente Integral 150 mm + 5th castor one antistatic castor	CA:75	
	4 castors Tente Integral 150 mm, one antiscatic castor, i-Drive Power (motorized patient in-bed transport) <i>Compatibility with the servers an overbed tables shall be checked individually prior to order! One antistatic castor.</i>	CA:84	
	4 castors Tente Integral double 150 mm, one antiscatic castor, i-Drive Power (motorized patient in-bed transport) <i>Compatibility with the servers an overbed tables shall be checked individually prior to order! One antistatic castor.</i>	CA:85	
Brake system			
	Centrally operated castors, 2 brake levers, foot side operated + i brake	BR:27	
	Centrally operated castors e with 4 pedals - head sided operation + i-brake <i>Compatible with CA84 and CA85 only</i>	BR:35	
	Centrally operated castors with 4 pedals - head sided operation <i>Compatible with CA84 and CA85 only</i>	BR:36	
Foot control			
	foot controls (height adjustment, examination position)	FC:04	
Get up handles			
	Mobi Lift® handles; 1 pair with bed height adjustment buttons	MH:06	
Accessory holders			
	Pair of lateral accessory rails with removable plastic hooks	AH:01	
Nominal voltage			
	power supply 100 V	VO:02	
	power supply 110 V	VO:03	

punktas 1.29.3

Priedų bėgelių(pora) su nuimamais plastikiniais kabliukais