



RVTM2

Tube of 10 single use antibacterial filters

Single-use antibacterial plastic filter up-front

Non sterile filter

Product class I

Dimensions:

Diameter: 50 mm

Thickness: 0.47 mm

Performance:

Density: 100 g/m²

DEHS retention: 5.3 cm/s

Particle diameter: 0.3 mm

Efficiency: > 99.97 %

Material:

Filter: borosilicate fiberglass + laminated polyester

Plastic housing: polypropylene

Packing:

Tube of 10 single-use antibacterial plastic filters

Change of the filter:

Change the antibacterial filter after each patient or, for a longstay patient, after inspection of the filter clogging level.

Fixing of the filter:

Push until click-lock adjustment

Removing of the filter:

Tightly pull the filter while operating a rotational movement. Throw the filter away respecting the usual precautions.