

# ES-5400 / ES-5200

## 1.1. Multi-Channel Electrotherapy Unit 2.1.

Daugiakanalis elektroterapijos prietaisas



ES-5400



ES-5200



Electrotherapy is a treatment that promotes or represses vital reactions through the application of electric stimulation to the neuromuscular system or sensory nerves. Intended use is to reduce pain, improve muscle strength, increase joint mobility, suppress spasms, promote blood circulation and reduce the severity of edema. 28 pre-programmed therapy parameters are installed in ES-5400 and ES-5200 for over 20 types of typical pathologies. You can simply choose the area to be treated from Human Body Diagram.

### Features 4 arba 2 nepriklausomi kanalai

- 2.1.1. 1.1.1. 4 or 2 independent electrotherapy channels
- 1.1.6. 1.1.6. 28 effective clinical programs with Human Body Diagram 2.1.2.
- 1.1.2. 1.1.2. 13 types of frequently used current modes installed
- 1.1.7. 1.1.7. Sequential mode to apply two designated therapeutic courses without having to modify unit parameters 2.1.7.

- Exclusive vacuum unit available (optional)
- Constant current / constant voltage modes
- Full-color, 7" LCD touch screen 1.1.4. 2.1.4. Pilnų spalvų 7" LDC jutimasis ekranas
- 120 memorable program numbers
- 9 languages available (English, German, French, Italian, Spanish, Portuguese, Turkish, Vietnamese & Chinese)

1.1.6./2.1.6. 28 efektyvios klinikinės programos su žmogaus kūno diagrama

1.1.2./2.1.2. Instaliuota 13 dažniausiai naudojamų srovių tipų

1.1.7./2.1.7. Sekos režimas taikantis 2 pasirinktus kursus, nekeičiant prietaiso parametru

**Current Details**

1.1.2. 2.1.2.

Current Mode	Waveform image	Output mode	Maximum amplitude (peak)	Frequency
IF-4		Constant, Sweep	100 mA	1–250 Hz
IF-2		Constant, Sweep	100 mA	1–250 Hz
EMS		Surge Independent, Surge Co-Cont, Surge Alternate	100 mA	20–250 Hz
Russian		Surge Independent, Surge Co-Cont, Surge Alternate	100 mA	2.5 kHz
High Voltage (HV)		Constant, Sweep, Burst, Surge Independent, Surge Co-Cont, Surge Alternate	2.1.3. 1.1.3. <u>600 mA</u>	0.5–200 Hz
TENS		Constant, Sweep, Burst, Surge Independent, Surge Co-Cont, Surge Alternate	100 mA	0.5–250 Hz
Microcurrent (MCR)		Constant	750 µA	0.2–400 Hz
Galvanic		Continuous, Interrupted	20 mA	0.95–15.8 Hz
Faradic		Surge Independent	70 mA	20–250 Hz
Diadynamic		DF, MF, CP, CP-ISO, LP, RS	70 mA	50 Hz or 100 Hz
Traebert		Constant	70 mA	142 Hz
I/T Curve • AQ		Constant	70 mA	I/T Curve: 0.33–0.5 Hz AQ: 0.33 Hz

## Specifications

1.1.11./2.1.10.

	ES-5400	EU-5200
Power supply	AC 100–240 V, 50/60 Hz	
Power consumption	150 VA	100 VA
Number of channels	4	2
Display	Full-color, 7" LCD touch screen	
Safety class according to IEC 60601-1	Class I, Type BF	
Dimensions	345 (W) × 270 (D) × 145 (H) mm	
Weight	2.5 kg	2.3 kg
IF carrier frequency	2 kHz, 4 kHz, 5 kHz, 8 kHz, 10 kHz	
Vector sweep (for IF-4)	0°, 15°, 30°, 45°	
User-defined protocol memories	120	
Pre-programmed clinical protocols	28	
Timer	Max. 60 min.	

## Ordering Data

Standard Kit *Standard kit includes main unit.		ES-5400	ES-5200
B180534	Electrode Cable (Brown)	1x	1x
B180535	Electrode Cable (Red)	1x	1x
B180536	Electrode Cable (Orange)	1x	—
B180537	Electrode Cable (Yellow)	1x	—
011151	Rubber Electrode (M), 60 × 50 mm, 2pcs/pack	4x	2x
011148	Electrode Sponge A (M), 80 × 65 mm, 2pcs/pack	4x	2x
011655	Strap (L), 80 × 1200 mm	4x	2x
011654	Strap (S), 80 × 600 mm	4x	2x
180562	Power Supply Cord (220–240 V, Type F) or	1x	1x
180559	Power Supply Cord (110–120 V, Type A)		



Optional Accessories *Available to both ES-5400 and ES-5200	
011152	Rubber Electrode (L), 100 × 60 mm, 2pcs/pack
011150	Rubber Electrode (S), 50 × 30 mm, 2pcs/pack
011149	Electrode Sponge A (L), 120 × 80 mm, 2pcs/pack
011147	Electrode Sponge A (S), 70 × 45 mm, 2pcs/pack
010306	Self-adhesive Electrode, 49 × 49 mm, 4pcs/pack
010747	Self-adhesive Electrode, 89 × 51 mm, 4pcs/pack
010889	Self-adhesive Electrode, ø32, 4pcs/pack
011356	HV/DC Probe
011172	MCR (Microcurrent) Probe, 2pcs/pack
—	Vacuum Unit (Please see next page.)



# SU-540 / SU-520<BK>

## Blow-out System 4 or 2 Channels Vacuum Unit

Prapūtimo sistemos 4 kanalų vakuuminis prietaisas



1. SU-540 with ES-5400



SU-520<BK> with ES-5200

SU-540 for ES-5400 and SU-520<BK> for ES-5200, both models adopt blow-out vacuum system, and therefore no water reservoir is needed. It is virtually maintenance-free, since dust or moisture does not adhere inside the tube or on the connections, resulting in preventing blockages and oxidation.

### Specifications

1.1.11. / 2.1.10.

	SU-540	SU-520 <BK>
Power supply	<u>AC 100–240 V, 50/60 Hz</u>	AC 110, 120 or 220–240 V, 50/60 Hz
Power consumption	150 VA	110 and 120 V: 70 VA 220–240 V: 90 VA
Number of channels	4 independent	2 independent
Dimensions	345 (W) × 315 (D) × 115 (H) mm	345 (W) × 267 (D) × 116 (H) mm
Weight	8.5 kg	8 kg
Suction pressure	0 to -100 mmHg	

### Ordering Data

Standard Kit *Standard kit includes main unit.		SU-540	SU-520 <BK>	
011771	Electrode Hose (Blue)	2x	1x	
011772	Electrode Hose (Gray)	2x	1x	
012399	Vacuum Electrode B (S), ø80 (Gray)	4x	2x	
012400	Vacuum Electrode B (S), ø80 (Blue)	4x	2x	
011276	Electrode Sponge B (S), ø70, 4pcs/pack	4x	2x	
180562	Power Supply Cord (220–240 V, Type F) or	1x	1x	
180559	Power Supply Cord (110–120 V, Type A)			

### Optional Accessories \*Available to both SU-540 and SU-520 <BK>

012401	Vacuum Electrode B (L), ø100 (Gray)	012401	012402	011277	220279	
012402	Vacuum Electrode B (L), ø100 (Blue)					
120974	Suction Cup (L), ø100 (Gray) (cup only)					
120975	Suction Cup (L), ø100 (Blue) (cup only)					
120976	Suction Cup (S), ø80 (Gray) (cup only)					
120977	Suction Cup (S), ø80 (Blue) (cup only)					
151157	Suction Steel Plate (L), ø55			151157	151083	
151083	Suction Steel Plate (S), ø45					
121359	Suction Head			121359	121360	
121360	Suction Head Cup					
011277	Electrode Sponge B (L), ø90, 4pcs/pack					
220278	Paper Disk (L), ø90, 100pcs/pack					
220279	Paper Disk (S), ø70, 100pcs/pack					